

Procool Thumb Restriction Splint

Whilst providing full finger function, this very low profile wrap around support provides adjustable support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Osteoarthritis
- + Rheumatoid Arthritis
- + Basal joint arthritis
- + CMC Joint Arthritis
- + Juvenile Chronic Arthritis
- + CMC joint instability
- + Gamekeeper's thumb
- + Inflamed/swollen tendons
- + Overuse injuries
- + Post carpal/metacarpal surgery
- + Repetitive strain injury
- + Rhizarthrosis
- + Scapholunate instability
- + Scapholunate tear
- + Skier's thumb
- + Soft tissue injury
- + Sprained wrist
- + Stener lesion
- + Tendonitis
- + Thumb basal joint instability
- + Thumb basal joint/saddle joint irritation
- + UCL injury
- + UCL lesion
- + Weakness of the thumb joint
- + Worn cartilage





Made with tiny perforations in the Neoprene for comfortable wear and a terry cloth lining that feels good against the skin.



Available in black or beige.



Procool offers you just the fit you want - excellent stability without excessive bulk. This brace can be trimmed if required.

Code (Left Black)	Code (Right Black)	Code (Left Beige)	Code (Right Beige)	Size
PTRS/11	PTRS/21	-	-	Infant
PTRS/12	PTRS/22	-	-	Child
PTRS/13	PTRS/23	PTRSB/13	PTRSB/23	S
PTRS/13+	PTRS/23+	PTRSB/13+	PTRSB/23+	S+
PTRS/15	PTRS/25	PTRSB/15	PTRSB/25	M
PTRS/15+	PTRS/25+	PTRSB/15+	PTRSB/25+	M+
PTRS/17	PTRS/27	PTRSB/17	PTRSB/27	L
PTRS/17+	PTRS/27+	PTRSB/17+	PTRSB/27+	L+
PTRS/18	PTRS/28	PTRSB/18	PTRSB/28	XL