

Standard Wrist Brace

An 18cm wrist brace made of quality woven elastic, with easy open overstitched hook/loop tabs designed for the patient with limited dexterity.

The unique overlapped edges offer excellent comfort especially to minimise chafing at the thumb web and palmar crease. The aluminium palmar bar is adjustable to any wrist position. Also available in black.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Osteoarthritis
- + Rheumatoid arthritis
- + CMC joint arthritis
- + Basal joint arthritis
- + Juvenile chronic arthritis
- + Carpal tunnel syndrome
- + De Quervain's syndrome
- + Gamekeeper's thumb
- + Guyon's canal syndrome
- + Inflamed tendons
- + Luxation
- + Medial subluxation
- + Repetitive strain injury
- + Skier's thumb
- + Sprained wrist
- + Stener lesion
- + Thumb basal/saddle joint irritation
- + Ulnar tunnel syndrome
- + Weakness of the thumb joint
- + Worn cartilage





Quality woven elastic, with easy open overstitched hook/loop tabs designed for the patient with limited dexterity. Available in black or beige colour options.



Convenient hook / loop closures make it easy to adjust the support, and are ideal for users with limited dexterity.



The malleable palmar bar has a spoon shaped end which conforms to the palm of the hand for a comfortable fit.

Code (Left Beige)	Code (Right Beige)	Code (Left Black)	Code (Right Black)	Size	Wrist Circumference
SWB/12	SWB/22	SWB/12/BK	SWB/22/BK	XS	11.5 - 14cm
SWB/13	SWB/23	SWB/13/BK	SWB/23/BK	S	15 - 17cm
SWB/15	SWB/25	SWB/15/BK	SWB/25/BK	M	17 - 19cm
SWB/17	SWB/27	SWB/17/BK	SWB/27/BK	L	19 - 22cm
SWB/18	SWB/28	SWB/18/BK	SWB/28/BK	XL	23 - 25cm
SWB/19	SWB/29	SWB/19/BK	SWB/29/BK	XXL	25cm >