

# Sling Strap




## NEW SUSTAINABLE PACKAGING

The universal, adjustable non-stretch Sling Strap is designed to provide support and immobilisation to the upper limb. It is padded with 5cm wide material that is combined with hook Y-tabs to provide an adjustable and versatile sling system.

Sling Strap can be cut to the required length to produce a simple collar and cuff, or similar type of sling. This reduces waste and the adjustability of the material ensures that the sling will fit and can be easily adjusted.

Now available in new **sustainable, eco-friendly packaging** using no harmful dyes or inks, reducing waste and promoting a greener, more sustainable future.

- + Dislocation
- + Strains and sprains
- + Pre and post operatively
- + Rehab stages (post fractures or soft tissue injuries)
- + Brachial plexus or peripheral nerve



| Code | Length x Width | Quantity | Y-Tabs | NHSSC Code |
|------|----------------|----------|--------|------------|
| SS   | 12m x 5cm      | 2 Rolls  | 28     | GRZ554     |