

Wrist/Ulnar Deviation Support

A unique functional splint to correct early MCP Ulnar Deviation (U.D.) deformity and provide wrist support for patients with Rheumatoid arthritis. Made in ventilated Neoprene with a towelling to reduce heat build up and reduce perspiration. Incorporating a malleable aluminium palmar bar which can be shaped to fit the wrist.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + MCP Ulnar Deviation
- + Rheumatoid Osteoarthritis





Breathable fabric helps reduce heat build-up, the towelling liner provides a fantastic level of comfort



This wrist support features a malleable aluminium palmar bar that provides rigidity.



The easy to use hook and loop fastenings ensure a secure, comfortable fit.

Code (Left)	Code (Right)	Size
UDSA/12	UDSA/22	XS
UDSA/13	UDSA/23	S
UDSA/15	UDSA/25	M
UDSA/17	UDSA/27	L
UDSA/18	UDSA/28	XL