


Wrist Thumb Brace

The Wrist Thumb Brace holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint, which is crucial for alleviating pain and assisting the healing process. This thumb support also incorporates a padded aluminium thumb stay and is made from a gently elasticated material that provides compression to the injured area, reducing swelling.


Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management
- + Inflamed tendons
- + Worn cartilage
- + Tendinosis
- + Skier's thumb
- + Wrist trauma & many more





Made from Neoprene material for ultimate comfort. Compression applies gentle pressure to the area around the injury, this helps to improve circulation which prevents excess swelling.



Incorporates a padded aluminium thumb stay for great protection and support, and the aluminium bar through the palm keeps the wrist stable.

Code (Left)	Code (Right)	Size
WTN/12	WTN/22	XS
WTN/13	WTN/23	S
WTN/15	WTN/25	M
WTN/17	WTN/27	L
WTN/18	WTN/28	XL