

# Rigid Post-Op Shoe

A padded durable nylon upper with return loop buckles. The upper will accommodate swelling and dressings.

Padded insole with fully rigid lightweight composite outer sole unit. The sole unit is rigid to prevent any flexion following foot surgery that requires immobilisation to aid healing and rehabilitation.

- + Accommodate bandages or swelling
- + Foot Surgery
- + Broken toe management
- + Metatarsal injuries
- + Soft tissue procedures
- + Management of diabetic foot
- + Post-operative recovery



**Durable  
hard-wearing  
sole unit**



**Repositionable  
hook tabs**

Size	Code - Men's	Sole Length	Sole Width	Code - Women's	Sole Length	Sole Width
Paediatric	-	-	-	<b>RPO/PAED</b>	20cm	8cm
S	<b>RPO/MS</b>	28cm	10cm	<b>RPO/WS</b>	24cm	8cm
M	<b>RPO/MM</b>	29cm	10cm	<b>RPO/MM</b>	25cm	9cm
L	<b>RPO/ML</b>	30cm	10cm	<b>RPO/WL</b>	27cm	9cm
XL	<b>RPO/MXL</b>	32cm	11cm	-	-	-