

Reversible Wrist Brace

A universal left or right design. Easily converted by moving the aluminium palmar bar. The unique elastic design offers excellent compression and conforms well around the hand and wrist.

The additional elastic wrist wrap offers improved compression and wrist stabilisation.

Provides wrist support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Post trauma
- + Post fracture management



The universal design allows for use on left or right wrists

Code	Size	Wrist Circ.
RWB/02	XS	> 15cm
RWB/03	S	15cm - 17cm
RWB/05	M	17cm - 19cm
RWB/07	L	19cm - 22cm
RWB/08	XL	23cm - 25cm