

Standard Neobrace Wrist Brace

Made in Neoprene which retains heat to soothe painful joints, this brace conforms well and can be cut to fit swollen, arthritic joints without fraying. Easy open Wrioughton tabs designed for patients with limited dexterity.

Available in short (18cm) or long (23cm) versions.

- + Arthritis
- + Sprains & strains
- + Soft tissue injuries & conditions
- + Instability
- + Post fracture management



Size	Left (Short)	Right (Short)	Left (Long)	Right (Long)	Wrist Circumference
XS	SNB/12/BX5	SNB/22/BX5	LNB/12	LNB/22	11.5cm - 13cm
S	SNB/13/BX5	SNB/23/BX5	LNB/13	LNB/23	14cm - 15cm
M	SNB/15/BX5	SNB/25/BX5	LNB/15	LNB/25	16.5cm - 18cm
L	SNB/17/BX5	SNB/27/BX5	LNB/17	LNB/27	20cm - 21.5cm
XL	SNB/18/BX5	SNB/28/BX5	LNB/18	LNB/28	23cm - 24cm
XXL	SNB/19/BX5	SNB/29/BX5	-	-	25 - 28cm