

Standard Wrist Brace

An 18cm wrist brace made of quality woven elastic, with easy open overstitched hook/loop tabs designed for patients with limited dexterity. The unique overlapped edges offer excellent comfort especially to minimise chafing at the thumb web and palmar crease. The aluminium palmar bar is adjustable to any wrist position.

Available in beige or black.

- + Arthritis
- + Sprains & strains
- + Soft tissue injuries & conditions
- + Instability
- + Post fracture management



Size	Left (Beige)	Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	SWB/12/BX6	SWB/22/BX6	SWB/12/BK	SWB/22/BK	11.5cm - 13cm
S	SWB/13/BX6	SWB/23/BX6	SWB/13/BK	SWB/23/BK	14cm - 15cm
M	SWB/15/BX6	SWB/25/BX6	SWB/15/BK	SWB/25/BK	16.5cm - 18cm
L	SWB/17/BX6	SWB/27/BX6	SWB/17/BK	SWB/27/BK	20cm - 21.5cm
XL	SWB/18/BX6	SWB/28/BX6	SWB/18/BK	SWB/28/BK	23cm - 24cm