

Ulnar Deviation Splint

A unique functional splint to correct early MCP Ulnar Deviation (U.D.) deformity and provide wrist support for patients with Rheumatoid arthritis. Made in ventilated Neoprene with a towelling to reduce heat build up and reduce perspiration. Incorporating a malleable aluminium palmar bar which can be shaped to fit the wrist.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + MCP Ulnar Deviation
- + Rheumatoid Osteoarthritis



Size	Left	Right	Wrist Circumference
XS	UDSA/12	UDSA/22	11.5cm - 13cm
S	UDSA/13	UDSA/23	14cm - 15cm
M	UDSA/15	UDSA/25	16.5cm - 18cm
L	UDSA/17	UDSA/27	20cm - 21.5cm
XL	UDSA/18	UDSA/28	23cm - 24cm