

# Universal Toe Adjustment Splint

The Promedics Toe Adjustment Splint is designed to help maintain correct toe positioning following surgery or when managing conditions such as hallux valgus, hammer toe, tailor's bunion and capsulitis.

It features an elasticated metatarsal band and soft toe straps that gently guide the toe into alignment. The low-profile design fits comfortably inside most shoes, and the non-latex material helps reduce the risk of skin irritation. Suitable for use on either foot, it is supplied as a single, universal size.

For treatment of:

- + Hammer Toe
- + Tailor's Bunion
- + Hallux Valgus
- + Capsulitis



Size	Code
Universal	UTA/10