




Symmetrical Back Panel

The Vista II collar features an innovative symmetrical back panel that integrates with the VRTX CTLSO for a cervical to sacrum solution.



Vista II Cervical Collar


Incorporates the trusted design of the original Vista Cervical Collar, including a one size adjustable fit and dial height adjustment.



Vista II Cervical Collar & Vista II CTLSO

The Vista II Cervical Collar integrates the proven technology of the original Vista Cervical Collar with the reversible back panel of the VRTX Collar, ensuring precise and consistent application. This not only optimizes application times but also proves invaluable in critical situations, promoting efficiency and reliability in clinical settings.

This comprehensive off the-shelf cervical to sacrum solution enhances the overall efficiency of hospital operations.



Symmetrical Back Panel

The Vista II collar features an innovative symmetrical back panel that streamlines fitting and reduces the risk of incorrect patient application. The symmetrical posterior panel serves as a visual cue, ensuring that the collar is applied accurately and securely.



Connecting Straps

Adjustable straps easily connects the Vista II collar with the Vista II CTLSO upgrade kit for multilevel spinal support.



Cervical and Sternal Adjustability

Easy multi step cervical and sternal adjustability allows clinician to provide the care and comfort needed for the patient.


VISTA II COLLAR

Vista II Collar	984300
Vista II Collar Set with Replacement Pads	984302
Vista II Collar Replacement Pads	984320

VISTA II CTLSO STANDALONE

	Waist Circumference	
Vista II CTLSO - Standard	31 - 60 in / 78.7 – 152.4 cm	984680
Vista II CTLSO - XS	26 - 36 in / 66 - 91.5 cm	984690
Vista II CTLSO - XL	51 - 70 in / 129.5 - 177.8 cm	984700

Anterior Panel Height



Posterior Panel Height

