

# Elastic Wrist Thumb Brace

An elastic wrist brace incorporating a padded malleable aluminium thumb stay. Holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint.

- + Arthritis
- + Sprains & strains
- + Soft tissue injuries & conditions
- + Instability
- + Post fracture management



Size	Left	Right	Wrist Circumference
XS	<b>WTE/12</b>	<b>WTE/22</b>	11.5cm - 13cm
S	<b>WTE/13</b>	<b>WTE/23</b>	14cm - 15cm
M	<b>WTE/15</b>	<b>WTE/25</b>	16.5cm - 18cm
L	<b>WTE/17</b>	<b>WTE/27</b>	20cm - 21.5cm
XL	<b>WTE/18</b>	<b>WTE/28</b>	23cm - 24cm