

Neoprene Wrist Thumb Brace

Features a malleable aluminium palmar bar that provides fantastic support and stabilisation for the wrist and thumb. This brace also features a specially formulated Neoprene fabric that is ideal for orthotic braces because it provides the perfect combination of insulation and compression, offering relief for arthritic conditions and reducing swelling by improving circulation.

Provides wrist and thumb support for; an array of conditions and injuries, including;

- + Arthritis
- + CMC joint arthritis
- + Carpal tunnel syndrome
- + De Quervain's syndrome
- + Degenerative irritation
- + Inflamed tendons
- + Soft tissue injury / wrist trauma



Supports both the thumb MPC joint and restricts movement at the CMC joint, allowing you to go about your day in comfort and reduced pain.

Wrist Circ.	Code - Left	Code - Right	Wrist Circ.
XS	WTO/12	WTO/22	11.5cm - 13cm
S	WTO/13	WTO/23	14cm - 15cm
M	WTO/15	WTO/25	16.5cm - 18cm
L	WTO/17	WTO/27	20cm - 21.5cm
XL	WTO/18	WTO/28	23cm - 24cm
XXL	WTO/19	WTO/29	24cm >