

# Wrist Thumb Wrap

Holds the thumb lightly in abduction and stabilises the wrist with a double wrist wrap in Neoprene. Universal design fits left and right.

Provides wrist and thumb support for; an array of conditions and injuries, including;

- + Arthritis
- + Sprains & strains
- + Soft tissue injuries & conditions
- + Instability
- + Post fracture management



Size	Code	Wrist Circ.
S	<b>WTW/03</b>	14cm - 15cm
M	<b>WTW/05</b>	16.5cm - 18cm
L	<b>WTW/07</b>	20cm - 22cm
XL	<b>WTW/08</b>	22cm - 24cm