

Product Catalogue

A wide range of orthopaedic products

www.promedics.co.uk



Ordering Information

By Telephone:

01475 746 400

Our customer service department is open from 8.30am to 5.00pm Monday to Friday. There is a voicemail service available at all other times for orders and enquiries.

By Post:

Post details of your order, including product codes and a delivery and invoice address.

Promedics Orthopaedics Ltd
Port Glasgow Industrial Estate
Gareloch Road, Port Glasgow
Scotland, PA14 5XH

By Email:

Orders can be emailed to:
orders@promedics.co.uk

Online:

See online for our complete range of products, training courses, videos and web shop.

www.promedics.co.uk

Enquiries:

If you have an order or product enquiry, or would like to order a product catalogue, telephone our customer service department on:

01475 746 400

Delivery (UK Mainland Only)

We endeavour to deliver all orders within 48 hours, with the majority being despatched on the same day that purchase orders are received.

Delivery Dates for made-to-measure or specially requested products will be advised on receipt of an official purchase order. Please note, made to measure orders must be received in writing by email or post and once acknowledged cannot be cancelled and are subject to additional charges.

Promedics reserve the right to make partial deliveries unless otherwise instructed. All orders are subject to a standard post and packaging charge of only £5.95. Saturday delivery and timed delivery incur additional charges.

Payment


All prices are exclusive of VAT. If you have an account, payment should be made within 30 days of the date of invoice. Payment can be made by credit card, most major cards are accepted. Payment can also be made by Debit card. Cheques should be made payable to 'Promedics Orthopaedics Ltd'

Returns Policy

All goods for return must have authorisation. Please call or email our customer service department for a returns note: all returns must be notified within 28 days of receipt of goods. A restocking fee may be applied. All returned goods must be as new and in original packaging.



Contents

Introduction	2-6
Hand Therapy	7-38
Spinal Support 	39-80
Upper Limb & Shoulder	81-91
Lower Limb	92-103
Ankle & Footwear	104-127
• NEW Carbon Fibre AFO's	112-113
Maternity Support	128-131
Fracture Bracing & Splinting	132-141
Hot & Cold Therapy	142-144
Paediatric	145-154
Index	155-161
Terms & Conditions	162-166



Featuring **NEW** Carbon Fibre AFO's

Find them on page 112-113





Welcome

to the Promedics Orthopaedics Product Catalogue

Promedics are now an Aspen Medical Products company and the exclusive provider of Aspen Products in the U.K. We are very proud to present to you our catalogue which includes the Aspen product range and is packed with new and innovative products. The majority have been designed and manufactured by us in our 40,000 square foot facility in the UK.

Promedics remain absolutely committed to using British employees to manufacture our products. We employ over 150 people in the UK to ensure we continue to deliver to you, high quality products that have been designed specifically for the UK medical market and your patients needs.

The vast majority of our products are despatched on the same day that we receive your order. Our commitment to manufacturing our goods in the UK means that we continue to offer unrivalled levels of service to our customers in the UK and are consistently a top supplier to the NHS.


Your business is very important to us, and for over 35 years we have remained committed to offering exceptional service, quality and educational support to healthcare professionals. Again this year Promedics will run a full training and education workshop programme for events throughout the UK on a wide range of our orthopaedic, orthotic and therapy products.

If you would like any further information or help, then please get in touch with our Customer Service Team, your local Commercial Manager or contact us through our website at www.promedics.co.uk





Manufacturing, Partnership & Development



Promedics constantly listen to the needs of our customers, not only in the UK but also internationally where we sell to over 40 countries. For over 30 years we have invested in British manufacturing, design and production and our UK facility continues to grow and be a leader in European Orthotics and Orthopaedics.

Manufacturing in the UK not only gives us the ability to react quickly to the changing demands and needs of our customers, but it also allows us to control the raw materials and final quality of our products to ensure that we adhere to the highest possible standards.

Our facility is fully accredited to the Medical Device Directive ISO 13485:2016 and is fully audited by British Standards Institute.

We are committed to Research & Development in the UK. Our in-house team work independently and in-association with our partners and the NHS, to constantly strive to integrate the latest textiles, designs and benefits to ensure our products are designed specifically to the needs of our customers in each country we work in. Everyone is different and it is our job to meet your needs.

Promedics are very proud to represent on an exclusive basis, the world class companies of Aspen, North Coast Medical and Orfit.

 **ascender[®]**



NEW Telescopic Range of Motion Knee Brace

Find on page 93



Customer Service



Promedics operate throughout the UK and are able to supply every hospital, department and Trust in the UK on a next day basis if required. We are a key supplier to NHS Supply Chains in England, Scotland, Northern Ireland and Wales. So whether you prefer to order direct or via contract, we are there to help.

We have a full dedicated Customer Service department to ensure we provide you with outstanding, first class service every time. Whatever your needs our team are available to help, answer your call, and to ensure you get the product you are looking for. We have over 15,000 product lines in stock and our team are fully trained and waiting to help you.

We also have a full national team of local Commercial Managers throughout the UK. To get in touch with yours please contact our Customer Service Team.



Export

Promedics export to over 40 countries and we can offer you a full design and manufacture service to give you a range of products that work for your country.

We recognise that each country has its own unique demands and for many years have been a leading supplier throughout USA, Germany, Scandinavia and Western Europe.

We would be happy to work with you in Great Britain, to design and manufacture the Orthopaedic range that's right for you.



Contact Information



Tel: 01475 746 400



Email: enquiries@promedics.co.uk



Fax: 01475 746 410



@promedics



Hand Therapy



Standard Wrist Brace

Product Code: SWB

An 18cm wrist brace made of quality woven elastic, with easy open overstitched hook/loop tabs designed for the patient with limited dexterity. The unique overlocked edges offer excellent comfort especially to minimise chafing at the thumb web and palmar crease.

The aluminium palmar bar is adjustable to any wrist position.

Also available in black.

Contains Latex



Size	Left (Beige)	NHSCC Code	Right (Beige)	NHSCC Code	Left (Black)	Right (Black)	Wrist Circumference
XS	*12/BX6	GQV2226	*22/BX6	GQV2224	*12/BK	*22/BK	11.5 – 13 cm
S	*13/BX6	GQV2227	*23/BX6	GQV2225	*13/BK	*23/BK	14 – 15 cm
M	*15/BX6	GQV2221	*25/BX6	GQV2229	*15/BK	*25/BK	16.5 – 18 cm
L	*17/BX6	GQV2222	*27/BX6	GQV2230	*17/BK	*27/BK	20 – 21.5 cm
XL	*18/BX6	GQV2223	*28/BX6	GQV2228	*18/BK	*28/BK	23 – 24 cm

* Insert product code

(Supplied in boxes of 6)



Beta Wrist Brace

Product Code: BWB

An 18cm wrist brace made of a dual compression elastic which conforms to the wrist. The vinyl pocket encloses an aluminium palmar bar which can be shaped to the wrist. Lycra binding finishes the edges for comfort and prevents chafing at the thumb web.

Also available in beige.

Size	Left (Beige)	NHSCC Code (Beige)	Right (Beige)	NHSCC Code (Beige)	Left (Black)	Right (Black)	Left (Webbing)	NHSCC Code	Right (Webbing)	NHSCC Code	Wrist Circumference
XS	*12/BX6	GQV2195	*22/BX6	GQV2197	*12/BK	*22/BK	*12/W	GRZ2314	*22/W	-	11.5 – 13 cm
S	*13/BX6	GQB032	*23/BX6	GQB019	*13/BK	*23/BK	*13/W	GRZ2315	*23/W	GRZ2320	14 – 15 cm
M	*15/BX6	GQB016	*25/BX6	GQB015	*15/BK	*25/BK	*15/W	GRZ2316	*25/W	GRZ2321	16.5 – 18 cm
L	*17/BX6	GQB018	*27/BX6	GQB017	*17/BK	*27/BK	*17/W	GRZ2317	*27/W	GRZ2322	20 – 21.5 cm
XL	*18/BX6	GQV2196	*28/BX6	GQV2198	*18/BK	*28/BK	*18/W	GRZ2318	*28/W	GRZ2323	23 – 24 cm

* Insert product code

(Supplied in boxes of 6)



Nuform Wrist Brace

Product Code: NUF

An 18cm wrist brace made of a ventilated elastic for cool comfort. The vinyl pocket encloses an aluminium palmar bar which can be shaped to the wrist.

Lycra binding finishes the splint edges for a comfortable fit and prevents chafing at the thumb web.

Also available in black.

Size	Left (Beige)	NHSCC Code Left (Beige)	Right (Beige)	NHSCC Code Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*12/BX6	GQV105	*22/BX6	GQV107	*12/BK	*22/BK	11.5 – 13 cm
S	*13/BX6	GQB060	*23/BX6	GQB057	*13/BK	*23/BK	14 – 15 cm
M	*15/BX6	GQB045	*25/BX6	GQB033	*15/BK	*25/BK	16.5 – 18 cm
L	*17/BX6	GQB059	*27/BX6	GQB052	*17/BK	*27/BK	20 – 21.5 cm
XL	*18/BX6	GQV106	*28/BX6	GQV108	*18/BK	*28/BK	23 – 24 cm

* Insert product code

(Supplied in boxes of 6)



Deltaform Wrist Brace

Product Code: DFS

An 18cm woven elastic/cotton support, with an aluminium stay inside a vinyl pocket and a longitudinal seam for excellent conformity to the wrist and forearm.

Also available in beige.

Contains Latex



Size	Left (Beige)	NHSCC Code (Beige)	Right (Beige)	NHSCC Code (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*12/BX6	GQV109	*22/BX6	GQV112	*12/BK	*22/BK	11.5 – 13 cm
S	*13/BX6	GQB066	*23/BX6	GQB067	*13/BK	*23/BK	14 – 15 cm
M	*15/BX6	GQB050	*25/BX6	GQB048	*15/BK	*25/BK	16.5 – 18 cm
L	*17/BX6	GQB058	*27/BX6	GQB055	*17/BK	*27/BK	20 – 21.5 cm
XL	*18/BX6	GQV111	*28/BX6	GQV113	*18/BK	*28/BK	23 – 24 cm

* Insert product code

(Supplied in boxes of 6)



Deltaform PLUS Wrist Brace

Product Code: DFP

A longer, 23cm woven elastic/cotton support, with an aluminium stay inside a vinyl pocket and a longitudinal seam for excellent conformity to the wrist and forearm.

Extends to 2/3 of the forearm for radial/ulnar stability.

Also available in black.

Contains Latex

Size	Left (Beige)	Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*12	*22	*12/BK	*22/BK	11.5 – 13 cm
S	*13	*23	*13/BK	*23/BK	14 – 15 cm
M	*15	*25	*15/BK	*25/BK	16.5 – 18 cm
L	*17	*27	*17/BK	*27/BK	20 – 21.5 cm
XL	*18	*28	*18/BK	*28/BK	23 – 24 cm

* Insert product code



Freedom Wrist Brace

Product Code: FWB

An 18cm long brace with cut-away design to allow flexion at the palmar crease and opposition of the thumb to the fingers. The soft Lycra binding prevents chafing of the thumb web.

Longitudinal seam gives good conformity and support to the wrist and forearm.

Also available in black.

Contains Latex

Size	Left (Beige)	Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*12/BX6	*22/BX6	*12/BK	*22/BK	11.5 – 13 cm
S	*13/BX6	*23/BX6	*13/BK	*23/BK	14 – 15 cm
M	*15/BX6	*25/BX6	*15/BK	*25/BK	16.5 – 18 cm
L	*17/BX6	*27/BX6	*17/BK	*27/BK	20 – 21.5 cm
XL	*18/BX6	*28/BX6	*18/BK	*28/BK	23 – 24 cm

* Insert product code

(Supplied in boxes of 6)

 **Freedom**



Freedom SX Wrist Brace

Product Code: FSX

Using latex free elastic, the Freedom SX brace is manufactured with a full length seam assembly that shapes the brace to fit the wrist and forearm perfectly. This newly developed elastic feels soft and kind to the skin. 18cm long with reinforced casing and aluminium palmar bar.

Available in grey with black (Pictured) Also available in black.

Size	Left (Grey)	Right (Grey)	Left (Black)	Right (Black)	Wrist Circumference
XS	*12/GB	*22/GB	*12/BK	*22/BK	11.5 – 13 cm
S	*13/GB	*23/GB	*13/BK	*23/BK	14 – 15 cm
M	*15/GB	*25/GB	*15/BK	*25/BK	16.5 – 18 cm
L	*17/GB	*27/GB	*17/BK	*27/BK	20 – 21.5 cm
XL	*18/GB	*28/GB	*18/BK	*28/BK	23 – 24 cm

* Insert product code

(Supplied in boxes of 6)



Economy Standard Neobrace

Product Code: ESNB (Thin) ESNBW (3mm Thick)

An 18cm wrist brace made in Neoprene which retains heat to soothe painful joints, the brace conforms well and can be cut to fit swollen arthritic joints without fraying. Made in thinner Neoprene for light support.

Size	Left	NHSCC Code	Right	NHSCC Code	Wrist Circ.
XS	ESNB/12/BX5	GRZ852	ESNB/22/BX5	GRZ857	11.5 – 13 cm
S	ESNB/13/BX5	GRZ853	ESNB/23/BX5	GRZ858	14 – 15 cm
M	ESNB/15/BX5	GRZ854	ESNB/25/BX5	GRZ859	16.5 – 18 cm
L	ESNB/17/BX5	GRZ855	ESNB/27/BX5	GRZ860	20 – 21.5 cm
XL	ESNB/18/BX5	GRZ856	ESNB/28/BX5	GRZ861	23 – 24 cm

(Supplied in boxes of 5)

Size	Left	Right	Wrist Circumference
XS	ESNBW/12	ESNBW/22	11.5 – 13 cm
S	ESNBW/13	ESNBW/23	14 – 15 cm
M	ESNBW/15	ESNBW/25	16.5 – 18 cm
L	ESNBW/17	ESNBW/27	20 – 21.5 cm
XL	ESNBW/18	ESNBW/28	23 – 24 cm

poroflex

Poroflex is a new three dimensional fabric offering excellent properties in the design of orthopaedic braces. Poroflex offers improved patient comfort due to breathability and temperature control.



Poroflex Pull

Product Code: PFWB

- Simple "slip-on" design with two strap D ring closures. Excellent for patients with limited dexterity
- Adjustable aluminium palmar bar with flexible dorsal stays
- Excellent wrist compression and stability

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code




Poroflex Pull & Thumb

Product Code: PFTB

- Simple "slip-on" design with two strap D ring closures. Excellent for patients with limited dexterity
- Thumb extension for MCP & CMC control
- Adjustable aluminium thumb stay
- Adjustable aluminium palmar bar with flexible dorsal stays
- Excellent wrist compression and stability

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



Lycrafleece™

Lycrafleece™ is a Lycra® Spandex™ 3 ply laminate, designed for the ultimate in patient comfort. The unique construction makes this material perfect for orthopaedic bracing. It is our most comfortable and durable material which makes this range of products the elite in medical bracing.

Lycrafleece™ has unique properties which keep the patient warm when it's cold and cool when it's warm for optimum comfort without sacrificing performance.

Lycra® Spandex™ gives the material unique stretch and recovery properties whilst a polyester lining material wicks moisture away from the skin to allow greater patient compliance.

Lycrafleece Thumb Spica

Product Code: BLTS

Supports the wrist and thumb, leaving the palm completely free. An 11cm removable flexible stay restricts movement at the CMC and MCP joints, leaving the IP free. The angled hook/loop attachment at the wrist and thumb ensures a perfect fit. Use for patients with Rheumatoid and Osteoarthritis.

Lycrafleece™

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
XS	*/12	-	*/22	-	11.5 – 13 cm
S	*/13	-	*/23	-	14 – 15 cm
M	*/15	GRZ1776	*/25	GRZ1781	16.5 – 18 cm
L	*/17	-	*/27	-	20 – 21.5 cm
XL	*/18	-	*/28	-	23 – 24 cm

* Insert product code

Lycrafleece Wrist Brace

Product Code: BLWB

An 18cm three-piece wrist brace featuring longitudinal seams to give greater conformity at the wrist. The malleable palmar bar has a spoon shaped end which conforms to the palm of the hand for a comfortable fit.

Lycrafleece™

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



Pro-rheuma Wrist Brace

Product Code: PRHW (Beige) PRHWB (Black)

Made from a breathable high compression material with microfleece lining. Easy 'slip-on' design with D-ring closure straps. This brace is very comfortable and excellent for patients with limited dexterity, such as osteoarthritis and wrist instability. Supplied with medium stability plastic palmar support which can be heat moulded.

Available in black & beige. (NHSSC Beige only)

Size	Left	NHSCC Code	Right	NHSCC Code	Wrist Circumference
XS	*/12	-	*/22	-	13 – 14 cm
S	*/13	-	*/23	-	14 – 16 cm
M	*/15	GRZ5368	*/25	GRZ5341	16 – 19 cm
L	*/17	-	*/27	GRZ5342	19 – 20 cm
XL	*/18	-	*/28	-	20 – 22 cm
XXL	*/19	-	*/29	-	22 – 24 cm

* Insert product code

Pro-rheuma Wrist Thumb Brace

Product Code: **PRHT** (Beige) **PRHTB** (Black)

Made from a breathable high compression material with microfleece lining. Easy 'slip-on' design with D-ring closure straps. This brace is very comfortable and excellent for patients with limited dexterity, such as osteoarthritis and wrist instability. Supplied with medium stability plastic palmar support which can be heat moulded. Semi-flexible stay provides thumb abduction and support to CMC and MCP joints. Available in black & beige.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	13 – 14 cm
S	*/13	*/23	14 – 16 cm
M	*/15	*/25	16 – 19 cm
L	*/17	*/27	19 – 20 cm
XL	*/18	*/28	20 – 22 cm
XXL	*/19	*/29	22 – 24 cm

* Insert product code



Jura Latex Free Wrist Brace

Product Code: **JLFW**

A latex free elastic brace with horizontal wrist seam. This, together with tri-weave compression elastic, provides improved fit and compression. 18cm long with vinyl casing and aluminium palmar bar.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



Neobrace Wrist Brace

Product Code: **SNB** (Standard) **LNB** (Long)

Made in Neoprene which retains heat to soothe painful joints, this brace conforms well and can be cut to fit swollen, arthritic joints without fraying. Standard: 18cm length. Long: 23cm length.

Size	Left	NHSCC Code	Right	NHSCC Code	Wrist Circ.
XS	SNB/12/BX5	GQV2207	SNB/22/BX5	GQV2219	11.5 – 13 cm
S	SNB/13/BX5	GQV2208	SNB/23/BX5	GQV2220	14 – 15 cm
M	SNB/15/BX5	GQV2209	SNB/25/BX5	GQV2214	16.5 – 18 cm
L	SNB/17/BX5	GQV2210	SNB/27/BX5	GQV2215	20 – 21.5 cm
XL	SNB/18/BX5	GQV2211	SNB/28/BX5	GQV2216	23 – 24 cm
XXL	SNB/19/BX5	GQV2218	SNB/29/BX5	GQV2217	25 – 28 cm

(SNB Supplied in boxes of 5)

Size	Left	NHSCC Code	Right	NHSCC Codes	Wrist Circumference
XS	LNB/12	GQV557	LNB/22	GQV559	11.5 – 13 cm
S	LNB/13	GQV2205	LNB/23	GQV2203	14 – 15 cm
M	LNB/15	GQV2206	LNB/25	GQV2204	16.5 – 18 cm
L	LNB/17	GQV2200	LNB/27	GQV2212	20 – 21.5 cm
XL	LNB/18	GQV558	LNB/28	GQV560	23 – 24 cm





Air X Wrist Thumb

Product Code: AWT

A two-piece wrist thumb brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Aluminium palmar bar and thumb stay provides excellent support to wrist and first CMC & MCP joints. The thumb stay can be easily adjusted for the required degree of abduction and the thumb closure is easily adjusted for a comfortable fit

- Breathable
- Vapour transfer for reduced moisture
- Skin Friendly
- Comfort, compression & Conforms

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



Air X Thumb Restriction Splint


Product Code: ATRS

A low profile AirX wrap that allows uninterrupted finger function and provides adjustable support to the thumb CMC joint. Excellent in the treatment of osteoarthritis and medial subluxation. The adjustable strap can be wrapped around the first CMC and MCP joints to provide targeted joint support.

Size	Left	Right	Wrist Circumference
Infant	*/11	*/21	11 – 14 cm
Child	*/12	*/22	14 – 17 cm
S	*/13	*/23	15 – 18 cm
S+	*/13+	*/23+	17 – 18 cm
M	*/15	*/25	18 – 20 cm
M+	*/15+	*/25+	20 – 21 cm
L	*/17	*/27	20 – 23 cm
L+	*/17+	*/27+	21 – 24 cm
XL	*/18	*/28	23 – 25 cm

* Insert product code

- Breathable
- Comfortable
- Skin Friendly



Moisture-transporting material

The unique AirX™ fabric has the ability to transport moisture efficiently away from the body. The layering of fibre structures promotes a moisture wicking process, which results in a dry comfortable feeling to the surface of the skin so that body always feels fresh and dry. AirX™ gives the body optimal freedom to breathe – reduces the risk of skin maceration which results in improved patient compliance.

A pleasant surprise

It is hard not to be impressed. AirX™ is made of a soft, lightweight material that perfectly follows the contours of the body. This means that AirX™ adapts itself perfectly to the body so that it feels extremely comfortable to wear – even for long periods at a time – a feature, which is often a problem with orthopaedic and sport medical braces and supports. In this way, AirX™ is a comfortable material that helps to increased patient compliance without you even knowing!

AirX™
Let your body breathe

Air X One Piece Wrist

Product Code: AWB

A one piece wrist brace utilising the excellent properties of AirX to provide a comfortable well fitting brace. The textile is very durable and, together with the reinforced bar casing, makes this brace an excellent choice for patients using braces long-term.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



Air X Two Piece Wrist

Product Code: AFS

This two piece brace is shaped to achieve excellent fit, compression and stability around the hand, wrist and forearm. The 3D structure of the AirX textile gives a high degree of cushioning and increases comfort and compression. It also allows air transfer making the brace comfortable to wear for extended periods

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



Air X Thumb Spica

Product Code: ATS

Adjustable thumb Spica offering support around the wrist, MCP and CMC joints. AirX offers improved properties

- Breathable
- Vapour transfer for reduced moisture
- Skin Friendly
- Comfort, compression & Conforms

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
S	*/13	GRZ5313	*/23	GRZ5315	14 – 15 cm
M	*/15	GRZ5362	*/25	GRZ5316	15 – 16.5 cm
L	*/17	GRZ5314	*/27	GRZ5317	16.5 – 18 cm
XL	*/18	-	*/28	-	18 – 19 cm

* Insert product code





Jura Wrist Brace

Product Code: JWB

Dual compression elastic for improved support and fit. A two-piece brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Hard wearing webbing bar casing with aluminium palmar bar. Hook/loop fastening

Size	Left	Right	Wrist Circumference
XS	*12	*22	11.5 – 13 cm
S	*13	*23	14 – 15 cm
M	*15	*25	16.5 – 18 cm
L	*17	*27	20 – 21.5 cm
XL	*18	*28	23 – 24 cm

* Insert product code



Jura Wrist Thumb

Product Code: JWT

Dual compression elastic for improved support and fit. A two-piece wrist thumb brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Aluminium palmar bar and thumb stay provides excellent support to wrist and first CMC & MCP joints. The thumb stay can be easily adjusted for the required degree of abduction and the thumb closure is easily adjusted for a comfortable fit.

Size	Left	Right	Wrist Circumference
XS	*12	*22	11.5 – 13 cm
S	*13	*23	14 – 15 cm
M	*15	*25	16.5 – 18 cm
L	*17	*27	20 – 21.5 cm
XL	*18	*28	23 – 24 cm

* Insert product code



Jura Wrist Thumb (Long)

Product Code: JWTL

Similar design and function as the standard Jura Wrist Thumb Brace. Increasing length to 23cm offers improved wrist stability and immobilisation.

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
XS	*12		*22	-	11.5 – 13 cm
S	*13	GQT396	*23	-	14 – 15 cm
M	*15	GQT397	*25	GQT402	16.5 – 18 cm
L	*17	GQT398	*27	GQT403	20 – 21.5 cm
XL	*18	GQT399	*28	GQT404	23 – 24 cm

* Insert product code



Thumbs Up

Product Code: TUJ

Lightweight structure and comfortable fitting thumb spica. Innovative design with breathable holes. Protection glove and contact closure make this thumb spica an easy to use, well fitting solution for thumb care. Anatomically shaped, provides the right immobilisation to relieve pain during everyday activities.

Contains Latex

Size	Left	NHSCC Code	Right	NHSCC Code	Wrist Circumference
S	*/13	GRZ1981	*/23	GRZ1871	13 – 17 cm
M	*/15	GRZ1982	*/25	GRZ1872	17 – 19 cm
L	*/17	GRZ1983	*/27	GRZ1873	19 – 23 cm
XL	*/18	GRZ1870	*/28	GRZ1874	23 – 28 cm

* Insert product code



Jura Palm Protector

Product Code: JPC

Fleece lined lightweight palm cover. Designed to allow full finger movement whilst reducing finger contractures and skin damage/breakdown in the palmar region. The Jura Palm Cover is easily laundered and is easily applied using hook/loop fastening. One size fits all.

Code (Left)	NHSCC Codes	Code (Right)	NHSCC Codes
JPC/L	GQT198	JPC/R	GQT199

(Supplied in packs of 3)



Intrinsic Wrist Brace

Product Code: IWB

Padded elastic wrist brace with a full-length adjustable aluminium resting splint. The splint is wide enough to contain all four fingers and can be easily adjusted to hold the wrist and fingers in the desired position. There is a wide, padded adjustable dorsal closure strap that can be used to hold the fingers in position.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



Wrightington Neobrace

Product Code: WNB

An 18cm wrist brace with easy-open, over-stitched hook/loop tabs designed for the patient with limited dexterity.

Made in Neoprene which retains heat to soothe painful joints, it conforms well and can be cut to fit swollen arthritic joints without fraying.

Size	Left	Right	Wrist Circumference
XS	*/12/BX5	*/22/BX5	11.5 – 13 cm
S	*/13/BX5	*/23/BX5	14 – 15 cm
M	*/15/BX5	*/25/BX5	16.5 – 18 cm
L	*/17/BX5	*/27/BX5	20 – 21.5 cm
XL	*/18/BX5	*/28/BX5	23 – 24 cm

* Insert product code

(Supplied in boxes of 5)



Ventoprene Neobrace

Product Code: VNB

An 18cm wrist brace made in ventilated Neoprene with a towelling lining to reduce perspiration build-up. Incorporating a malleable aluminium palmar bar which can be shaped to fit the wrist.

Size	Left (Beige)	NHSCC Code	Right (Beige)	NHSCC Code	Wrist Circumference
XS	*/12	GRZ2359	*/22	-	11.5 – 13 cm
S	*/13	-	*/23	-	14 – 15 cm
M	*/15	GRZ2361	*/25	GRZ2366	16.5 – 18 cm
L	*/17	-	*/27	-	20 – 21.5 cm
XL	*/18	-	*/28	-	23 – 24 cm

* Insert product code

Size	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12/BK	*/22/BK	11.5 – 13 cm
S	*/13/BK	*/23/BK	14 – 15 cm
M	*/15/BK	*/25/BK	16.5 – 18 cm
L	*/17/BK	*/27/BK	20 – 21.5 cm
XL	*/18/BK	*/28/BK	23 – 24 cm

* Insert product code



Ventoprene Wrist Brace

Product Code: VWNB (Beige) R/VWNB (Blue)

An 18cm wrist brace with Wrightington easy-open, over-stitched hook/loop tabs designed for the patient with limited dexterity.

Made in beige or blue ventilated Neoprene with a towelling lining to reduce heat build-up and perspiration. Incorporating a malleable aluminium palmar bar which can be shaped to fit the wrist.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code

Wrist/Ulnar Deviation Support

Product Code: UD^{SA}

A unique functional splint to correct early MCP Ulnar Deviation (U.D.) deformity and provide wrist support for patients with Rheumatoid arthritis. Made in ventilated Neoprene with a towelling to reduce heat build up and reduce perspiration. Incorporating a malleable aluminium palmar bar which can be shaped to fit the wrist.

Size	Left	NHSCC Code	Right	NHSCC Code	Wrist Circumference
XS	*/12	GQT1187	*/22	GQT1087	11.5 – 13 cm
S	*/13	GQT1188	*/23	GQT1082	14 – 15 cm
M	*/15	GQT1189	*/25	GQT1083	16.5 – 18 cm
L	*/17	GQT1191	*/27	GQT1084	20 – 21.5 cm
XL	*/18	GQT1086	*/28	GQT1085	23 – 24 cm

* Insert product code



Wristlet

Product Code: WL

A 15cm wrist support made in Neoprene which retains heat to soothe painful joints. It conforms well giving light compression and can be cut without fraying, to fit swollen joints.

Size	Left	Right	Wrist Circumference
S-M	*/14	*/24	14 – 18 cm
M-L	*/16	*/26	17 – 22 cm
L-XL	*/LXL/L	*/LXL/R	23 – 24 cm

* Insert product code



Elastic Wrist Thumb Brace

Product Code: WTE

An elastic wrist brace incorporating a padded malleable aluminium thumb stay.

Holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint.

Contains Latex

Size	Left	NHSCC Code	Right	NHSCC Code	Wrist Circumference
XS	*/12	GQV071	*/22	GQV101	11.5 – 13 cm
S	*/13	GQV072	*/23	GQV102	14 – 15 cm
M	*/15	GQV073	*/25	GQV070	16.5 – 18 cm
L	*/17	GQV098	*/27	GQV103	20 – 21.5 cm
XL	*/18	GQV099	*/28	GQV104	23 – 24 cm

* Insert product code





Neoprene Wrist Thumb Brace

Product Code: WTO

A malleable aluminium thumb stay holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint.

Size	Left	NHSCC Code	Right	NHSCC Code	Wrist Circumference
XS	*/12	GQT1126	*/22	GQT1137	11.5 – 13 cm
S	*/13	GQT1127	*/23	GQT1138	14 – 15 cm
M	*/15	GQT1128	*/25	GQT1133	16.5 – 18 cm
L	*/17	GQT1129	*/27	GQT1134	20 – 21.5 cm
XL	*/18	GQT1132	*/28	GQT1135	23 – 24 cm
XXL	*/19	GQT000	*/29	GQT1136	24 cm +

* Insert product code



Wrist Thumb Brace

Product Code: WTN

A Neoprene wrist brace incorporating a padded, malleable aluminium thumb stay. Holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint. The thumb section opens fully for easy application and to accommodate swelling.

Size	Left	NHSCC Code	Right	NHSCC Code	Wrist Circumference
XS	*/12	GQT1122	*/22	GQT1121	11.5 – 13 cm
S	*/13	GQT1123	*/23	GQT1124	14 – 15 cm
M	*/15	GQT1118	*/25	GQT1125	16.5 – 18 cm
L	*/17	GQT1119	*/27	GQT1130	20 – 21.5 cm
XL	*/18	GQT1120	*/28	GQT1131	23 – 24 cm

* Insert product code



Bedford Wrist Thumb Brace

Product Code: BWTB

A unique wrist immobilisation brace which can be used to support undisplaced fractures of the scaphoid, trapezium and thumb metacarpal. It features a mouldable bar which holds the thumb firmly in abduction and has been improved to reduce pressure points. It also now features an open thumb seam for ease of application and comfort.

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
XS	*/12	GQT1114	*/22	GQT1115	11.5 – 13 cm
S	*/13	GQV062	*/23	GQV067	14 – 15 cm
M	*/15	GQV063	*/25	GQV060	16.5 – 18 cm
L	*/17	GQV064	*/27	GQV068	20 – 21.5 cm
XL	*/18	GQV065	*/28	GQV069	23 – 24 cm

* Insert product code

Abducted Thumb Wrap

Product Code: ATW

A Neoprene wrist/thumb wrap with a hook/loop closure, incorporating a flexible stay to hold the thumb firmly in abduction.

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
XS	*/12	GQT1098	*/22	-	11.5 – 13 cm
S	*/13	GQT1099	*/23	GQT1111	14 – 15 cm
M	*/15	GQT1100	*/25	GQT1112	16.5 – 18 cm
L	*/17	GQT1108	*/27	GQT1113	20 – 22 cm
XL	*/18	GQT1116	-	-	23 – 24 cm

* Insert product code



Wrist Thumb Wrap

Product Code: WTW

Holds the thumb lightly in abduction and stabilises the wrist with a double wrist wrap in Neoprene.

Universal design fits left and right.

Size	Code	NHSCC Code	Wrist Circumference
S	*/03	GRZ2475	14 – 15 cm
M	*/05	GRZ2476	16.5 – 18 cm
L	*/07	GRZ2477	20 – 22 cm
XL	*/08	GRZ2478	22 – 24 cm

* Insert product code



Wrist Wrap

Product Code: WW

The double wrist wrap in Neoprene provides warmth and compression to soothe arthritic joints.

Universal design fits left and right.

Size	Code	Wrist Circumference
One Size	WW/10	13 – 21.5 cm





Elastic Thumb Spica

Product Code: ETSLF

Dual compression beige elastic supports the wrist and thumb, leaving the palm free. The 11cm removable flexible thumb stay restricts movement at the CMC and MCP joints, leaving the IP joint free. Latex Free.

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
XS	ETSLF/12	GQT1088	ETSLF/22	GQT1093	11 – 13 cm
S	ETSLF/13	GQT050	ETSLF/23	GQT054	14 – 15 cm
M	ETSLF/15	GQT051	ETSLF/25	GQT049	15 – 16.5 cm
L	ETSLF/17	GQT052	ETSLF/27	GQT055	16.5 – 18 cm
XL	ETSLF/18	GQT053	ETSLF/28	GQT056	18 – 19 cm



Elastic Thumb Spica – Modified

Product Code: ETSMLF

One-way stretch elastic supports the wrist and thumb, leaving the palm completely free. An 10cm removable flexible thumb stay restricts movement at the CMC and MCP joints, leaving the IP joint free. Angled hook/loop attachment at wrist and thumb ensures perfect fit. Latex Free.

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
XS	*/12	GQT1094	*/22	GQT1095	11 – 13 cm
S	*/13	GQT1089	*/23	GQT1101	14 – 15 cm
M	*/15	GQT1090	*/25	GQT1102	16 – 17 cm
L	*/17	GQT1091	*/27	GQT1096	17 – 18 cm
XL	*/18	GQT1092	*/28	GQT1097	18 – 19 cm

* Insert product code



Promedics Thumb Shield

Product Code: TS

A Neoprene wrist/thumb wrap with adjustable cinch strap and removable thermoplastic insert. The cinch strap allows greater compression at the site of the pain, without complete immobilisation.

The thermoplastic insert, made from Orfit NS 2.0mm, increases the support at the MCP joint and can be removed and remoulded as required.

Size	Code	NHSCC Code	Wrist Circumference
XS	*/02	GRZ1976	11 – 13 cm
S	*/03	GRZ1977	14 – 16 cm
M	*/05	GRZ1978	17 – 19 cm
L	*/07	GRZ1979	20 – 22 cm
XL	*/08	GRZ1980	23 – 25 cm

* Insert product code

Pollo Form Spica – Standard

Product Code: SPF

The unique two-way stretch material gives the Pollo Form excellent properties for supporting the wrist and thumb. Outer material is durable and supports an inner pocket with a pre-cut low temperature thermoplastic insert. The insert is activated by pouring warm water into the pocket. After removing the water the brace is then fitted to the patient whilst forming the insert to hold the thumb in the required position.

The insert can be reactivated to allow serial adjustments or removed later in the treatment programme. The Long version gives greater control around the wrist.

Size	Left	Right	Wrist Circumference
XS	*12	*22	Up to 14 cm
S	*13	*23	14 – 16 cm
M	*15	*25	16 – 19 cm
L	*17	*27	19 – 21 cm
XL	*18	*28	21 – 24 cm

* Insert product code



Pollo Form Spica – Long

Product Code: LPF

Details above.

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
XS	*12	GRZ2457	*22	GRZ2462	Up to 14 cm
S	*13	GRZ2458	*23	GRZ2463	14 – 16 cm
M	*15	GRZ2459	*25	GRZ2464	16 – 19 cm
L	*17	GRZ2460	*27	-	19 – 21 cm
XL	*18	-	*28	-	21 – 24 cm

* Insert product code



Pro-rheuma Thumb Spica

Product Code: PRHS

Made from a breathable high compression material with microfleece lining. Removable semi-flexible thumb stay helps stabilise CMC & MCP joints. Simple hook and loop closures.

Size	Left	Right	Wrist Circumference
XS	*12	*22	13 – 14 cm
S	*13	*23	14 – 16 cm
M	*15	*25	16 – 18 cm
L	*17	*27	18 – 20 cm
XL	*18	*28	20 – 22 cm

* Insert product code





Jura Ulnar Deviation Support

Product Code: JUDS

Helps support the MCP joints and reduce ulnar drift whilst allowing functional use of the fingers. Manufactured from neoprene to provide warmth and support

Size	Left	NHSCC Codes	Right	NHSCC Codes	Width at MCP Joints
S – M	*/14	GRZ1736	*/24	GRZ1737	6 – 9 cm
M – L	*/16	GRZ1738	*/26	GRZ1739	9 – 11 cm

* Insert product code

Jura Ulnar Deviation Support Extended

Product Code: JUDSE

The extended version of the popular Jura Ulnar Deviation Splint wraps around the thumb, providing a secure and beneficial fit. It supports the MCP joints and helps to reduce ulnar drift, whilst allowing functional use of the fingers. The splint is manufactured from neoprene to provide warmth and support. Supports and maintains neutral alignment of the MCP joints in rheumatoid arthritis.

Product Features

- Finger tabs provide ulnar-to-radial pull
- Easy to apply and remove
- Low profile to provide patient comfort

Size	Left	Right	Width at MCP Joints
S – M	*/14	*/24	6 – 9 cm
M – L	*/16	*/26	9 – 11 cm

* Insert product code

Juraprene Wrist Thumb Wrap

Product Code: JTWS

A latex free neoprene with terry plush lining makes this a very comfortable support. Full finger flexion whilst giving good circumferential compression and support to the wrist and thumb.

The wrap is easily fastened using hook and loop closures.

Fits either left or right hand and can be trimmed in areas for a custom fit. Available in Beige, Black and Blue.

Size	Colour	Code	Wrist Circumference
S	Beige	*/03/BE	Up to 17 cm
M	Beige	*/05/BE	17 – 19 cm
L	Beige	*/07/BE	19 cm +
S	Blue	*/03/BL	Up to 17 cm
M	Blue	*/05/BL	17 – 19 cm
L	Blue	*/07/BL	19 cm +
S	Black	*/03/BK	Up to 17 cm
M	Black	*/05/BK	17 – 19 cm
L	Black	*/07/BK	19 cm +

* Insert product code



Fits either left or right hand and can be trimmed in areas for a custom fit. Available in Beige, Black and Blue .

Juraprene Wrist Thumb Wrap – Long

Product Code: JTWWL

A latex free neoprene with terry plush lining makes this a very comfortable support. Full finger flexion whilst giving good circumferential compression and support to the wrist and thumb. This longer wrist wrap adds to the compression and firmness around the wrist and is easily fastened using hook and loop closures.

Size	Colour	Code	Wrist Circumference
S	Beige	*/03/BE	Up to 17 cm
M	Beige	*/05/BE	17 – 19 cm
L	Beige	*/07/BE	19 cm +
S	Blue	*/03/BL	Up to 17 cm
M	Blue	*/05/BL	17 – 19 cm
L	Blue	*/07/BL	19 cm +
S	Black	*/03/BK	Up to 17 cm
M	Black	*/05/BK	17 – 19 cm
L	Black	*/07/BK	19 cm +

* Insert product code



Juraprene Short Wrist Wrap – Universal

Product Code: JSWWU

A latex free neoprene with terry plush lining makes this a very comfortable universal support. Full finger and thumb flexion whilst giving good circumferential compression and support to the wrist.

Size	Colour	Code	NHSCC Codes
Universal	Beige	*/10/BE	-
Universal	Blue	*/10/BL	-
Universal	Black	*/10/BK	GRZ2235

* Insert product code



Juraprene Long Wrist Wrap – Sizes


Product Code: JLWW

A latex free neoprene with terry plush lining makes this a very comfortable support. Full finger and thumb flexion whilst wider reinforced hook and loop closure gives improved circumferential compression and support to the wrist.

Available in Beige, Black and Blue.

Size	Colour	Code	Wrist Circumference
S	Beige	*/03/BE	Up to 17 cm
M	Beige	*/05/BE	17 – 19 cm
L	Beige	*/07/BE	19 cm +
S	Blue	*/03/BL	Up to 17 cm
M	Blue	*/05/BL	17 – 19 cm
L	Blue	*/07/BL	19 cm +
S	Black	*/03/BK	Up to 17 cm
M	Black	*/05/BK	17 – 19 cm
L	Black	*/07/BK	19 cm +

* Insert product code



Procool

Procool is a range of lightweight and breathable splints producing excellent compression and support for a wide range of conditions and injuries.

Made with tiny perforations in the Neoprene for comfortable wear and a terry cloth lining that feels good against the skin.

Procool offers you just the fit you want – excellent stability without excessive bulk.

Procool Thumb Restriction Splint

Product Code: PTRS (Black) PTRSB (Beige)

Whilst providing full finger function, this very low profile wrap around support provides adjustable support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support.

An extensive size range covers paediatric to adult sizes and plus sizes are a full size bigger in thumb circumference to fit oedematous thumb conditions. Function and fit without bulk.

Size	Left	Right	Wrist Circumference
Infant	*/11	*/21	11 – 14 cm
Child	*/12	*/22	14 – 17 cm
S	*/13	*/23	15 – 18 cm
S+	*/13+	*/23+	17 – 18 cm
M	*/15	*/25	18 – 20 cm
M+	*/15+	*/25+	20 – 21 cm
L	*/17	*/27	20 – 23 cm
L+	*/17+	*/27+	21 – 24 cm
XL	*/18	*/28	23 – 25 cm

Infant & child sizes
not available in
beige

* Insert product code

Procool

This brace can be trimmed if required. Available in black and beige.

Indications Include:

- CMC Joint Arthritis
- CMC Instability
- Tendinitis, Pain associated with repetitive motion.



Procool

Indications Include:

- Arthritic conditions
- Post reconstructive surgery
- Sprains / strains
- Tendinitis
- Post wrist fracture

Procool Wrist Support

Product Code: PWS

This support provides effective circumferential compression and support to the wrist whilst allowing full thumb and finger motion.

Useful as a step-down splint in later stage rehab or to encourage function with patients who struggle with transition from immobilisation to movement due to pain or anxiety.

The brace can be trimmed if required.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



This brace can be trimmed if required.

Indications Include:

- Sprains
- Arthritis
- Carpal Tunnel Syndrome
- Strains



Indications Include:

- Post scaphoid fracture
- De Quervian's
- Post surgical repair
- Carpal Tunnel Syndrome
- Skier's Thumb
- Tendinitis



Indications Include:

- Post scaphoid fracture
- De Quervian's
- Post surgical repair
- Carpal Tunnel Syndrome
- Skier's Thumb
- Tendinitis

Procool D-ring Wrist Splint

Product Code: **PDR** (Standard) **PDRF** (Firm)

Low profile support with circumferential slip-on design and D-Ring fastening to allow easier one handed application. Aluminium palmar bar can be removed for further adjustment. Dorsal stockinette lining facilitates application and cushions strap closures. Design of the Standard and Firm braces allow full finger and thumb function whilst limiting wrist flexion.

The Firm splint has the addition of adjustable dorsal stays to reduce unwanted flexion and extension and improve overall immobilisation of acute soft tissue injuries.

Size	Left (Std)	NHSCC Code	Right (Std)	NHSCC	Left (Firm)	Right (Firm)	Wrist Circ.
XS	*12	GRZ2257	*22	GRZ2262	*12	*22	> 15 cm
S	*13	GRZ2258	*23	GRZ2263	*13	*23	15 – 17 cm
M	*15	GRZ2259	*25	GRZ2264	*15	*25	17 – 19 cm
L	*17	GRZ2260	*27	GRZ2265	*17	*27	19 – 22 cm
XL	*18	GRZ2261	*28	GRZ2266	*18	*28	22 cm +

* Insert product code

Procool Deluxe Wrist Thumb Brace

Product Code: **PDWT**

A uniquely shaped wrist thumb brace offering excellent support and compression to the wrist, CMC and MCP joints.

Malleable palmar bar positions wrist in functional extension and flexible stay provides abduction and limited motion to support and protect CMC & MCP joints. Adjustable thumb closure allows easy application of the brace.

Size	Left	Right	Wrist Circumference
XS	*12	*22	12 – 14 cm
S	*13	*23	14 – 16 cm
M	*15	*25	16 – 18 cm
L	*17	*27	18 – 20 cm
XL	*18	*28	20 – 22 cm

* Insert product code

Procool Wrist Thumb Restriction

Product Code: **PWTRS**

This extended version of the thumb restriction splint provides additional wrist support for a combination of wrist and thumb conditions. A wide 'pull through' wrist strap allows easier brace application, adjustment and fastening, offering extra compression and support at the wrist. Function and fit without bulk.

Size	Left	Right	Wrist Circumference
XS	*12	*22	12 – 14 cm
S	*13	*23	15 – 18 cm
M	*15	*25	18 – 20 cm
L	*17	*27	20 – 24 cm
XL	*18	*28	23 – 25 cm

* Insert product code


Jura Black Wrist Brace

Product Code: **JRWB** (Regular) **JRWBL** (Long Version)

The Jura Black Wrist Brace immobilises and protects the wrist following tendonitis and other soft tissue injuries such as sprains/strains, and post-trauma healing. Lightweight, breathable and contoured for a comfortable fit. Slide-on design with D-ring closures allow easy application.

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
XS	*/12	GRZ891	*/22	GRZ896	> 15 cm
S	*/13	GRZ892	*/23	GRZ897	15 – 17 cm
M	*/15	GRZ893	*/25	GRZ898	17 – 19 cm
L	*/17	GRZ894	*/27	GRZ899	19 – 22 cm
XL	*/18	GRZ895	*/28	GRZ900	23 – 25 cm

* Insert product code



Long version available (JRWBL)



Jura Black Wrist / Thumb Brace

Product Code: **JRWTB**

The Jura Black Wrist/Thumb Brace immobilises and protects the wrist and thumb following tendonitis and other soft tissue injuries such as sprains/strains, and post-trauma healing. Lightweight, breathable and contoured for a comfortable fit. Slide-on design with D-ring closures allow easy application.

Size	Left	NHSCC Codes	Right	NHSCC Codes	Wrist Circumference
XS	*/12	GMB095	*/22	GMB100	> 15 cm
S	*/13	GMB096	*/23	GMB101	15 – 17 cm
M	*/15	GMB097	*/25	GMB102	17 – 19 cm
L	*/17	GMB098	*/27	GMB103	19 – 22 cm
XL	*/18	GMB099	*/28	GMB104	23 – 25 cm

* Insert product code





Ventilated Thumb Restriction Splint

Product Code: TRSV

A ventilated, low profile wrap around support providing adjustable support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support.

An extensive size range covers paediatric to adult sizes and plus sizes are a full size bigger in thumb circumference to fit oedematous thumb conditions. Function and fit without bulk.

Size	Left	Right	Wrist Circumference
XXS	*/11	*/21	11 – 14 cm
XS	*/12	*/22	14 – 17 cm
S	*/13	*/23	15 – 18 cm
S+	*/13+	*/23+	17 – 18 cm
M	*/15	*/25	18 – 20 cm
M+	*/15+	*/25+	20 – 21 cm
L	*/17	*/27	20 – 23 cm
L+	*/17+	*/27+	21 – 24 cm
XL	*/18	*/28	23 – 25 cm

* Insert product code

Wrist & Thumb Restriction Splint Ventilated

Product Code: WTRSV

A ventilated, low profile wrap around support providing adjustable support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 - 13 cm
S	*/13	*/23	14 - 15 cm
M	*/15	*/25	16.5 - 18 cm
L	*/17	*/27	20 - 21.5 cm
XL	*/18	*/28	23 - 24 cm

* Insert product code



Reversible Wrist

Product Code: RWB

A universal left or right design. Easily converted by moving the aluminium palmar bar. The unique elastic design offers excellent compression and conforms well around the hand and wrist. The additional elastic wrist wrap offers improved compression and wrist stabilisation.

Contains Latex

Size	Code	NHSCC Code	Wrist Circumference
XS	*/02	GQV517	11.5 – 13 cm
S	*/03	GQV518	14 – 15 cm
M	*/05	GQV519	16.5 – 18 cm
L	*/07	GQV520	20 – 21.5 cm
XL	*/08	GQV521	23 – 24 cm

* Insert product code



Jura Universal Wrist Brace

Product Code: JUWB (Standard) JUWBU (Unbound)

A universal support that is easily adjusted for a good fit. Made from a comfortable ventilated material that offers compression and support.

- Sprains / strains
- Carpal Tunnel
- Post Trauma
- Arthritic pain

Size	Code	NHSCC Code (Standard Only)	Code
Left	JUWB/L	GRZ901	JUWBU/L
Right	JUWB/R	GRZ902	JUWBU/R

Mini Thumb Wrap

Product Code: MTW

A one piece low profile wrap offering mild compression and support to the thumb CMC and MCP.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16 – 17 cm
L	*/17	*/27	17 – 18 cm
XL	*/18	*/28	18 – 19 cm

* Insert product code



Two Piece Neoprene Spica

Product Code: TPNS

A unique two piece construction offers excellent support and fit around the wrist and thumb.

Mouldable stay provides support to the thumb CMC and MCP.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16 – 17 cm
L	*/17	*/27	17 – 18 cm
XL	*/18	*/28	18 – 19 cm

* Insert product code



Neobrace PLUS Wrist Brace

Product Code: NNWB

This enhanced version of the Standard Neobrace offers easy grip-tabs to the hook closures. Bound edges and a hard wearing palmar bar casing ensures long term wear. Adjustable aluminium palmar bar.

Also available in black.

Size	Left (Beige)	Right (Beige)	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12	*/22	*/12/BK	*/22/BK	11 – 13 cm
S	*/13	*/23	*/13/BK	*/23/BK	14 – 15 cm
M	*/15	*/25	*/15/BK	*/15/BK	16 – 17 cm
L	*/17	*/27	*/17/BK	*/17/BK	17 – 18 cm
XL	*/18	*/28	*/18/BK	*/18/BK	18 – 19 cm

* Insert product code





Air X Low Profile Spica

Product Code: LPS

A low profile thumb spica offering targeted support to the thumb MCP and CMC joints.


The unique construction of the closure area offers excellent adjustment and stability around the joint.

The 3D structure of the AirX textile gives a high degree of cushioning and increases comfort and compression.

- Breathable
- Vapour transfer for reduced moisture
- Skin Friendly
- Comfort, compression & Conforms

Size	Left	Right	Wrist Circumference
XS	*12	*22	11 – 13 cm
S	*13	*23	14 – 15 cm
M	*15	*25	16 – 17 cm
L	*17	*27	17 – 18 cm
XL	*18	*28	18 – 19 cm

* Insert product code



Rainbow Putty

Rainbow Putty is a cost effective option for hand rehabilitation and treatment of a wide range of conditions and injuries.

It comes in 5 distinctive colours to indicate levels of resistance from extra soft to firm. The putty can be stretched, squeezed and pinched to provide various treatment modalities. It can also be combined to create a customised resistance.

Resistance	Product Code	Colour
X-Soft	JP01/*weight code	Orange
Soft	JP02/* weight code	Yellow
Soft-Medium	JP03/* weight code	Pink
Medium	JP04/* weight code	Green
Firm	JP05/* weight code	Blue

Weight	Weight Code
2oz (57g)	02
3oz (85g)	03
4oz (113g)	04
1lb (453g)	1
5lb (2.3kg)	5

*When ordering please quote your desired tub weight code with your order

Selected sizes available NHS Supply Chain



Jura Basic Wrist Thumb Neoprene

Product Code: JBWTN

A malleable aluminium thumb stay holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint. Length 18cm. Made from neoprene which retains heat. Can be cut without fraying.

Size	Left (Beige)	NHSCC Code	Right (Beige)	NHSCC Code	Wrist Circumference
XS	*12	-	*22	-	11.5 - 13 cm
S	*13	GRZ2554	*23	GRZ2558	14 - 15 cm
M	*15	GRZ2555	*25	GRZ2559	16.5 - 18 cm
L	*17	GRZ2556	*27	GRZ2560	20 - 21.5 cm
XL	*18	-	*28	GRZ2561	23 - 24 cm

Size	Left (Black)	Right (Black)	Wrist Circumference
XS	*12/BLK	*22/BLK	11.5 - 13 cm
S	*13/BLK	*23/BLK	14 - 15 cm
M	*15/BLK	*25/BLK	16.5 - 18 cm
L	*17/BLK	*27/BLK	20 - 21.5 cm
XL	*18/BLK	*28/BLK	23 - 24 cm

* Insert product code



Procool Wrist Thumb Support

Product Code: PWTS

An adjustable wrist thumb wrap offering good support around the wrist, MCP and CMC joints. Adjustable thumb closure with excellent shaping around the palm.

Indications Include:

- Basal thumb arthritis
- Thumb CMC & MCP arthritis
- Post UCL injury
- Degenerative joints
- Rheumatoid arthritis
- Osteo-arthritis

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	12 – 14 cm
S	*/13	*/23	14 – 16 cm
M	*/15	*/25	16 – 18 cm
L	*/17	*/27	18 – 20 cm
XL	*/18	*/28	20 – 22 cm

* Insert product code



Carbon SX Wrist Brace

Product Code: CSX

A newly developed ultra thin material that allows the CARBON SX brace to conform exactly to the contours of the wrist and hand.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



Black Seal SX Wrist Brace

Product Code: SSX

A four way stretch high compression Lycra is combined with a fully breathable microfleece to provide a material that offers great wrist support and compression. Lycrafleece is a comfortable and durable material and has unique properties that keep the patient's wrist and hand at the optimal temperature. The inner microfleece also helps wick away moisture from the skin to allow greater patient comfort.

Use CSX size chart above



Air SX Wrist Brace

Product Code: ASX

Using a 3D four way stretch fabric the Air brace is designed to offer compression, support and comfort. The Air brace offers excellent breathability and temperature control due to the open 3D construction of the textile. It is well suited to those patient's that prefer the security of a thicker more cushioned material.

The Air brace utilises the 3D material to offer a comfortable, well fitting brace with excellent breathability and overall temperature control.

Size	Left	Right	Wrist Circumference
XS	*/12	*/22	11.5 – 13 cm
S	*/13	*/23	14 – 15 cm
M	*/15	*/25	16.5 – 18 cm
L	*/17	*/27	20 – 21.5 cm
XL	*/18	*/28	23 – 24 cm

* Insert product code



Padded Liner

- **Easy Liner** – Preformed plastic shell with a towel liner that is attached using hook/loop tabs. The shell is supplied with two towel lined padded liners that can be easily removed and replaced. Removable straps can be positioned as required. This option is more suited to conditions that benefit from easily changeable liners such as burns or trauma and is also a lower cost option
- **Thermoplastic Insert** – The full padded cover functional resting splint, described as above with a shaped piece of Orfit Eco 3.2mm non perforated thermoplastic insert, which is perfect for splinting to the patients needs and can be used as a serial splint also. Splint can be shaped and at a later date taken out of the cover back into a waterbath and reshaped to match changing condition.

Functional Resting Splint

Product Code: **JRS (Padded)** **JRSE (Easy Liner)**
JRS/TP (Thermoplastic Insert)

A functional resting splint that is available in three options:

- **Full Padded Cover** – Preformed plastic shell with fully removable padded cover and comfortable wide closure straps. This is excellent for patients using the splint on a longer term basis and the fully padded cover offers excellent comfort and compliance to those patients that are very sensitive.

Size	Left (Beige)	NHSCC Codes (JRS Only)	Right (Beige)	NHSCC Codes (JRS Only)	Wrist Circumference
XS	*/12	-	*/22	-	11.5 – 13 cm
S	*/13	GRZ2173	*/23	GRZ2159	14 – 15 cm
M	*/15	GRZ2174	*/25	GRZ2160	16.5 – 18 cm
L	*/17	GRZ2175	*/27	GRZ2161	20 – 21.5 cm
XL	*/18	-	*/28	-	23 – 24 cm

* Insert product code

Jura Universal Wrist Thumb Brace

Product Code: **JUWBT**

A universal support that is easily adjusted for a good fit. Made from a comfortable ventilated material that offers compression and support.

Size	Left	NHSCC Code	Right	NHSCC Code
Universal	JUWBT/L	GRZ903	JUWBT/R	GRZ904



Jura Basic Wrist Brace Neoprene

Product Code: **JBNWB**

Made in Neoprene which retains heat to soothe painful joints, this brace conforms well and can be cut to fit swollen, arthritic joints without fraying.

Also available in black.

Size	Left (Beige)	NHSCC Code	Right (Beige)	NHSCC Code	Wrist Circumference
XS	*/12	-	*/22	-	11.5 - 13 cm
S	*/13	GRZ905	*/23	GRZ908	14 - 15 cm
M	*/15	GRZ906	*/25	GRZ2213	16.5 - 18 cm
L	*/17	GRZ907	*/27	GRZ2214	20 - 21.5 cm
XL	*/18	GRZ2211	*/28	GRZ2215	23 - 24 cm

Size	Left (Black)	Right (Black)	Wrist Circumference
XS	*/12/BLK	*/22/BLK	11.5 - 13 cm
S	*/13/BLK	*/23/BLK	14 - 15 cm
M	*/15/BLK	*/25/BLK	16.5 - 18 cm
L	*/17/BLK	*/27/BLK	20 - 21.5 cm
XL	*/18/BLK	*/28/BLK	23 - 24 cm

* Insert product code





Jura Basic Wrist Brace Elastic

Product Code: JBWB

Dual compression elastic wrist brace. Palmar bar and hook and loop closures.

Also available in black.

Size	Left (Beige)	NHSCC Code	Right (Beige)	NHSCC Code	Wrist Circumference
XS	*12	GRZ2216	*22	GRZ2221	11.5 - 13 cm
S	*13	GRZ2217	*23	GRZ2222	14 - 15 cm
M	*15	GRZ2218	*25	GRZ2223	16.5 - 18 cm
L	*17	GRZ2219	*27	GRZ2224	20 - 21.5 cm
XL	*18	GRZ2220	*28	GRZ2225	23 - 24 cm

Size	Left (Black)	Right (Black)	Wrist Circumference
XS	*12/BLK	*22/BLK	11.5 - 13 cm
S	*13/BLK	*23/BLK	14 - 15 cm
M	*15/BLK	*25/BLK	16.5 - 18 cm
L	*17/BLK	*27/BLK	20 - 21.5 cm
XL	*18/BLK	*28/BLK	23 - 24 cm

* Insert product code

Neoprene Thumb Spica

Product Code: NTS

A Neoprene thumb spica featuring a removable flexible stay which restricts movement at the CMC and MCP joints, leaving the IP joint free.

Size	Left	NHSCC Code	Right	NHSCC Code	Wrist Circumference
S	*13	GRZ1802	*23	GRZ1806	14 - 15 cm
M	*15	GRZ1803	*25	GRZ1807	15 - 16.5 cm
L	*17	GRZ1804	*27	GRZ1808	16.5 - 18 cm
XL	*18	GRZ1805	*28	GRZ1809	18 - 19 cm

* Insert product code



Carpotec Carpal Brace

Product Code: COT

Features

Malleable aluminium frame covered with perforated breathable material. The orthosis stabilises the wrist, hand and fingers and can easily be adapted to the size of the wrist and thumb of the patient. Its anatomical shape provides comfortable and effective therapy.

Size	Left	Right	Wrist Width
XS	*12	*22	< 7 cm
S	*13	*23	7 - 8 cm
M	*15	*25	8 - 9 cm
L	*17	*27	9 - 10 cm

* Insert product code



Soft, padded
hook / loop straps

Clinical Applications

Post-operative treatment, sprains, arthritis, osteoarthritis of articulation of the wrist: radial nerve palsy, tendonitis, to support the hand and wrist in cases of neurological deficits.



Breathable Material



Comfortable, removable
finger spaces



Jura CMC Thumb Brace

Product Code: JURACMC

- Ergonomic frame provides superior support and excellent comfort of CMC joint
- Inner mouldable aluminium stay offers extra support and stability of thumb
- Easy to apply and adjust with simple straps

Size	Left	NHSCC Code	Right	NHSCC Code	Circumference
S	*13	GQT481	*23	GQT484	14 – 15 cm
M	*15	GQT482	*25	GQT485	16.5 – 18 cm
L	*17	GQT483	*27	GQT486	20 – 21.5 cm

* Insert product code

Push CMC

The Push CMC consists of a cylindrical section around the ball of the thumb and a fixed connection running along the palm to the outside of the hand. This can be closed at the back of the hand with two nonelastic bands. The cylindrical section consists of two components: an aluminium strip and a flexible, synthetic part. The aluminium strip can be shaped and individually adapted for maximum stabilisation of the base of the thumb.

The selected materials and rounded edges guarantee maximum comfort. The simple closing system enables patients to fit the brace themselves without any effort.

The Push CMC is a slim and compact product. The brace can be used in the water. Hygiene is important around the hand. Therefore the synthetic material is antibacterial and the brace can be washed in a washing machine at 40°.



Judy C. Colditz

OT/L, CHT, FAOTA International hand therapy expert, President IFSHT

The Push CMC is an ingenious design that solves the problem of pain for those with CMC osteoarthritis, including me! I am eager to wear it because of its comfort, lack of functional impediment, and durability. In my experience it is the only commercially produced thumb brace that is superior to what I could custom mold for a patient.

Size	Code (Left)	NHSCC Code	Circumference
1	310111	GRZ2299	16 – 19.5 cm
2	310112	GRZ2300	19.5 – 22.5 cm
3	310113	GRZ2301	22.5 – 26 cm

Size	Code (Right)	NHSCC Code	Circumference
1	310121	GRZ2302	16 – 19.5 cm
2	310122	GRZ2303	19.5 – 22.5 cm
3	310123	GRZ2304	22.5 – 26 cm



Circumference





Single Armchair

Product Code: CAP/2

A dynamic splint with coil design for PIP and DIP joint flexion contractures and Boutonniere's deformity. Maintains assisted extension at the PIP joint while allowing active flexion and can be used for soft tissue injury affecting the PIP joint with a possible flexion deformity.

Size	Code	NHSCC Codes	Circumference at PIP
XS	*XS	GRZ814	4 – 5 cm
S	*S	GRZ812	5.5 – 6 cm
M	*M	-	6 – 6.5 cm
L	*L	GRZ811	6.5 – 7 cm
XL	*XL	GRZ813	7 – 7.5 cm

* Insert product code



Double Armchair

Product Code: CAP/3

A double armchair design for index/middle or middle/ring fingers which maintains assisted extension at the PIP joint while allowing active flexion. A dynamic splint with coil design for PIP and DIP joint flexion contractures and Boutonniere deformity. Can be used for soft tissue injury affecting the PIP joint with a possible flexion deformity.

Size	Code	Circumference at PIP
S	*S	5.5 – 6 cm
M	*M	6 – 6.5 cm
L	*L	6.5 – 7 cm
XL	*XL	7 – 7.5 cm

* Insert product code



Triple Armchair

Product Code: CAP/8

A dynamic splint with coil design for PIP and DIP joint flexion contractures which allows DIP joint to be free for active range of movement. The shorter, stronger spring supports PIP, leaving the DIP joint free to flex. Can be used as a progressive splint after intensive use of the Single Armchair Splint CAP/2.

Size	Code	Circumference at PIP
XS	*XS	5 – 5.5 cm
S	*S	5.5 – 6 cm
M	*M	6 – 6.5 cm
L	*L	6.5 – 7 cm
XL	*XL	7 – 7.5 cm
XXL	*XXL	7.5 – 8 cm

* Insert product code



Cambridge Finger Splint

Product Code: CAP/5

A dynamic splint with coil design for PIP joint flexion contractures, which allows DIP joint to be free for active range of movement. Offers the same level of support as the Simpson Armchair Splint CAP/4, providing additional support to the middle phalanx by the application of a padded strap at MP.

Size	Code	Circumference at PIP
S	*/S	6 cm
M	*/M	7 cm
L	*/L	8 cm

* Insert product code



Simpson Modified Armchair

Product Code: CAP/4

A dynamic splint with coil design for PIP and DIP joint flexion contractures which allows DIP joint to be free for active range of movement. The shorter, stronger spring supports PIP, leaving the DIP joint free to flex. Can be used as a progressive splint after intensive use of the Single Armchair Splint CAP/2.

Size	Code	NHSCC Codes	Circumference at PIP
XS	*XS	GRZ820	5 – 5.5 cm
S	*S	GRZ817	5.5 – 6 cm
M	*M	GRZ818	6 – 6.5 cm
L	*L	GRZ816	6.5 – 7 cm
XL	*XL	GRZ819	7 – 7.5 cm

* Insert product code



Plastic Mallet Splint

Product Code: MFS

Supports the DIP joint in neutral, to help correct mallet deformity, while protecting the distal phalanx and extending the PIP joint. Perforated design prevents skin maceration.

Size	Code	NHSCC Code	Circumference at PIP
S	*1/PK10	GQT1044	4.5 – 5 cm
S-M	*2/PK10	GQT1045	5 – 5.5 cm
M	*3/PK10	GQT1051	6 – 6.5 cm
M-L	*4/PK10	GQT1052	6.5 – 7 cm
L	*5/PK10	GQT1046	7 – 7.5 cm
L-XL	*5.5/PK10	GQT1047	7.5 – 8 cm
XL	*6/PK10	GQT1048	8 – 8.5 cm

* Insert product code

(Supplied in packs of 10)



Double Bedford Finger Stall

Product Code: DBF

Cotton/Elastic fabric conforms to the size of the finger. The stalls support both the PIP and DIP joints.

Contains Latex

Size	Code	NHSCC Codes	Width	Length
Paediatric	*01	GQT1041	4.5 cm	5 cm
Child	*02	GQT1042	4.5 cm	7 cm
Small	*03	GQT1161	6.25 cm	7 cm
Medium	*05	GQT1162	6.25 cm	9 cm
Large	*07	GQT1156	6.25 cm	10 cm

* Insert product code

(Supplied in packs of 10)



Whitechapel Finger Splints

Product Code: WFS

To increase range of movement by securing the injured finger next to a stronger finger ('buddy strap' effect). Use for soft tissue injury to PIP joint. Made in 12.5cm elasticated loop and plain hook.

Contains Latex

Size	Code	Length
S	*15	15 cm
M	*30	30 cm
L	*50	50 cm

* Insert product code

(Supplied in packs of 50)



Jura Basic Mallet Finger Splints

Product Code: JBPM

Supports the DIP joint in neutral, to help correct mallet deformity, while protecting the distal phalanx and extending the PIP joint.

Size	Code	NHSCC Code	Circumference at PIP
S	*1	GRZ1785	4.5 – 5 cm
S-M	*2	GRZ1786	5 – 5.5 cm
M	*3	GRZ1787	6 – 6.5 cm
M-L	*4	GRZ1788	6.5 – 7 cm
L	*5	GRZ1789	7 – 7.5 cm
L-XL	*5.5	GRZ1790	7.5 – 8 cm
XL	*6	GRZ1791	8 – 8.5 cm

* Insert product code

(Supplied in packs of 10)



Swan Neck

Oval – 8 Finger Splint

Complete Kit Code: P1008-K

Sizing Set Code: P1008-S

Made of 1.6mm high-temperature polypropylene, Oval-8 splints can be adjusted with a heat gun for a small angle alteration. With a simple turn of the splint, each Oval-8 Finger Splint can be used to correct Swan Neck, Boutonniere and Mallet deformities, reduce Lateral Deviation and rest Trigger Finger.

Oval-8 splints feature a wide, flat band for comfortable pressure distribution and control. Available in 12 sizes, each splint fits both a full and half ring size. Sizes correspond to jeweller's ring sizes, though correlation may vary due to swelling or joint changes.



Boutonniere



Complete Kit

Refills (Pack of 5)	NHSCC Code	Size	Individual	NHSCC Code	Size
P1008-2-5	GRZ1835	2	P1008-2	-	2
P1008-3-5	GRZ1836	3	P1008-3	GQT1168	3
P1008-4-5	GRZ1837	4	P1008-4	GQT1169	4
P1008-5-5	GRZ1838	5	P1008-5	GQT1163	5
P1008-6-5	GRZ1839	6	P1008-6	GQT1164	6
P1008-7-5	GQV327	7	P1008-7	GQT1165	7
P1008-8-5	GQV328	8	P1008-8	GQT1166	8
P1008-9-5	GRZ1840	9	P1008-9	GQT1167	9
P1008-10-5	GQV326	10	P1008-10	GQT1197	10
P1008-11-5	GRZ1830	11	P1008-11	GQT1157	11
P1008-12-5	GRZ1831	12	P1008-12	GQT1158	12
P1008-13-5	GRZ1832	13	P1008-13	GQT1159	13
P1008-14-5	GRZ1833	14	P1008-14	GQT1160	14

The economical Oval-8 Kit includes 42 splints in sizes 3-14, a carrying case, a ring gauge, a pen light and instructions. Packaged with most splints in the more commonly used mid-range sizes.

The Oval-8 Sizing Set includes one each of splint size 3-14 and the ring gauge.

Aluminium Finger Splints

These padded aluminium finger splints help provide protection, support and promotes healing comfortably for injuries such as trauma, breaks, cuts, dislocations and many others.

Easy to apply, simply position the splint on the finger, bend into position, trim with heavy duty scissors if required and secure with microporous tape.

Product Code	NHSCC Code	Splint Size
100302P	GRZ1771	3/4 inches
100303P	GRZ1772	1/2 inches

Supplied in boxes of 12





Spinal Support

Spine Applications – Upper Spine



C1-C7

Decompressive procedures
Fracture management
Cervicogenic headache
Cervical disc syndrome
Postsurgical stabilization
Musculoskeletal pain
Trauma / whiplash
Sprains / strains
Radiculopathy

Cervicogenic headaches
Loss of range of motion
General neck pain
Muscle spasms
Radiculopathy
Thoracic pain
Neuropathy
Whiplash
Stiff neck

C1-T2

Postsurgical stabilization
Cervicothoracic instability
Fracture management



Vista Collar



Vista TX Collar




Vista MultiPost

Aspen Collar
Aspen Pediatric CollarSierra Universal Collar
Sierra Pediatric Collar

Vista MultiPost Therapy Collar

Vista CTO4
Vista CTOAspen CTO
Aspen Pediatric CTO

Spine Applications – Lower Spine



SCOLIOSIS

Adult scoliosis with postural, pain or quality of life issues

T8-S1

Postoperative support / pain relief
Kyphosis secondary to osteoporosis
Decompressive procedures
Degenerative disc disease
Bulging or herniated disc

Fracture management
Spondylolisthesis
Spinal stenosis
Spondylosis
Severe DJD

T9-S1

Postoperative support / pain relief
Deconditioned trunk musculature
Post laminectomy syndrome
Decompressive procedures
Degenerative disc disease
Bulging or herniated disc
Insidious low back pain
Chronic low back pain
Fracture management
Spondylolisthesis

IDET procedure
Spinal stenosis
Sprains / strains
Facet syndrome
SI dysfunction
Radiculopathy
Spondylosis
Severe DJD
Sacrolitis
Injections

L1-L5

Postoperative support / pain relief
Deconditioned trunk musculature
Post laminectomy syndrome
Decompressive procedures
Degenerative disc disease
Bulging or herniated disc
Insidious low back pain
Chronic low back pain
Fracture management
Spondylolisthesis

IDET procedure
Spinal stenosis
Sprains / strains
Facet syndrome
SI dysfunction
Radiculopathy
Spondylosis
Severe DJD
Sacrolitis
Injections

S1-S5

Symphysis pubis dysfunction
Sacroiliac joint dysfunction
Chronic sacroiliac support
Low back sprains / strains
Sacroiliac joint fusion
Fracture management

Piriformis syndrome
Pelvic muscle pain
Sacroiliac pain
Low back pain
Arthritis
Sciatica



Aspen Tri-Point FSO



Vista 464 TLSO



Summit 456 TLSO



Horizon 456 TLSO



Contour TLSO



Sierra TLSO



Aspen TLSO



Summit 637 LSO



Horizon 637 LSO



Contour LSO



Aspen Sierra 637 LSO

Aspen LSO
Aspen LSO LoPro

QuikDraw PRO



Horizon 627 Lumbar



Aspen Lumbar Support



Evergreen 621 SI Belt

Aspen Vista Upper Spine



Vista Cervical Collar

Six sizes in one collar
Dial height adjustment
Effective motion restriction
Reduce inventory and waste
Improved skin care



Vista TX Collar

Thoracic extension
Extended support



Vista MultiPost Collar

Pivoting occipital panels
Designed to reduce pressure
Superior motion restriction
Fully adjustable front and back panels



Vista MultiPost Therapy Collar


Fully adjustable front and back panels
Pneumatic compression bladder
Integrated air pump
Hot/cold therapy



Vista ICU Back Panel

Extended back panel height
Dual-density foam
Pressure reducing pad
Effective motion restriction
Open by design

Vista Collar	984000
Vista Collar Set with Replacement Pads	984002
Vista Collar Replacement Pads	984020
Vista Collar Large Back Panel Replacement Pad	984028
Vista TX Collar	984100
Vista TX Collar Set with Replacement Pads	984102
Vista TX Collar Replacement Pads	984120
Vista MultiPost Collar	984200
Vista MultiPost Collar Set with Replacement Pads	984202
Vista MultiPost Collar Replacement Pads	984220
Vista MultiPost Therapy Collar	984250
Vista MultiPost Therapy Collar Replacement Pad Set	984252
Vista ICU Back Panel	984600
Vista ICU Back Panel Replacement Pad	984620
Large Back Panel	984008



C1-C7

Vista Collar • Vista TX Vista MultiPost

Decompressive procedures
Fracture management
Cervicogenic headache
Cervical disc syndrome
Postsurgical stabilization
Musculoskeletal pain
Trauma / whiplash
Sprains / strains
Radiculopathy

Vista MultiPost Therapy Collar

Cervicogenic headaches
Loss of range of motion
General neck pain
Muscle spasms
Radiculopathy
Thoracic pain
Neuropathy
Stiff neck
Whiplash

Aspen Vista Upper Spine



Vista CTO4

Maximum motion restriction
Comfort without compromising care
Supports progression of care



Vista CTO

Maximum motion restriction
Multiple adjustment points
Accommodates anatomical changes

Vista CTO4

984550

Vista CTO

984510

Vista CTO4 Upgrade Kit • Posterior


984554

Vista CTO Upgrade Kit • Anterior

984504

Vista CTO4 / Vista CTO Replacement Pads

984520




Vista CTO • Vista CTO4
Postsurgical stabilization
Cervicothoracic instability
Fracture management

Aspen Cervical Collar



Effective motion restriction
Minimize skin breakdown
Better patient outcomes

Short	983108
Regular	983110
Tall	983112
X-Tall	983114
Short Collar Set with Replacement Pads	983128
Regular Collar Set with Replacement Pads	983130
Tall Collar Set with Replacement Pads	983132
X-Tall Collar Set with Replacement Pads	983134
Small Back Panel — 30-51 cm circumference	983118
Standard Back Panel — 33-53 cm circumference	983120
Large Back Panel — 38-64 cm circumference	983119
Aspen Collar Replacement Pads	983159
Aspen Collar Small Back Panel Replacement Pad	983161



Aspen Cervical Collar
 Decompressive procedures
 Fracture management
 Cervicogenic headache
 Cervical disc syndrome
 Postsurgical stabilization
 Musculoskeletal pain
 Trauma / whiplash
 Sprains / strains
 Radiculopathy

Aspen Pediatric Collar



Effective motion restriction
Safe for delicate skin
Gentle by design
Custom fit

PD1

1-18 months
53-84 cm length
5 - 13 kg



PD2

9-24 months
74-94 cm in length
10 - 15 kg



PD3

1-3 years
84-102 cm in length
11 - 16 kg



PD4


2-5 years
89-114 cm in length
12 - 19 kg



PD5

3-6 years
94-122 cm in length
12 - 24 kg

PD1	983190
PD2	983195
PD3	983106
PD4	983102
PD5	983104
PD1 Set with Replacement Pads	983192
PD2 Set with Replacement Pads	983197
PD3 Set with Replacement Pads	983126
PD4 Set with Replacement Pads	983122
PD5 Set with Replacement Pads	983124
PD3, PD4, PD5 Back Panel — 22-33 cm circumference	983116
Aspen Pediatric Collar Replacement Pads	983199




Aspen Pediatric Collar
 Decompressive procedures
 Fracture management
 Cervicogenic headache
 Cervical disc syndrome
 Postsurgical stabilization
 Musculoskeletal pain
 Trauma / whiplash
 Sprains / strains
 Radiculopathy

Aspen CTO



66-135 cm circumference
Proven performance
Multiple orthotic options
Available in two or four post option
Available in a pediatric CTO

Short	983308
Regular	983310
Tall	983312
X-Tall	983314
PD CTO with PD4 and PD5 Front	983301
Short Front Panel	983309
Regular Front Panel	983311
Tall Front Panel	983313
X-Tall Front Panel	983315
Small Back Panel	983318
Standard Back Panel	983319
Large Back Panel	983321
Aspen CTO Replacement Pads	983320
PD CTO Replacement Pads	983302
Control Clamps (2)	983330



Aspen CTO

Postsurgical stabilization
Cervicothoracic instability
Fracture management

Aspen Sierra Universal Collar




- Universal adult size
- Universal pediatric size
- Reduce inventory
- Eliminate SKUs
- Excellent shower collar
- Available in a pediatric collar

Adult

983125

Pediatric

983150



Aspen Sierra Collar

Decompressive procedures
Fracture management
Cervicogenic headache
Cervical disc syndrome
Postsurgical stabilization
Musculoskeletal pain
Trauma / whiplash
Sprains / strains
Radiculopathy

Aspen Vista Lower Spine




Vista 464 TLSO

Three shoulder strap options
Height adjustment technology



One Size Adjustable — Fits waists ranging from 66-152 cm and up to 178 cm with the extension panel

Vista 464 TLSO	993640
Vista 464 TLSO Upgrade Kit — 35-55 cm thoracic bar height	993644
Aspen Extension Panel	992505

**T8-S1 — Vista 464**

Decompressive procedures - epidurals
Postoperative support / pain relief
Kyphosis secondary to osteoporosis
Degenerative disc disease
Bulging or herniated disc
Fracture management
Spondylolisthesis
Spinal stenosis
Spondylosis
Severe DJD

Aspen Summit



Summit 456 TLSO

Versatile adjustable posterior support
Comfortable — easy to don and doff
Thoracic to lumbar




Summit 637 LSO

Additional coverage
Anterior support
Conforms to different anatomies

One Size Adjustable — Fits waists ranging from 66-152 cm and up to 178 cm with the extension panel

Summit 456 TLSO	One Size Adjustable	992710
Summit 637 LSO	One Size Adjustable	992310
Aspen Extension Panel		992505



T8-S1 — Summit 456

Decompressive procedures - epidurals
 Postoperative support / pain relief
 Kyphosis secondary to osteoporosis
 Degenerative disc disease
 Bulging or herniated disc
 Fracture management
 Spondylolisthesis
 Spinal stenosis
 Spondylosis
 Severe DJD

T9-S1 — Summit 637

Postoperative support / pain relief
 Deconditioned trunk musculature
 Post laminectomy syndrome
 Decompressive procedures
 Degenerative disc disease
 Bulging or herniated disc
 Insidious low back pain
 Chronic low back pain
 Fracture management
 Spondylolisthesis
 IDET procedure
 Spinal stenosis
 Sprains / strains
 Facet syndrome
 SI dysfunction
 Radiculopathy
 Spondylosis
 Severe DJD
 Sacrolitis
 Injections

Aspen Horizon



Horizon 456 TLSO

Height adjustment technology
Wide circumference range



Horizon 637 LSO

Flexion and extension support
with lateral control
Integrated anterior support




Horizon 627 Lumbar

Rapid return to activity
Lowest profile design

One Size Adjustable — Fits waists ranging from 61-127 cm and up to 178 cm with the extension panel

Horizon 456 TLSO	993740
Horizon 637 LSO	993730
Horizon 627 Lumbar	993710
Horizon Extension Panels	993745



T8-S1 — Horizon 456

Decompressive procedures - epidurals
 Postoperative support / pain relief
 Kyphosis secondary to osteoporosis
 Degenerative disc disease
 Bulging or herniated disc
 Fracture management
 Spondylolisthesis
 Spinal stenosis
 Spondylosis
 Severe DJD

T9-S1 — Horizon 637

L1-L5 — Horizon 627

Postoperative support / pain relief
 Deconditioned trunk musculature
 Post laminectomy syndrome
 Decompressive procedures
 Degenerative disc disease
 Bulging or herniated disc
 Insidious low back pain
 Chronic low back pain
 Fracture management
 Spondylolisthesis
 IDET procedure
 Spinal stenosis
 Sprains / strains
 Facet syndrome
 SI dysfunction
 Radiculopathy
 Spondylosis
 Severe DJD
 Sacrolitis
 Injections

Aspen Evergreen



Evergreen 621 SI Belt
Sacral stabilization
Bi-lateral pull
Non-slip silicone strips

Small — 64-104 cm

Medium — 89-127 cm

Large — 114-152 cm



Evergreen 621 SI Belt	Small	36-2102	Medium	36-2104	Large	36-2106
-----------------------	-------	---------	--------	---------	-------	---------



S1-S5 — Evergreen 621

Symphysis pubis dysfunction
 Sacroiliac joint dysfunction
 Chronic sacroiliac support
 Low back sprains / strains
 Sacroiliac joint fusion
 Fracture management
 Piriformis syndrome
 Pelvic muscle pain
 Sacroiliac pain
 Low back pain
 Arthritis
 Sciatica

S1-S5

Aspen Contour Bracing System



Contour LSO

QuikDraw PRO with Rigid Anterior Panel (RAP)
Contour Back Panel



Contour SPK Complete

QuikDraw PRO with Rigid Anterior Panel (RAP)
Contour Back Panel
Sternal Pad Kit (SPK) Option



Contour TLSO

QuikDraw PRO with Rigid Anterior Panel (RAP)
Contour Back Panel
Shoulder Straps
Chest Panel



QuikDraw PRO


Easy effective compression
Independent upper and lower tightening

Custom Builds Available

Upon request Promedics offer a custom built Contour SPK Complete. Please provide us with the patients measurements (see "How to Measure" guide) and we will deliver a fully built custom Contour SPK Complete ready to fit your patient.

All components are available separately and are interchangeable please contact Customer Service for further details on **01475 746 400** or email enquiries@promedics.co.uk

QuikDraw PRO	Black X-Small	993081
	Black Small	993011
	Black Medium	993021
	Black Large	993031
	Black X-Large	993041
	Black XX-Large	993091
Black Regular RAP (Rigid Anterior Panel)		993200
Black Large RAP (Rigid Anterior Panel)		993201
Black Extension Panel		993061




QuikDraw PRO


Postoperative support / pain relief
 Deconditioned trunk musculature
 Post laminectomy syndrome
 Decompressive procedures
 Degenerative disc disease
 Bulging or herniated disc
 Insidious low back pain
 Chronic low back pain
 Fracture management
 Spondylolisthesis
 IDET procedure
 Spinal stenosis
 Sprains / strains
 Facet syndrome
 SI dysfunction
 Radiculopathy
 Spondylosis
 Severe DJD
 Sacrolitis
 Injections

How to Measure

- 1 Lower rib circumference = cm
- 2 Hip circumference = cm
- 3 T9 to Sacrococcygeal junction = cm
- 4 Length from Symphysis Pubis to the Sternal Notch = cm
- A Use waist circumference (average of 1 and 2 = cm)



Contour TLSO SPK Complete • Contour Back Panel • QuikDraw RAP • SPK Option	Medium – 79 - 94cm	994351	X-Large – 104 - 119cm	994353
	Large – 91 - 107cm	994352	XX-Large – 117 - 145cm	994354
Contour TLSO SPK Upgrades	Small			994455
	Medium			994460
	Large			994465
Contour TLSO Upgrades	Medium Chest Panel			994260
	Short Shoulder Strap			994151
	Tall Shoulder Strap			994152



T8-S1 – Contour TLSO


Decompressive procedures - epidurals
Postoperative support / pain relief
Kyphosis secondary to osteoporosis
Degenerative disc disease
Bulging or herniated disc
Fracture management
Spondylolisthesis
Spinal stenosis
Spondylosis
Severe DJD



Contour Back Panel	12" Back Panel	992212
	15" Back Panel	992215






Determining Sizes

Use "How to Measure" guide on pg. 65

1. Determine Quikdraw™ Size		
Size		Code
<input type="checkbox"/> X-Small 53 - 69 cm		993081
<input type="checkbox"/> Small 66 - 81 cm		993011
<input type="checkbox"/> Medium 79 - 94 cm		993021
<input type="checkbox"/> Large 91 - 107 cm		993031
<input type="checkbox"/> X-Large 104 - 119 cm		993041
<input type="checkbox"/> XX-Large 117 - 145 cm		993091

2. Determine RAP Size		
Size		Code
<input type="checkbox"/> Regular		993200
<input type="checkbox"/> Large		993201

3. Determine Chest Panel Size		
Size		Code
<input type="checkbox"/> Small < 33 cm		994255
<input type="checkbox"/> Medium 28 - 41 cm		994260
<input type="checkbox"/> Large > 33 cm		994265

4. Determine Contour™ Back Panel Size		
Size		Code
<input type="checkbox"/> Contour < 38 cm		992212
<input type="checkbox"/> Contour < 47 cm		992215

Aspen Sierra



One Size Adjustable — Fits waists ranging from 71-137 cm and up to 188 cm with the extension panel

Aspen Sierra Adjustable TLSO
71-137 cm

501221

Aspen Sierra Adjustable LSO
71-137 cm


501111

Aspen Sierra Adjustable Extension Panel
51 cm

501322

Aspen Sierra Adjustable Upgrade Kit
Standard t-bar: 33-46 cm
Tall t-bar: 41-53 cm

501321



T8-S1 — Aspen Sierra Adjustable TLSO

Decompressive procedures - epidurals
Postoperative support / pain relief
Kyphosis secondary to osteoporosis
Degenerative disc disease
Bulging or herniated disc
Fracture management
Spondylolisthesis
Spinal stenosis
Spondylosis
Severe DJD

T9-S1 — Aspen Sierra Adjustable LSO

Postoperative support / pain relief
Deconditioned trunk musculature
Post laminectomy syndrome
Decompressive procedures
Degenerative disc disease
Bulging or herniated disc
Insidious low back pain
Chronic low back pain
Fracture management
Spondylolisthesis
IDET procedure
Spinal stenosis
Sprains / strains
Facet syndrome
SI dysfunction
Radiculopathy
Spondylosis
Severe DJD
Sacrolitis
Injections

Aspen TLSO • Aspen LSO



Aspen TLSO

Additional options:

- Full TLSO
 - Over the shoulder
 - Hyperextension
- Multiple orthotic options

Short / Small — 66-99 cm

33 cm Back Panel, 20 cm Cut-Out

Short / Large — 91 - 124 cm

33 cm Back Panel, 20 cm Cut-Out

Tall / Small — 81-119 cm

40 cm Back Panel, 25 cm Cut-Out

Tall / Large — 112-150 cm

40 cm Back Panel, 25 cm Cut-Out



Aspen LSO

Superior motion restriction

Customized support

Large lumbar window

Easy LeverLatch™ Closure

Rigid layer of support

Short / Small — 66-99 cm

33 cm Back Panel, 20 cm Cut-Out

Short / Large — 91-124 cm

33 cm Back Panel, 20 cm Cut-Out

Tall / Small — 81-119 cm

40 cm Back Panel, 25 cm Cut-Out

Tall / Large — 112-150 cm

40 cm Back Panel, 25 cm Cut-Out



Aspen LSO LoPro

Low-profile lateral panels

Small — 66-99 cm

40 cm Back Panel, 25 cm Cut-Out


Large — 91-124 cm

40 cm Back Panel, 25 cm Cut-Out

X-Large — 117-158 cm

40 cm Back Panel, 25 cm Cut-Out

Aspen TLSO	Short / Small	991030	Tall / Small	991070
	Short / Large	991031	Tall / Large	991071
Aspen LSO	Short / Small	990040	Tall / Small	990060
	Short / Large	990041	Tall / Large	990061
Aspen LSO LoPro	Small	990020	X-Large	990022
	Large	990021		
Aspen LSO / TLSO Replacement Pads				991035
Aspen LSO LoPro Replacement Pads				990025
Aspen LSO LoPro X-Large Replacement Pads				990027



T8-S1 — Aspen TLSO

Decompressive procedures - epidurals
 Postoperative support / pain relief
 Kyphosis secondary to osteoporosis
 Degenerative disc disease
 Bulging or herniated disc
 Fracture management
 Spondylolisthesis
 Spinal stenosis
 Spondylosis
 Severe DJD

T9-S1 — Aspen LSO • Aspen LSO LoPro

Postoperative support / pain relief
 Deconditioned trunk musculature
 Post laminectomy syndrome
 Decompressive procedures
 Degenerative disc disease
 Bulging or herniated disc
 Insidious low back pain
 Chronic low back pain
 Fracture management
 Spondylolisthesis
 IDET procedure
 Spinal stenosis
 Sprains / strains
 Facet syndrome
 SI dysfunction
 Radiculopathy
 Spondylosis
 Severe DJD
 Sacrolitis
 Injections

Aspen Lumbar Support




Simple, elegant design
Two-stage relief
Essential trunk stability
Independent Compression
Tailored, comfortable fit
Easy to locate pull tabs

X-Small — 53-69 cm
Small — 66-79 cm

Medium — 79-91 cm
Large — 91-107 cm

X-Large — 104-122 cm
XX-Large — 109-132 cm

X-Small	580101
Small	580102
Medium	580103
Large	580104
X-Large	580105
XX-Large	580106



L1-L5

L1-L5 — Aspen Lumbar Support

Postoperative support / pain relief
 Deconditioned trunk musculature
 Post laminectomy syndrome
 Decompressive procedures
 Degenerative disc disease
 Bulging or herniated disc
 Insidious low back pain
 Chronic low back pain
 Fracture management
 Spondylolisthesis
 IDET procedure
 Spinal stenosis
 Sprains / strains
 Facet syndrome
 SI dysfunction
 Radiculopathy
 Spondylosis
 Severe DJD
 Sacrolitis
 Injections

Aspen Tri-Point FSO



Multiple Options from one brace
Anterior and posterior tension straps
Ergonomically designed molded pull tabs
Adjustable thoracic pad
Adjustable trochanter pad
Malleable aluminum struts
Optional chest strut

One Size Adjustable — Fits waists ranging from 61-114 cm and up to 139 cm with the extension panel

One Size Adjustable

591911

Dual Strut Configuration Kit

993905

Aspen Extension Panel

992505



Aspen Tri-Point FSO

Adult Degenerative Scoliosis
Adult Idiopathic Scoliosis
Spinal Cord Injury
Cerebral Palsy
Parkinson's Disease
Multiple Sclerosis
Brown-Séquard Syndrome
Hemiplegia/Post CVA (Stroke)





MP45 Plus

3-Point Hyperextension Brace

Product Code: TO1303N

Lightweight aluminium construction with unique spring-loaded 'floating' pelvic band. Thoracolumbar pad can be positioned horizontally or vertically. Unique 'Plus-lock' closure ensures controlled posterior loading.

Sternal component can be adjusted without bending, using serrated disc for angular adjustments in 5° settings.

Size	Code
Small	TO1303N/S
Medium	TO1303N/M
Large	TO1303N/L
Extra Large	TO1303N/XL

* Insert product code



Nhyco Cruciform

Hyperextension Orthosis

Product Code: TO410

A strong, lightweight aluminium alloy frame brace, weighing less than 700g. The replaceable EVA padding is non-toxic, waterproof and latex-free.

Features new 'Plus' safety lock closure system, easy for patient to apply/remove brace and ensures the same lumbar pressure each time.

Size	Code
Short	TO410CS
Regular	TO410CR
Long	TO410CL

* Insert product code



Chairback Brace

Product Code: CBB

A lightweight, spinal orthosis to maintain the spine in neutral and provide sacrolumbar stability. Made in heat-mouldable Kydex for a custom fit and lined with closed cell foam, this brace is lighter and thinner than conventional metal and leather braces, and ideal for elderly and osteoporotic patients who cannot tolerate heavier braces.

Size	Code	Hip Circumference	Brace Height (Centre back)
S	*/03	74 – 81 cm	29 cm
M	*/05	84 – 91 cm	32 cm
M-L	*/06	94 – 102 cm	33 cm
L	*/07	104 – 112 cm	34 cm

* Insert product code



Pro Tem Belt

Product Code: **TB**

Simple but effective, the high quality elastic and two flexible posterior stays make this our most popular belt.

The edges have quality Lycra binding and a back depth of 20cm cut away to 18cm at the front for increased comfort when sitting.

Size	Code	NHSCC Code	Waist Measurement
S	*/03	GRZ1252	66 – 74 cm
M	*/05	GRZ1253	76 – 84 cm
L	*/07	GRZ1254	86 – 94 cm
XL	*/08	GRZ1255	96 – 104 cm
XXL	*/09	GRZ1256	107 – 114 cm

* Insert product code



Proform Belt

Product Code: **PFM**

An elastic front with hook/loop closure and a cotton back section, with detachable comfort pad and two malleable steels which can be shaped to fit the spinal curve.

A webbing strap which can be positioned in the lumbar or sacral area, for extra compression on the spinal curve. Back depth is 28cm to support the thoracic area reducing to 20cm at the front for comfort.

Size	Code	Waist Measurement
S	*/83	58 – 66 cm
M	*/85	68 – 76 cm
L	*/87	79 – 86 cm
XL	*/88	89 – 97 cm
XXL	*/89	99 – 107 cm

* Insert product code

Contains Latex



The Proflex Belt

Features Include:

- Hook/loop fastening elastic belt, with an elastic back pocket to accommodate a panel.
- Available with a choice of back panels.
- A cotton fulcrum strap which can be positioned in the lumbar or sacral area from extra compression on the spine.
- The front depth is 20cm.



Proflex Belt (No panel)

Product Code: **PFN/8** (20cm back depth)

Size	Code	Waist Measurement
S	*3	64 – 76 cm
M	*5	79 – 89 cm
L	*7	89 – 99 cm
XL	*8	99 – 109 cm
XXL	*9	109 – 119 cm

* Insert product code

Contains Latex

Proflex Belt (Flexible Panel)

Product Code: **PF/8** (20cm back depth)

The back pocket encloses a panel of polythene lined with closed cell foam, which becomes flexible with body heat.

Size	Code	Waist Measurement
S	*3	64 – 76 cm
M	*5	79 – 89 cm
L	*7	89 – 99 cm
XL	*8	99 – 109 cm
XXL	*9	109 – 119 cm

* Insert product code

Contains Latex



The Lumbarform Belt

Features Include:

- Triple elastic fulcrum straps at the back which attach to one simple strap at the front, giving extra compression on the spinal curve.
- Two flexible stays which keep the belts smooth at the front for comfort when sitting.
- Ventilated elastic which has Lycra binding to prevent chafing.
- The Lumbarform belt range is available with a choice of back panels and stays.

Standard Lumbarform Belt (Flexible Panel)

Product Code: **LBF/9** (23cm back depth)

The double-layer back pocket encloses a panel of polythene lined with closed cell foam, which becomes flexible with body heat.

Size	Code	Waist Measurement
S	*3	69 – 76 cm
M	*5	79 – 86 cm
L	*7	89 – 97 cm
XL	*8	99 – 107 cm

* Insert product code

Contains Latex

Lumbarform Belt with Stays

Product Code: **LBS/9** (23cm back depth)

Two flexible posterior stays and a double-layer ventilated elastic back pocket stabilise the spine.

Size	Code	Waist Measurement
S	*3	69 – 76 cm
M	*5	79 – 86 cm
L	*7	89 – 97 cm
XL	*8	99 – 107 cm

* Insert product code

Contains Latex

Deep Lumbarform Belt (Flexible Panel)

Product Code: **LBF/12** (31cm back depth)

The double-layer back pocket encloses a panel of polythene lined with closed cell foam, which becomes flexible with body heat.

Size	Code	Waist Measurement
S	*3	69 – 76 cm
M	*5	79 – 86 cm
L	*7	89 – 97 cm
XL	*8	99 – 107 cm

* Insert product code

Contains Latex



Sacro Lumbar Support

Product Code: SLS

To accommodate the waist/hip differential, two layers of elastic are sewn at an angle giving this belt increased strength and good conformity. The back depth is 23cm, front depth is 18cm.

The posterior pocket is lined with Neoprene for extra comfort and warmth and a flexible panel (FPP/8) or a thermoplastic panel (TPP/8) can be supplied separately if required.

Size	Code	Waist Measurement
S	*/03	64 – 74 cm
M	*/05	76 – 84 cm
L	*/07	86 – 94 cm
XL	*/08	96 – 104 cm

* Insert product code

Contains Latex



TLSO Back Brace

The TLSO back brace is a unique prefabricated orthosis for post-operative, post traumatic, or chronic spinal pathology. The orthosis provides full circumferential rigid support from T-2 through S-5 and restricts forward, lateral and posterior motion and encourages a hyperextension posture.

The anterior and posterior panels are removable and heat formable to customize fit and functionality of the brace when worn. The anterior sternal assembly can be locked to prevent tampering and assure proper fit when donning and doffing the assembly. Elastic shoulder straps adjust easily and are wide to minimise any pressure. Compression adjustments are made by pulling each of the lateral lacer assemblies handles to the desired level of comfort for the patient. The components are washable. The TLSO is made of a breathable foam and mesh laminate that allows for optimal comfort to the patient."

Code	Size	Circumference Waist
60002	S	25" – 30"
60003	M	30" – 35"
60004	L	35" – 42"
60005	XL	42" – 48"
60006	XXL	48" – 57"

Universal Abdominal Binder

Product Code: ABW

Elastic and soft abdominal binder with universal sizing, ideal for compressing and supporting the abdomen

Size	Code	NHSCC Code
Universal	ABW/10	GP2717





Universal Two Piece Collar

Product Code: UTP/10

A universal two piece height adjustable collar that offers good support for the upper cervical spine.

- Preformed polyethylene shells with comfort liners
- Multiple height adjustment
- Large anterior aperture to allow airway management
- Lateral hook/loop closures with slim posterior panel for ease of fitting

Size
Universal



Budget Collar

Product Code: BC

Made in fully washable, medium density closed cell foam, which will not collapse with use. Stabilises the neck after trauma, or for long term use as a night resting collar.

Size	Code	NHSCC Codes	Depth x Length
S	*1	GQB693	7.5 x 43 cm
M	*2	GQB354	9 x 47 cm
L	*3	GQB694	9 x 51 cm
L/Deep	*4	GQB806	10 x 51 cm

* Insert product code



Philly Collar

Product Code: PCC

A lightweight, two-piece collar in Plastazote, with rigid occipital and mandibular posts. Strong comfortable support for the head and neck, incorporating a tracheal access hole. A Philly Collar Stabiliser (Product code: PCS/10), one size, is available on request; its anterior and posterior attachments further restrict cervical movement.

Code	NHSCC Code	Neck Circumference	Collar Depth
*33	GQB825	25.5 x 30.5 cm	8.5 cm
*53	GQB822	33 x 38 cm	8.5 cm
*73	GQB827	41 x 48 cm	8.5 cm
*34	GQB819	25.5 x 30.5 cm	11 cm
*35	GQB820	-	-
*54	GQB826	33 x 38 cm	11 cm
*55	GQB830	-	-
*74	GQB828	41 x 48 cm	11 cm
*75	GQB829	-	-

* Insert product code (5.5cm & 13.5cm depths available on request)



Halton Collar

Product Code: FC

Available in 8 sizes to fit most patients. Made in fully washable medium density closed cell foam, which will not collapse with use.

Size	Code	NHSCC Codes	Depth x Length
Large	FC/1	GQB1710	9 x 51 cm
Medium	FC/2	GQB1711	9 x 46 cm
Small	FC/3	GQB1712	9 x 42 cm
Long/Medium	FC/4	-	7.5 x 51 cm
Medium/Medium	FC/5	GQB1719	7.5 x 46 cm
Long/Shallow	FC/6	GQB1714	7 x 48.5 cm
Medium/Shallow	FC/7	GQB1715	7 x 44.5 cm
Extra Small	FC/8	-	7 x 39.5 cm



Upper Limb & Shoulder

Universal Lancaster Shoulder Sling

Product Code: LSS/10 (Single) LSS/25 (25 Pack)

- Padded, Non-stretch straps
- Fast release buckles
- 'Y'-Tab Adjustments
- Fold lines

Size	Code	NHSCC Codes	Dimensions
One Size	LSS/10	GRW004	13 x 41 cm
One Size	LSS/25	GRZ851	13 x 41 cm



Adjust sling length
using fold lines and circular
hook pads



Fast release buckles for easy
patient removal & re-application



Removable Hook/Loop
'Y' – Tabs



Scan QR Code
to watch video

www.youtube.com/promedicsorthopaedic



Easily adjust body and
shoulder strap length





Deluxe Shoulder Immobiliser

Product Code: DSI

A fully-adjustable body binder and humeral strap attachment holds the arm firmly to the body. The wrist sling is made in padded foam for comfort, and is attached to a strong webbing strap to position the wrist and stabilise the arm and shoulder.

Universal design fits left and right.

Size	Code		Dimensions
One Size	DSI/10	Body Binder	15 x 119 cm
		Wrist Cuff	15 x 18 cm
		Humeral Strap	8 x 51 cm



Arm Elevator Sling

Product Code: AES

Made in coutil cotton, this sling supports the forearm and upper arm close to the body, using the shoulder strap and body strap.

The strap extending from the wrist section includes a D-ring which can be attached to a drip stand to allow arm elevation.

Universal design fits right and left arm.

Size	Code	NHSCC Code	Dimensions
One Size	AES/10	GRZ2107	28 x 43 cm



Bradford Arm Elevator Sling

Product Code: AESB

Made in hook/loop attractive foam, this sling supports the forearm and upper arm close to the body, using the shoulder strap and body strap.

The strap extending from the wrist section includes a D-ring which can be attached to a drip stand to allow arm elevation.

Universal design fits left and right arm.

Size	Code	Dimensions
One Size	AESB/10	28 x 43 cm



Economy Shoulder Immobiliser

Product Code: ESI

An arm sling made in coutil cotton with an adjustable shoulder strap and body strap, to stabilise the arm close to the body.

Universal design fits right and left arm.

Size	Code	Dimensions
S	*/03	14 x 38 cm
M	*/05	17 x 43 cm
L	*/07	19 x 48 cm

* Insert product code



Sling & Swathe

Product Code: **SSW**

A simple padded sling supports the arm, and webbing shoulder straps cross at the back to hold the arm in position. An adjustable body binder holds the upper arm close to the body for increased shoulder stability.

Universal design fits right and left arm.

Size	Code	Body Binder	Padded Sling
S	*/03	13 x 41 cm	20 x 28 cm
M	*/05	13 x 94 cm	23 x 31 cm
L	*/07	15 x 119 cm	25 x 33 cm

* Insert product code



Universal Pro Sling

Product Code: **UPS**

A universal collar and cuff type sling that is fully adjustable to fit all sizes from child to adult. Hook/loop fastening with adjustable loop buckle to secure sling around the wrist.

Universal design fits right and left arm.

Size	Code	Dimensions
Standard	UPS/10	5 x 112 cm

(Supplied in packs of 10)



Boscombe Sling

Product Code: **BOS**

A development of the Pro Sling design, incorporating an adjustable double sling which supports shoulder, elbow, forearm and wrist.

Size	Code	Dimensions
S	*/03	5 x 102 cm
M	*/05	5 x 117 cm
L	*/07	5 x 132 cm

* Insert product code




Trauma Sling

Product Code: TAU/10

The Trauma Sling is a high quality universal sling that has multiple uses within the healthcare field. The Trauma Sling can be used as a traditional sling or for arm elevation following trauma, burns or post-operative management of the shoulder or upper arm.

The design and quality of materials used offer protection, comfort and ease of access to the damaged limb.

The Trauma Sling is easy to apply and is supplied vacuum packed to reduce storage volume.



Adjustment at elbow allows for varying degree of elbow flexion/extension.

Size	NHSCC Code
Universal	EKM000



Paediatric Trauma Sling

Product Code: TAU/P

Same excellent features as Trauma Sling, but with closed seam behind elbow.

The design and quality of materials used offer protection, comfort and ease of access to the damaged limb. The Trauma Sling is easy to apply and is supplied vacuum packed to reduce storage volume.



Size	NHSCC Code
Universal	EKM008

- Comfortable
- Protective
- Universal



Jura Abduction Sling PLUS

Product Code: JASSP/10

The Jura Plus has a three point off-loading harness that offers excellent support and comfort. The contoured abduction pad is profiled for improved patient compliance and holds the arm in 45° of abduction. A body swathe combined with upper/lower arm wraps offer good stability and shoulder fixation. Hand resting splint supplied.

- Rotator cuff repairs
- Immobilisation following joint replacement
- Glenohumeral dislocations / subluxation
- Posterior/ Anterior capsule repairs
- Total shoulder reconstructions
- Soft tissue repairs / strains

Size
Universal



Jura Shoulder Immobiliser

Product Code: JSI

The Jura Shoulder Immobiliser can be used as an arm sling or as a shoulder immobiliser and it is designed to be used on both left and right arms.

Strong cotton sling with fully adjustable padded shoulder and body straps.

Size	NHSCC Codes	Forearm Length
XS	GRZ549	< 33 cm
S	GRZ550	33 – 38 cm
M	GRZ551	37 – 43 cm
L	GRZ552	42 – 48 cm
XL	-	48 cm >

* Insert product code



Jura Shoulder Lock

Product Code: JSL/10

A universal fast-fit low profile shoulder immobiliser. The Jura Shoulder Lock comes complete with a body wrap to which the upper arm wrap is secured to effectively immobilise the shoulder. The wrist wrap fixes to the anterior of the body wrap and this in turn is secured to the shoulder strap using a fast-fit snap buckle.

The shoulder strap can be adjusted posteriorly to achieve a comfortable fit for the patient and also to alter the brace to fit left or right. There is a large comfort pad to distribute pressure over the shoulder /clavicle.

This is a shoulder immobiliser that is easy to fit and the modular design greatly improves patient compliance.

Size
Universal



Universal Combo Sling

Product Code: UCS/10

A comfortable off-loading sling designed for patients recovering from shoulder surgery. Two wedges are included which can be applied in a choice of 3 different combinations: 45° Neutral, 90° Neutral or 15° Abduction. Lightweight and breathable 3D spacer fabric sling.

- Exercise ball
- Detachable waist strap
- Opening to relieve elbow pressure



Size
Universal

Fold to fit length



Evans Abduction Sling

Product Code: EAS

This is a low cost universal option that utilizes the popular Lancaster Shoulder Immobiliser. In addition to the comfortable sling with shoulder and body straps – the Evans Abduction Sling incorporates a 15° abduction pad that is easily secured to the inside of the sling.

The inner profile of the abduction pad is shaped to rest comfortably against the patient. The Evans Abduction sling is useful in the treatment of various shoulder conditions

- Rotator cuff repairs
- Glenohumeral dislocations / subluxation
- Posterior/ Anterior capsule repairs
- Soft tissue repairs / strains

Size
Universal



Jura Abduction Sling

Product Code: JASS

The Jura Abduction has a three-point off-loading harness that offers excellent support and comfort. The sling is made from a lightweight breathable 3D spacer fabric and attaches to the abduction pads. The sling is supplied with a contoured 15° pad that can be used on its own.

In addition to this there is an additional pad that increases abduction to 45° and positions the shoulder in a neutral position. A detachable exercise ball is also provided. The shoulder harness is fitted with easy release buckles to facilitate easy application and removal of the sling.

Size	Code	NHSCC Codes	Forearm Length	
S	*/03	-	38 cm	<ul style="list-style-type: none"> • Rotator cuff repairs • Glenohumeral dislocations/ subluxation • Posterior/ Anterior capsule repairs • Total shoulder reconstructions • Soft tissue repairs / strains
M	*/05	GRZ2125	40 cm	
L	*/07	GRZ2126	43 cm	
XL	*/08	-	46 cm	

* Insert product code



Jura Clavicle Splint

Product Code: JCS/10

Soft padded shoulder straps with return loops and padded back pad. Comfort pads can be positioned anywhere along the length of the straps – for example in the axilla area.

Straps can be easily trimmed to length and refastened using the hook Y-tab closures.

Size	NHSCC Codes
Universal	GP2719




Subluxation Cuff

Product Code: SCF

A humeral cuff with soft cotton lining supports the injured arm, while the padded webbing strap fits under the opposite axilla. Universal design fits right and left arm.

Size	Code	Mid Upper Arm Circumference
S	*/03	24 – 29 cm
M	*/05	29 – 34 cm
L	*/07	34 – 40 cm

* Insert product code



Clavicle Splint


Product Code: CLS

Two padded webbing straps adjust in a figure of eight design fixed to a pad over the scapulas. Designed to assist shoulder retraction.

The length of the shoulder strap is measured from top of back pad, under axilla, to bottom of back pad.

Size	Code	NHSCC Codes	Length of Shoulder Straps
XS	*/02	-	48 cm
S	*/03	-	56 cm
M	*/05	-	68 cm
L	*/07	GP1440	87 cm
XL	*/08	-	92 cm

* Insert product code



Neoprene Shoulder Brace

Product Code: NSB

Two adjustable upper arm cuffs fix with hook/loop at the back of the neck, to prevent shoulder subluxation and hold the head of humerus in the glenoid cavity.

Size	Code	NHSCC Code	Mid Upper Arm Circumference
S	*/03	-	23 – 28 cm
M	*/05	GR22129	28 – 33 cm
L	*/07	-	33 – 38 cm
XL	*/08	-	38 – 43 cm

* Insert product code



Scan QR Code
to watch video



Ascender Telescopic ROM Elbow Brace


Product Code: ASCE

Designed and manufactured in the UK by Promedics. The new Ascender telescopic ROM elbow brace is easily adjustable and secure. With additional features offering superior comfort and support to patients recovering from surgery, injuries or instabilities.

- Maximum strength, minimum weight
- High grade anodised aluminium
- 5 drop lock settings
- Quick-clip Range of Motion Settings

The range of motion settings can be easily adjusted via the spring loaded quick-clips by a simple 'pull and slide' motion. The design prevents unintentional changes to the flexion and extension settings. The drop lock is ergonomically designed to be easily locked and unlocked due to its curved concave profile.

Size	Product Code
Left	ASCE/L/REG
Right	ASCE/R/REG



Cubital Tunnel Splint

Product Code: CTS/10

A universal easy to fit splint that effectively restricts elbow flexion to help reduce the pain associated with Cubital Tunnel Syndrome. A soft padded wrap with an adjustable aluminium stay. The splint is easy to apply, lightweight and breathable.

Size
Universal





Hinged Elbow Wrap

Product Code: HEW

A front-opening Neoprene sleeve with hinges, which limit medial/lateral movement and elbow flexion. Hook/loop straps above and below the elbow ensure a good fit. Specifically designed for the patient with limited dexterity.

Size	Code	NHSCC Codes	Circumference 13cm Below Elbow Crease	Circumference 13cm Above Elbow Crease
S	*/03	GQT1154	13 – 18 cm	23 – 28 cm
M	*/05	GQT1148	18 – 23 cm	28 – 33 cm
L	*/07	-	23 – 28 cm	33 – 38 cm
XL	*/08	-	28 – 33 cm	38 – 43 cm

* Insert product code



Orthoprene Elbow Sleeve

Product Code: OES

A simple pull-on Neoprene sleeve to give light support and compression around the elbow area and ease the pain of inflammation.

Size	Code	NHSCC Code	Elbow Circumference
S	*/03	-	18 – 20 cm
M	*/05	GRZ1951	23 – 25 cm
L	*/07	GRZ1952	28 – 30 cm

* Insert product code



Tennis Elbow Strap

Product Code: TEP

A D-ring on strong cotton webbing adjusts tension to ease the pain of medial or lateral epicondylitis.

Size	Code	NHSCC Code	Strap Length
One Size	TEP/10	GRZ5300	5 x 46 cm

(Supplied in packs of 6)



Jura Clasp

Product Code: JUC

Designed to reduce the pain and discomfort associated with elbow tendinopathies, conditions more commonly referred to as epicondylitis and tennis or golfer's elbow.

Contains Latex

Size	Code	NHSCC Codes	Arm Circumference
XS	*/02	GQT117	22 – 23 cm
S	*/03	GQT118	24 – 25 cm
M	*/05	GQT119	26 – 27 cm
L	*/07	GQT120	28 – 30 cm
XL	*/08	GQT121	31 – 34 cm

* Insert product code



Vari Pad Tennis Elbow Support

Product Code: VPS

Made in Neoprene with a D-ring hook/loop strap and adjustable pressure pad, which is easily positioned on the insertion point to ease the pain of medial or lateral epicondylitis.

This support has a high level of patient compliance, as the stretch Neoprene accommodates the changing muscle bulk of the moving forearm.

Size	Code	NHSCC Code	Circumference 2.5cm Below Elbow Crease
XS	*/02	GRZ1927	18 – 20 cm
S	*/03	GQT1105	23 – 25 cm
M	*/05	GQT1106	28 – 30 cm
L	*/07	GQT1107	33 – 35 cm
XL	*/08	GQT1110	37 – 40 cm

* Insert product code



Replay Tennis Elbow Support

Product Code: RVP

A tennis elbow support in 'Ventoprene', a ventilated Neoprene providing compression and heat retention. The adjustable pressure pad can be positioned to ease the pain of medial or lateral epicondylitis.

Stretch 'Ventoprene' accommodates the changing muscle bulk of the moving forearm when playing racquet sports.

Size	Code	Circumference 2.5cm Below Elbow Crease
S	*/03	23 – 25 cm
M	*/05	28 – 30 cm
L	*/07	33 – 35 cm
XL	*/08	37 – 40 cm

* Insert product code



Sling strap

Sling Strap is a non-stretch padded 50mm wide material that is combined with Hook Y-Tabs to provide an adjustable and versatile sling system.

Sling Strap can be cut to the required length to produce a simple collar n' cuff or similar type of sling. This reduces waste and the adjustability of the material ensures that the sling will fit and can be easily adjusted. Sling Strap is padded for comfort around the neck and the modern black colour improves patient compliance and looks better for longer. Sling Strap is supplied in a pack that contains: 2 x 12m or 2 x 6m rolls and 20 Y-tab closures.

Product Code	Length	NHSCC Code	Width
SS	2 x 12 Metres	GRZ554	5 cm
SS/5.5	2 x 6 Metres	ESM137	5.5cm



Lower Limb

Ascender

Telescopic ROM knee Brace

Product Code: ASC/10

Designed and manufactured in the UK by Promedics. The new Ascender telescopic ROM knee brace is universal, easily adjustable and secure. With additional features offering superior comfort and support to patients recovering from knee surgery, injuries or instabilities.

- Universal sizing
- Maximum strength, minimum weight
- High grade anodised aluminium
- Length range 44 cm – 70 cm
- 5 drop lock settings
- Quick-clip Range of Motion Settings

The range of motion settings can be easily adjusted via the spring loaded quick-clips by a simple 'pull and slide' motion.

The design prevents unintentional changes to the flexion and extension settings.

The drop lock is ergonomically designed to be easily locked and unlocked due to its curved concave profile.

Size	Product Code
Universal	ASC/10

Ascender Telescopic ROM Knee Brace with Foam Pads

Product Code: ASCFP/10

Size	Product Code
Universal	ASCFP/10



Scan QR Code
to watch video

www.ascender.pro

ascender®

Range of Motion Control



5 DROP LOCK
SETTINGS

Universal Sizing
Length range 44cm – 70cm



Designed for Comfort
Optimised for Airflow
Silicone for Ultra Grip





Genurange Universal

Product Code: GN/0510

To reduce stock-holding of different circumferences, the one size Universal wraps can be cut, without fraying, to fit most sizes. Open-cell padding tibial and femoral liners are attached by four or six straps for an excellent fit.

Hinges, applied to the liners, align with the tibial condyles and the liners can be removed for replacement or laundering.

Available in four lengths. Replacement wraps and straps available, see below.

Size	Code	NHSCC Codes	Length
S	*/18	GRS368	45 cm
M	*/20	GRS369	50 cm
L	*/22	GRS370	55 cm
XL	*/24	GRS371	60 cm

* Insert product code



Genuranger

Product Code: GN/04

Towelling-lined tibial and femoral liners are attached by four straps; six straps on 24", for an excellent fit. The hinges are applied to the liners to align with the tibial condyles, and the liners can be removed for replacement or laundering.

Available in five circumferences and five lengths. Replacement wraps and straps available.

Replacement Liners and Straps

To order replacement liners and straps for different sizes of Genurange Universal and Genuranger, add ' W ' after ' GN ' in the code, e.g. **GNW/0510/18** or **GNW/0405/18**

Size	40 cm	45 cm	50 cm	55cm	60 cm	Knee Circumference
S	*03/16	-	*03/20	*03/22	*03/24	35 – 40 cm
M	*05/16	*05/18	*05/20	*05/22	*05/24	40 – 45 cm
L	*07/16	*07/18	*07/20	*07/22	*07/24	45 – 50 cm
XL	*08/16	*08/18	*08/20	*08/22	*08/24	55 – 60 cm
XXL	*09/16	*09/18	*09/20	*09/22	*09/24	60 – 65 cm

* Insert product code

Genurehab

Product Code: GN/0310

The 40cm Genurehab is a short version of the Genurange Universal, designed for long term support and motion control during the more active phase of recovery. The one-size Genurehab wraps can be cut without fraying, to fit most sizes.

To order replacement liners and straps, add 'W' after 'GN' in the code, e.g. GNW/0310

Size
Universal



Genurange Lite

Product Code: GN/08

A 'lite' wrap-around knee brace with ROM hinges. 3D spacer fabric with open mesh posterior panel makes this a supportive, comfortable well fitting brace. Suitable for the treatment various knee conditions including mild/moderate collateral ligament injuries and cruciate mild instability of the cruciate ligaments.

Size	Code	Circumference 15cm above mid-patella
XS	*02	32 – 38 cm
S	*03	38 – 44 cm
M	*05	44 – 51 cm
L	*07	51 – 58 cm
XL	*08	58 – 66 cm
XXL	*09	66 – 75 cm

* Insert product code



GENURANGE

- A geared, polycentric hinge which mimics the action of the knee.
- Four small pins which are inserted into apertures to set the hinge at the required angles of flexion and extension, between 10° and 110°.
- A plastic hinge cover, with screw closure, to prevent patients tampering with the hinge mechanism.
- Features plastic tags to prevent loss of pins.




Jura Knee Immobiliser

Product Code: BTPK

A universal three panel knee immobiliser manufactured from a quality foam laminate with terry towel lining. This brace is designed to offer a lower cost alternative to hospital departments with high volume usage.

Code	NHSCC Code	Splint Length
BTPK/12	-	30 cm
BTPK/14	GRA392	36 cm
BTPK/16	GRA393	41 cm
BTPK/18	GRA394	46 cm
BTPK/20	GRA395	51 cm
BTPK/22	GRA396	56 cm
BTPK/24	GRA397	60 cm
BTPK/26	-	66 cm
BTPK/28	-	71 cm

(Sold in packs of 5 via NHSCC)



Post-Op Knee Splint

Product Code: PKS/20

A three panel splint made in foam laminate with two detachable side panels which can be positioned to fit most leg circumferences. The four steels are removable for the patient to launder the splint.

The length of the splint is 51cm.

Size
Universal



Three Panel Knee Immobiliser

Product Code: TPK

A three panel splint made in high quality foam laminate with bound edges, designed for longer term use. The two detachable side panels can be positioned to fit most leg circumferences.

The four steels are removable for the patient to launder the splint.

Code	NHSCC Code	Splint Length	Code	NHSCC Code	Splint Length
*/12	GRS1180	30 cm	*/22	GRS1178	56 cm
*/14	GRS1181	36 cm	*/24	GRS1179	61 cm
*/16	GRS1175	41 cm	*/26	GRS1192	66 cm
*/18	GRS1176	46 cm	*/28	GRS1193	71 cm
*/20	GRS1177	51 cm			

* Insert product code



Jura One Piece Immobiliser

Product Code: OPK

A well fitting comfortable knee immobiliser designed to maintain knee extension following surgery, trauma and soft tissue injuries.

- One piece foam laminate wrap with towel lining
- Detachable hook / loop straps for easy positioning

Size	Code	Length	Circumference (Mid Thigh)
Paediatric	*/12	30 cm	Up to 32 cm
S	*/16	41 cm	Up to 41 cm
S	*/20SM	51 cm	Up to 52 cm
M-L	*/20ML	51 cm	47 – 66 cm
S	*/24SM	56 cm	Up to 52 cm
M-L	*/24ML	56 cm	47 – 66 cm
S	*/28SM	71 cm	Up to 52 cm
M-L	*/28ML	71 cm	47 – 66 cm

* Insert product code

Canvas Knee Immobiliser

Product Code: CNI

A one-piece knee immobiliser made in soft padded canvas with four removable steels which can be removed for the patient to launder the splint.

The length of the splint is 48cm.

Size	Code	Knee Circumference
S	*/03	28 – 34 cm
M	*/05	34 – 40 cm
L	*/07	40 – 46 cm
XL	*/08	46 – 56 cm

* Insert product code

Neosplint

Product Code: NS

A one-piece knee immobiliser in Neoprene, which conforms very well to the leg. Designed for long term wear, the four steels can be removed to launder the splint.

The length of the splint is 48cm.

Size	Code	Knee Circumference
S	*/03	28 – 34 cm
M	*/05	34 – 40 cm
L	*/07	40 – 46 cm
XL	*/08	46 – 56 cm

* Insert product code



Air X Polycentric Knee Brace

Product Code: AMP

A 30cm long wrap-around knee brace offering good medial and lateral support using quality polycentric hyperstop hinges. Hook/loop return closures make this brace easy to apply and the posterior straps can be adjusted quickly to allow accurate midline positioning of the polycentric joints.

Air X gives great compression and support and the properties of the textile offer a high degree of breathability and reduction in moisture. This makes the brace comfortable to wear for extended periods of time and improves patient compliance.

Size	Code	Knee Circumference
S	*/03	28 – 34 cm
M	*/05	34 – 40 cm
L	*/07	40 – 46 cm
XL	*/08	46 – 56 cm
XXL	*/09	56 – 61 cm

* Insert product code



Air X SK Knee Brace

Product Code: ASK

The SK Knee Brace has been the industry leading wrap-around knee brace for many years and thousands of patients have benefited from this easy to use comfortable brace in the treatment of mediolateral instability.

The SK Knee Brace is now available using the excellent properties of the Air X 3D spacer fabric. Air X gives great compression and support and the properties of the textile offer a high degree of breathability and reduction in moisture. This makes the brace even more comfortable to wear for extended periods of time and improves patient compliance.

Size	Code	Knee Circumference
XS	*/02	25 – 28 cm
S	*/03	28 – 34 cm
M	*/05	34 – 40 cm
L	*/07	40 – 46 cm
XL	*/08	46 – 56 cm
XXL	*/09	56 – 61 cm

* Insert product code



Air X Long ROM Knee Brace

Product Code: ARK

A 50cm long wrap-around knee brace with Range of Motion (ROM) Joints. Air X 3D spacer fabric offers excellent comfort, compression support and fit. The brace has a range of rigid and elastic straps that can be adjusted anteriorly and posteriorly to allow a perfect fit and joint alignment. The ROM joints have a wide range of flexion and extension settings and the brace is suitable for the treatment of various knee conditions including mild/moderate collateral ligament injuries and mild instability of cruciate ligaments.

Size	Code	Knee Circumference
XS	*/02	25 – 28 cm
S	*/03	28 – 34 cm
M	*/05	34 – 40 cm
L	*/07	40 – 46 cm
XL	*/08	46 – 56 cm
XXL	*/09	56 – 61 cm

* Insert product code



Air X VK Long Knee Brace

Product Code: AVK

The VK Knee Brace has been one of the industry leading long wrap-around knee braces for many years and thousands of patients have benefited from this easy to use comfortable brace in the treatment of mediolateral instability. The longer length (41cm) provides increased medial and lateral support for the unstable valgus or varus knee. The front opening design makes it easy to apply – even for those patients with limited dexterity

The VK Knee Brace is now available using the excellent properties of the Air X 3D spacer fabric.

Air X gives great compression and support and the properties of the textile offer a high degree of breathability and reduction in moisture. This makes the brace even more comfortable to wear for extended periods of time and improves patient compliance.

Size	Code (Left)	Code (Right)	Knee Circumference
XS	*/12	*/22	25 – 28 cm
S	*/13	*/23	28 – 34 cm
M	*/15	*/25	34 – 40 cm
L	*/17	*/27	40 – 46 cm
XL	*/18	*/28	46 – 56 cm
XXL	*/19	*/29	56 – 61 cm

* Insert product code



Jura 3 Panel Immobiliser

Product Code: JTI

- Universal three panel knee immobiliser
- Return hook/loop fastenings
- Removable aluminium stays
- Bound edges

Code	NHSCC Code	Splint Length
JTI/12	GRA5236	30 cm
JTI/14	GRA5237	36 cm
JTI/16	GRA5238	41 cm
JTI/18	GRA5239	46 cm
JTI/20	GRA5240	51 cm
JTI/22	GRA5241	56 cm
JTI/24	GRA5242	60 cm
JTI/26	GRA5243	66 cm
JTI/28	GRA5244	71 cm

(Sold in packs of 5 via NHSCC)



Ventilated SK Knee Plus

Product Code: VSKP

The SK Knee Brace has been the industry leading wrap-around knee brace for many years and thousands of patients have benefited from this easy to use comfortable brace in the treatment of mediolateral instability.

The Ventilated SK Plus takes the brace to the next level – using multi-perforated towel-lined latex free neoprene – the SK Plus provides excellent medial and lateral support and great patient comfort. The ventilated material reduces moisture build up whilst the towel lining improves comfort.

The brace and hinge casings are fully bound to enhance the overall look of the brace whilst the hook/loop straps have grip-taps to assist patients with limited dexterity.

Size	Code	Knee Circumference
XS	*/02	25 – 28 cm
S	*/03	28 – 34 cm
M	*/05	34 – 40 cm
L	*/07	40 – 46 cm
XL	*/08	46 – 56 cm
XXL	*/09	56 – 61 cm

* Insert product code



SK Knee Brace

Product Code: SK

The front opening design and easy-open tabs on the hook/loop straps make this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity. Hyperstop hinges provide medial and lateral support for the unstable knee, while allowing flexion.

Also available in black

Size	Code (Beige)	NHSCC Code	Code (Black)	Knee Circumference
XS	*/02	GRA413	*/02/BK	25 – 28 cm
S	*/03	GRA414	*/03/BK	28 – 34 cm
M	*/05	GRA415	*/05/BK	34 – 40 cm
L	*/07	GRA416	*/07/BK	40 – 46 cm
XL	*/08	GRA417	*/08/BK	46 – 56 cm
XXL	*/09	GRA418	*/09/BK	56 – 61 cm

* Insert product code



Nuffield SK Knee Brace

Product Code: NSK

The front opening design and easy-open tabs on the hook/loop straps make this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity.

The brace incorporates a polycentric hinge which follows the true motion of the knee. The top opening design allows easy removal of the hinges for laundering.

Size	Code	Knee Circumference
XS	*/02	25 – 28 cm
S	*/03	28 – 34 cm
M	*/05	34 – 40 cm
L	*/07	40 – 46 cm
XL	*/08	46 – 56 cm

* Insert product code



VK Knee Brace

Product Code: VK

The longer 41cm metal hyperstop hinges provide increased stability and medial and lateral support for the unstable valgus or varus knee. The front opening design and easy-open tabs on the hook/loopstraps make this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity.

The VK Knee brace is sided left and right for improved fit.

Size	Code (Left)	Code (Right)	Knee Circumference
XS	*/12	*/22	25 – 28 cm
S	*/13	*/23	28 – 34 cm
M	*/15	*/25	34 – 40 cm
L	*/17	*/27	40 – 46 cm
XL	*/18	*/28	46 – 56 cm
XXL	*/19	*/29	56 – 61 cm

* Insert product code




Lycrafleece Patella Lux

Product Code: BLPL

Provides medial or lateral stabilisation of the knee. The crescent shaped buttress pad can be placed on the medial or lateral side of the patella. A non-slip rubber lining is used to ensure the support does not move after fitting. The Patella Lux is ideal for lateral or medial subluxation of the patella. Universal design fits left and right.

Size	Code	NHSCC Codes	Knee Circumference
XS	*/02	-	25 – 28 cm
S	*/03	GRA344	28 – 36 cm
M	*/05	-	33 – 41 cm
L	*/07	GRA346	38 – 46 cm
XL	*/08	GRA347	45 – 55 cm

* Insert product code



Economy Knee Sleeve

Product Code: EN

A Neoprene sleeve with open patella providing warmth and light compressional support.

Size	Code	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/02	30 cm	30 – 33 cm	38 cm
S	*/03	33 cm	33 – 36 cm	46 cm
M	*/05	38 cm	36 – 38 cm	50 cm
L	*/07	41 cm	38 – 43 cm	53 cm
XL	*/08	43 cm	43 – 48 cm	58 cm
XXL	*/09	52 cm	48 – 56 cm	61 cm

* Insert product code



Standard Knee Sleeve – Open Patella

Product Code: SNO

A Neoprene sleeve with a reinforced open patella section providing warmth and light compressional support.

Size	Code	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/02	30 cm	30 – 33 cm	38 cm
S	*/03	33 cm	33 – 36 cm	46 cm
M	*/05	38 cm	36 – 38 cm	50 cm
L	*/07	41 cm	38 – 43 cm	53 cm
XL	*/08	43 cm	43 – 48 cm	58 cm
XXL	*/09	52 cm	48 – 56 cm	61 cm

* Insert product code



Spiral Knee Sleeve

Product Code: SP

A Neoprene sleeve with two flexible spirals, providing light compression and medial/lateral support.

Size	Code	NHSCC Code	Circumference 15cm Below Knee	Knee Circumfer- ence	Circumference 15cm Above Knee
XS	*/02	-	30 cm	30 – 33 cm	38 cm
S	*/03	-	33 cm	33 – 36 cm	46 cm
M	*/05	GRA386	38 cm	36 – 38 cm	50 cm
L	*/07	GRA387	41 cm	38 – 43 cm	53 cm
XL	*/08	GRA389	43 cm	43 – 48 cm	58 cm
XXL	*/09	-	52 cm	48 – 56 cm	61 cm

* Insert product code



Butress Knee Sleeve

Product Code: BN

A standard Neoprene knee sleeve with open patella, incorporating a buttress pad. The Butress Knee Sleeve can be selected with a buttress pad to oppose medial/lateral patella migration; caused by ligament deficiency.

The buttress pad positions are defined as viewed by the patient, for example a left buttress pad would give lateral support to the left leg and medial support to the right leg. A full buttress is also available.

Size	Code (Left)	NHSCC Codes	Code (Right)	NHSCC Codes	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/12	-	*/22	-	30 cm	30 – 33 cm	38 cm
S	*/13	-	*/23	-	33 cm	33 – 36 cm	46 cm
M	*/15	-	*/25	-	38 cm	36 – 38 cm	50 cm
L	*/17	GRA351	*/27	GRA357	41 cm	38 – 43 cm	53 cm
XL	*/18	GRA352	*/28	GRA358	43 cm	43 – 48 cm	58 cm
XXL	*/19	GRA353	*/29	GRA359	52 cm	48 – 56 cm	61 cm

* Insert product code



Stabilised Knee Sleeve

Product Code: ST

A pull-on Neoprene sleeve incorporating a buttress pad, spiral stays for medial and lateral support, and hook/loop straps above and below the knee for improved fit. The Stabilised Knee Sleeve can be selected with a buttress pad to oppose medial/lateral patella migration; caused by ligament deficiency.

The buttress pad positions are defined as viewed by the patient, for example a left buttress pad would give lateral support to the left leg and medial support to the right leg. A full buttress is also available.

Size	Code (Left)	Code (Right)	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/12	*/22	30 cm	30 – 33 cm	38 cm
S	*/13	*/23	33 cm	33 – 36 cm	46 cm
M	*/15	*/25	38 cm	36 – 38 cm	50 cm
L	*/17	*/27	41 cm	38 – 43 cm	53 cm
XL	*/18	*/28	43 cm	43 – 48 cm	58 cm
XXL	*/19	*/29	52 cm	48 – 56 cm	61 cm

* Insert product code



Hinged Knee Sleeve

Product Code: HN

A pull-on Neoprene sleeve, with strong metal hyperstop hinges to control medial and lateral instability.

Hook/loop straps above and below knee give improved fit.

Size	Code	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee
XS	*/02	30 cm	30 – 33 cm	38 cm
S	*/03	33 cm	33 – 36 cm	46 cm
M	*/05	38 cm	36 – 38 cm	50 cm
L	*/07	41 cm	38 – 43 cm	53 cm
XL	*/08	43 cm	43 – 48 cm	58 cm
XXL	*/09	52 cm	48 – 56 cm	61 cm

* Insert product code



Ankle & Footwear

Promedics Airstep™ Walker

Product Code: ASW (Tall) ASSW (Short)

Airstep™ is an innovative walker with a fully integrated inflation system that allows independent adjustment.

- Super lightweight – ventilated shell system
- Breathable fitted contoured liner
- Low profile sole unit – minimal height discrepancy
- Easy three strap closure
- Fast and easy to fit
- No separate pump required
- Excellent stability
- Achilles wedges available

Size	Code (Tall)	NHSCC Codes	Shoe Size
XS	ASW/02	GRA5126	1.5 – 3
S	ASW/03	GRA2847	3 – 5.5
M	ASW/05	GRA2848	6 – 9
L	ASW/07	GRA2849	9.5 – 11
XL	ASW/08	GRA2850	11 +
Size	Code (Short)	NHSCC Codes	Shoe Size
XS	ASSW/02	GRA5263	1.5 – 3
S	ASSW/03	GRA5253	3 – 5.5
M	ASSW/05	GRA5254	6 – 9
L	ASSW/07	GRA5255	9.5 – 11
XL	ASSW/08	GRA5264	11 +



Airstep Achilles Wedges

Product Code: ASAW

Sold in packs of 5.



Left	NHSCC Codes	Right	NHSCC Codes
ASAW/RL	ERD087	ASAW/RR	ERD091



Dual integrated pumps
for fast, easy inflation




Short version



Scan QR Code to watch the Airstep Walker Instructional video on our Youtube channel.



Low profile, hard wearing sole unit
Improved gait function



Adjustable / removable toe guard



**DUAL
AIRBAGS**
- INDEPENDENTLY
ADJUSTABLE



Jura Air Shield Walker

Product Code: JWA

Available with charcoal liners, this lightweight low profile walker allows good gait pattern and has independently adjustable dual airbags for improved stability and support.

Features a comfortable breathable liner with toe cover. Anterior & posterior shells offer improved stability and immobilisation with velour strap fastening. The shells can be removed for staged rehabilitation.

Replacement liners available in charcoal.

Size	Code	NHSCC Codes	Replacement Liners Code	Shoe Size
S	*03/C	GRA523	*03/L/C	3 – 5.5
M	*05/C	GRA525	*05/L/C	6 – 9
L	*07/C	GRA527	*07/L/C	9.5 – 11
XL	*08/C	GRA529	*08/L/C	11 +

* Insert product code

Jura Achilles Wedges

Product Code: JWW

Extension straps included.



Left	NHSCC Codes	Right	NHSCC Codes
JWW/RL	GRA219	JWW/RR	GRA220

Walker Boot Socks

- Designed for extra comfort whilst wearing walker boot
- Helps prevent chafing and absorb sweat
- Available in both Tall and Short to fit most walker boots
- Sold in pairs



Size	Code
Short	ASWSS
Tall	ASWS



Jura Air Shield Walker – Short

Product Code: JWA

Only available with charcoal liners.

Size	Code	NHSCC Codes	Replacement Liners Code	Shoe Size
S	*03/C	GRA535	*03/L/C	3 – 5.5
M	*05/C	GRA536	*05/L/C	6 – 9
L	*07/C	GRA537	*07/L/C	9.5 – 11
XL	*08/C	-	*08/L/C	11 +

* Insert product code

Jura Air Shield Walker

Extension straps

Supplied as standard with Jura Achilles Wedges.

Code
JWA-ES





Jura Fixed Walker

Product Code: JWF

One piece injection moulded walker with 'paddle' uprights to reduce pressure on lower limb and provide improved stability and immobilisation. Foot section can accommodate full dressings and padding.

New lightweight design sole unit to greatly improve shock absorption, gait pattern to reduce the overall weight of the walker to an absolute minimum.

Reduced strap fastening to simplify removal / application and improve patient acceptance. Comfortable detachable liners with replacements available.

Available with standard or Air Liner.

Fixed Walker

Size	Code	NHSCC Codes	Replacement Liners Code	Shoe Size
XS	*/02	GRA543	*/02/L	1.5 – 3
S	*/03	GRD814	*/03/L	3 – 5.5
M	*/05	GRD815	*/05/L	6 – 9
L	*/07	GRD816	*/07/L	9.5 – 11
XL	*/08	GRA550	*/08/L	11 +

* Insert product code

Fixed Walker with Air Liner

Size	Code	NHSCC Codes	Replacement Liners Code	Shoe Size
XS	*/02/A	GRA544	*/02/L/A	1.5 – 3
S	*/03/A	GRA545	*/03/L/A	3 – 5.5
M	*/05/A	GRA546	*/05/L/A	6 – 9
L	*/07/A	-	*/07/L/A	9.5 – 11
XL	*/08/A	-	*/08/L/A	11 +

* Insert product code



Jura Fixed Walker – Short

Product Code: JWFS

Short version of the Jura Fixed Walker.

Fixed Walker – Short

Size	Code	NHSCC Codes	Replacement Liners Code	Shoe Size
XS	*/02	-	*/02/L	3 – 5.5
S	*/03	GRA553	*/03/L	3 – 5.5
M	*/05	GRA555	*/05/L	6 – 9
L	*/07	GRA557	*/07/L	9.5 – 11
XL	*/08	-	*/08/L	11 +

* Insert product code

Fixed Walker Short with Air Liner



Size	Code	Replacement Liners Code	Shoe Size
XS	*/02/A	*/02/L/A	3 – 5.5
S	*/03/A	*/03/L/A	3 – 5.5
M	*/05/A	*/05/L/A	6 – 9
L	*/07/A	*/07/L/A	9.5 – 11
XL	*/08/A	*/08/L/A	11 +

* Insert product code



ROM Walker

Product Code: JWR

Fully moulded plastic foot section with rubber sole and padded insole. Easily adjustable 'Push-Dial-Set' range of motion settings. Chosen settings can be 'locked' in by addition of tie wrap through flexion/extension controls.

Plantarflexion Settings: 0, 10, 30, 45 degrees.
Dorsiflexion Settings: 0, 10, 30 degrees

The lightweight ROM walker features detachable foam liners and return-loop hook / loop Fastening. Replacement liners available.

Size	Code	Replacement Liners Code	Shoe Size
XS	*/02	*/02/L	1 – 3
S	*/03	*/03/L	3 – 5.5
M	*/05	*/05/L	6 – 9
L	*/07	*/07/L	9.5 – 11
XL	*/08	*/08/L	11 +

* Insert product code



Podus Boot

Product Code: PB

Complete pressure free heel suspension support, lined with Anti-Decubitus Kodel type liner to prevent pressure sores. The heel features a rotator bar which, when correctly positioned, corrects internal / external rotation of the hip, foot and ankle. The liner is removable for washing and can be replaced. Measure length of foot from back of calcaneus to the tip of the hallux.

Code	NHSCC Code	Item / Size
*/STDS	GSZ11731	Podus Boot Complete Standard, (With Straps)
*/XLS	GSZ11730	Podus Boot Complete Ex Large, (With Straps)
*/STDS-L	-	Replacement Liner for Standard Podus Boot
*/XLS-S	-	Replacement Liner for Ex Large Podus Boot

* Insert product code



Boot Foot Support

Product Code: BFS

A simple foot support to prevent drop foot. Fleece lined and snug fit for added comfort. The lining can be removed for washing or replaced. One size fits all. Only available in Grey.

Code	Item / Size
*G/10	Boot Foot Support Grey
*G/10-L	Replacement Liner Grey

* Insert product code



Universal Plantar Fasciitis Orthosis

Product Code: UFO

Holds foot in dorsi flexion and allows load application to soft tissues of plantar fascia and achilles tendon. Dorsi flexion easy to adjust using double hook/loop straps additional wedges (3 supplied)

Fits left and right. Replacement liners available.

Size	Code	NHSCC Code	Replacement Liners Code	Shoe Size
S	*/S	GSZ11733	*/S-L	3 – 5.5
M	*/M	GSZ657	*/M-L	6 – 9
L	*/L	GSZ658	*/L-L	9.5 – 11

* Insert product code

Heal-up® Protection Boot

Product Code: HUP/10

- Off-loads heel to prevent and manage ulcers
- Simple and effective closure straps
- Floats the heel
- Open heel provides visibility with additional opening for greater inspection visibility
- Universal size, left or right.

Size	Code	NHSCC Codes	Code
Universal	1	GRA5129	HUP/10
Universal	Box of 10	GRA5130	NH/HUP/10

* Insert product code



Treated with **Agion®** and **Silvadur™** antimicrobial to inhibit the growth of bacteria.



Numbered tabs for easy 1 – 2 – 3 application



heal-up.co.uk



Anatomical strapping
with configuration to
inhibit "drop-foot"

Friction free textile



Additional opening
for improved inspection

Floated heel

HEALUP®



Ankle Contracture AFO

Product Code: LEB

A quality AFO constructed from lightweight polycarbonate that provides heel suspension to assist in the treatment and prevention of pressure sores. The clear polycarbonate frame allows the heel to be monitored without the need to remove the orthosis. The AFO is also effective in the management of foot drop and non-fixed plantar flexion contractures.

The liner is treated to wick away moisture and provide a drier cleaner environment thus reducing the risk of bacterial build up and improve patient compliance. There is a padded anti-rotation bar for positioning and prevention of lateral/medial hip rotation. The wide design of the frame also provides lateral ankle protection if patients roll onto their side during sleep. Adjustable toe rest and removable anti-slip sole unit

Size	Code	NHSCC Codes	Size (Liner only)	Code	NHSCC Codes
Small	LEB/S	GRA4917	Small	LEB/S/L	GRA4918
Regular	LEB/R	GRA4934	Regular	LEB/R/L	GRA4916



Leaf Spring Orthosis

Product Code: 825

A lightweight, polypropylene splint which provides static dorsiflexion assistance and allows limited plantarflexion.

The low arch and open heel give a neat shape which fits into a shoe easily, the splint is also heat mouldable.

Size	Code (Left)	Code (Right)	Shoe Size
S	*0/0	*5/0	4 – 6
M	*1/0	*6/0	5 – 7
L	*2/0	*7/0	6 – 8
XL	*3/0	*8/0	7 – 9

* Insert product code



Swedish Ankle Foot Orthosis

Product Code: M.9050

Provides static dorsiflexion assistance and lateral stability for the ankle/foot area. Made in polythene, it is thicker on the vertical aspect for rigidity and thinner on the footplate to allow for trimming.

The low arch and open heel give a neat shape which fits into a shoe easily and the calf section is padded for comfort.

Can be trimmed to size with scissors and shaped with a heat gun if necessary.

Size	Code (Left)	Code (Right)	NHSCC Codes	Shoe Size
S	*/SL	*/SR	-	3 – 6
M	*/ML	*/MR	GRA242	7 – 8.5
L	*/LL	*/LR	-	7.5 – 10

* Insert product code



Trimmable sole plate



100% Pure Carbon Fibre

Size	Left	Right	Shoe Size (EUR)	Length Tolerance +/- 2mm
S	*/SL	*/SR	36 - 39	253 mm
M	*/ML	*/MR	39 - 42	275 mm
L	*/LL	*/LR	42 - 45	293 mm
XL	*/XLL	*/XLR	45 - 47	300 mm

* Insert product code

AFO-1 Classic Medial/KT AFO-1 Classic Medial EasyFit

Product Code: KTAFO1/M

Strong and lightweight solution to provide constant support to the foot. Maximum strength and durability. Sole can be easily adapted to different shoes. Velcro closure, with return loop.

Recommended for use:

- Foot drop can occur after a stroke or traumatic brain injury or as a result of multiple sclerosis, neuromuscular atrophy, or peroneal paralysis;
- Foot and calf extensors after operation;
- Ankle arthritis.

Maximum body weight tolerance is 100 kg.

Premium quality pre-preg carbon fibre 100%.

Available in a right and left.

*Easy Fit means product with trimmable sole part – what gives possibility to put it inside the shoe easier and to make corrections without any grinding.

Removable & replaceable padding



Trimmable sole plate



100% Pure Carbon Fibre

Size	Left	Right	Shoe Size (EUR)	Length Tolerance +/- 2mm
S	*/SL	*/SR	36 - 39	253 mm
M	*/ML	*/MR	39 - 42	275 mm
L	*/LL	*/LR	42 - 45	293 mm
XL	*/XLL	*/XLR	45 - 47	300 mm

* Insert product code

AFO-2 Classic Lateral/KT AFO-1 Classic Lateral EasyFit

Product Code: KTAFO1/L

Strong and lightweight solution to provide constant support to the foot. Maximum strength and durability. Sole can be easily adapted to different shoes. Velcro closure, with return loop.

Recommended for use:

- Foot drop can occur after a stroke or traumatic brain injury or as a result of multiple sclerosis, neuromuscular atrophy, or peroneal paralysis;
- Foot and calf extensors after operation;
- Ankle arthritis.

Maximum body weight tolerance is 100 kg.

Premium quality pre-preg carbon fibre 100%.

Available in a right and left.

*Easy Fit means product with trimmable sole part – what gives possibility to put it inside the shoe easier and to make corrections without any grinding.

Removable & replaceable padding





100% Pure Carbon Fibre

Size	Left	Right	Shoe Size (EUR)	Length Tolerance +/- 2mm
S	*/13	*/23	36 - 39	253 mm
M	*/15	*/25	39 - 42	275 mm
L	*/17	*/27	42 - 45	293 mm

* Insert product code

AFO-3 Flexible/KT AFO-3 Flexible EasyFit

Product Code: PRAFKT3

Spiral structure for resistance and elasticity; pre-loaded to provide constant support to the foot. Light. Maximum strength and durability. Sole can be easy adapted to different shoes. Velcro closure, with return loop.

Recommended for use:


- Neurological illnesses;
- Partial foot amputation;
- Limping foot;
- Multiple sclerosis;
- Ictus.

Maximum body weight tolerance is 100 kg.

Premium quality pre-preg carbon fibre 100%.

Available in a right and left.

*Easy Fit means product with trimmable sole part – what gives possibility to put it inside the shoe easier and to make corrections without any grinding.



100% Pure Carbon Fibre

Size	Left	Right	Shoe Size (EUR)	Length Tolerance +/- 2mm
S	*/SL	*/SR	36 - 39	253 mm
M	*/ML	*/MR	39 - 42	275 mm
L	*/LL	*/LR	42 - 45	293 mm

* Insert product code

AFO-4 Anterior Carbon AFO

Product Code: PREGFIX/AFO

Lightweight and low profile. Features an anterior moulded rigid shell with comfortable padding. Provides dorsiflexion support during swing phase.

Maximum body weight tolerance is 100 kg.

Premium quality pre-preg carbon fibre 100%.

Available in a right and left.

*Easy Fit means product with trimmable sole part – what gives possibility to put it inside the shoe easier and to make corrections without any grinding.





Fixed Ankle

Product Code: 821

A lightweight, polypropylene drop foot splint which provides dorsiflexion assistance and prevents plantarflexion.

The splint is heat mouldable.

Size	Code (Left)	Code (Right)	Shoe Size
S	*0/0	*5/0	4 – 6
M	*1/0	*6/0	5 – 7
L	*2/0	*7/0	6 – 8
XL	*3/0	*8/0	7 – 9

* Insert product code



Universal Foam Ankle

Product Code: UFA

Two plastic medial and lateral ankle supports, shaped to fit comfortably over the malleoli, are lined with closed-cell foam.

The plastic supports are attached to a 'stirrup' of closed cell foam, fixed with hook/loop and worn under the heel.

One size; universal design fits left and right.

Size	Code	NHSCC Codes
Universal	UFA/10	GRA4930



Air/Gel Ankle

Product Code: AGA/10

Two plastic medial and lateral ankle supports, shaped to fit comfortably over the malleoli, are lined with closed-cell foam. The two detachable air/gel bags line the splint and can be used for cold therapy.

The plastic supports are attached to a 'stirrup' of closed cell foam, fixed with hook/loop and worn under the heel.

One size; universal design fits left and right.

Size
Universal



Barnsley Ankle Brace

Product Code: EWB

The stitched forefoot section allows this brace to slip on like a sock, and the hook/loop closure applies variable tension for comfort.

Made in 10cm quality elastic, giving light ankle support.

Contains Latex

Size	Left	Right	Shoe Size
S	*/13	*/23	3 – 5
M	*/15	*/25	5 – 7
L	*/17	*/27	7 – 9

* Insert product code



Figure of 8 Ankle Brace

Product Code: F8LA

A lightweight ankle support which laces up and features an elastic figure of 8 strap for added support and compression. Fits either ankle.

Contains Latex

Size	Code	NHSCC Codes	Shoe Size
XS	*/02		1 – 3
S	*/03	GRA4914	3 – 5
M	*/05	GRA4932	6 – 8
L	*/07	GRA4915	9 – 12
XL	*/08		12 +

* Insert product code



Malleostrong

Product Code: M.A1028

Ankle brace for limiting pro – and supination. Thin breath-o-prene structure with padded medial and lateral shells. Additional rear padding for improved comfort. Circular strap for limiting talar shift (can be shortened to fit)

- Functional therapy of 2nd degree distortions
- Functional recovery following surgical/conservative treatment of ligament injuries (3rd degree distortions)
- Functional recovery of consequences of malleolar fractures (surgically/conservatively treated)

Contains Latex

Size	Code (Left)	NHSCC Codes	Code (Right)	NHSCC Codes	Shoe Size
XS	*/XSL	GRA5256	*/XSR	GRA5257	1 – 3
S	*/SL	GRA4926	*/SR	-	3 – 5
M	*/ML	-	*/MR	GRA4936	6 – 8
L	*/LL	GRA4924	*/LR	GRA253	9 – 12
XL	*/XLL	-	*/XLR	-	12 +

* Insert product code



Hallux Valgus Night Splint

Product Code: HVN

A plastic semi flexible splint with padded closed cell foam toe section. Features adjustable forefoot pad and adjustable hook/loop fastening.

Size	Code (Left)	NHSCC Codes	Code (Right)	NHSCC Codes
XS	*/12	-	*/22	-
S	*/13	GRA327	*/23	GRA330
M	*/15	-	*/25	GRA331
L	*/17	-	*/27	GRZ5344

* Insert product code



Universal Soft Hallux Valgus

Product Code: JSHV

A soft padded universal size splint that helps maintain the large toe in a neutral position. Easy to apply with adjustable hook/loop fastening, adjustable aluminium stay provides stability and assists in the positioning of the alignment of the large toe. The splint can be worn at night or during the day in footwear

Size	Code (Left)	Code (Right)	NHSCC Codes
Universal	JSHV/L	JSHV/R	ERP5492



Anklet

Product Code: ANK

A pull-on Neoprene support for light compression and ankle support.

Size	Code	Shoe Size
S	*/03	3 – 5
M	*/05	5 – 7
L	*/07	7 – 9
XL	*/08	9 – 11

* Insert product code



Lincoln Ankle Brace


Product Code: LAB

A zipped ankle support in Ventoprene, which is ventilated, towelling – lined Neoprene. The zip is worn on the lateral aspect of the ankle, and a flexible spiral stay gives good ankle support on the medial aspect of the ankle.

The 'figure of eight' strapping is designed to limit plantar flexion.

Size	Code (Left)	Code (Right)	Shoe Size
S	*/13	*/23	3 – 5
M	*/15	*/25	5 – 7
L	*/17	*/27	7 – 9

* Insert product code



The coloured graphic depicts the pressure distribution on the sole of the foot.

Jura Post-Op Shoe

Product Code: P-OP



- Square toe configuration offers protection and improved fit/ accommodation.
- Rocker-sole profile aids gait and reduces forces on foot.
- Outsole is non-slip and material composition provides reduced movement under the metatarsal heads.
- Simple closure system allows for additional dressings.
- Insole can be removed for alteration, improved accommodation or fitting 'Off-Loading Insole'
- Removable toe cover available if required

Mens Size	Code	NHSCC Code	Shoe Size	Ladies Size	Code	NHSCC Code	Shoe Size
S	*/MS	ERP2087	5 – 7	S	*/WS	ERP2108	3 – 5
M	*/MM	ERP2086	7.5 – 9	M	*/WM	ERP2107	5.5 – 7
L	*/ML	ERP2085	9.5 – 11	L	*/WL	ERP2088	7.5 – 9
XL	*/XL	ERP2091	11.5 – 13				

* Insert product code

** Insert size code from table above

Product	Code
Removable Toe Cover	P-OP/**/T



Jura High Post-Op Shoe


Product Code: HP-OP

- Square toe configuration offers protection and improved fit/ accommodation.
- Rocker-sole profile aids gait and reduces forces on foot.
- Outsole is non-slip and material composition provides reduced movement under the metatarsal heads.
- Simple closure system allows for additional dressings.
- Insole can be removed for alteration, improved accommodation or fitting 'Off-Loading Insole'
- Adjustable instep strap for improved seating of heel.
- Ankle high upper for improved fit and control for patients with compression dressings and total contact casts.

Mens Size	Code	NHSCC Codes	Shoe Size
S	*/MS	-	5 – 7
M	*/MM	ERP2106	7.5 – 9
L	*/ML	-	9.5 – 11
XL	*/XL		11.5 – 13

Ladies Size	Code	Shoe Size
S	*/WS	3 – 5
M	*/WM	5.5 – 7
L	*/WL	7.5 – 9

* Insert product code



The coloured graphic depicts the pressure distribution on the sole of the foot.






The coloured graphic depicts the pressure distribution on the sole of the foot.

Off-loader Heel Shoe

Product Code: OWH

- Hard wearing padded nylon upper.
- Adjustable hook / loop fastening.
- Hard-wearing sole unit shaped to prevent force on forefoot.
- Extendable forefoot closure panel to accommodate dressings padding.
- Lower profile wedge.
- Insole can be removed for alteration, improved accommodation or fitting 'Off-Loading Insole'

Size	Code	NHSCC Code	Shoe Size
XS	OWH/XS	ERP2026	2 – 4
S	OWH/S	ERP2040	5 – 6
M	OWH/M	ERP2039	7 – 8
L	OWH/L	ERP2038	9 – 11
XL	OWH/XL	ERP2041	12 – 14



Off-loading Insoles

- Easily remove sections of insole for custom off-loading in treatment of various foot pathologies.
- Greatly reduces pressure on plantar surface when patient is walking.
- Self-adhesive stabilising base board.
- Available for all Jura post op shoes.



Off-loading insoles for Post-op & High Post-op Shoes

Mens Size	Code	Shoe Size	Womans Size	Code	Shoe Size
S	P-OP/MSI	5 – 7	S	P-OP/WSI	3 – 5
M	P-OP/MMI	7.5 – 9	M	P-OP/WMI	5.5 – 7
L	P-OP/MLI	9.5 – 11	L	P-OP/WLI	7.5 – 9
XL	P-OP/XLI	11.5 – 13			

Off-loading insoles for Off-loader Heel Shoe

Size	Code	Shoe Size
XS	OWH/XSINS	2 – 4
S	OWH/SINS	5 – 6
M	OWH/MINS	7 – 8
L	OWH/LINS	9 – 11
XL	OWH/XLINS	12 – 14



Heel Relief Shoe

Product Code: HRS





- Alleviates pressure on the heel
- Comfortably padded
- Insole can be removed for alteration or to accommodate a custom orthotic
- Ankle strap reduces heel slippage
- Removable / adjustable closure to accommodate dressings etc

Common Examples of Use

- Trauma to the rear of the foot
- Wounds or ulcerations on the heel area
- Post surgical healing
- Pain or bruising to the heel

Size	Code	Shoe Size (Women's)	Shoe Size (Men's)
XS	*XS	4 – 7.5	-
S	*S	8 – 10	6 – 7.5
M	*M	10.5 – 13	8 – 10
L	*L	-	10.5 – 12
XL	*XL	-	12.5 +

* Insert product code





Hi-Sided Cast Sandal

Product Code: HSCS

Designed with hard wearing, non-stretch nylon faced padded upper for long term wear. Loop-lock D-ring fastening for secure and easy fastening. Non-slip rocker sole.

Also available in paediatric sizes.

Size	Code	NHSCC Codes	Length	Width
XS	*/AXS	-	23 cm	8 cm
S	*/AS	CVL219	24.5 cm	9 cm
M	*/AM	CVL218	26.5 cm	10 cm
L	*/AL	CVL217	29 cm	10.5 cm
XL	*/AXL	-	31 cm	11 cm

* Insert product code



Open Toe Cast Sandal with Buckles

Product Code: CS/30

Designed for long term wear, a strong moulded sole is stitched to a hard-wearing canvas upper. The return loop buckles ensure a good fit on the patient.

For use with patients in a weight bearing cast. Colours may vary.

Size	Code	Sole Length	Sole Width
XS	*XS	18 cm	9 cm
S	*S	23 cm	10 cm
M	*M	27 cm	11.5 cm
L	*L	31 cm	13 cm
XL	*XL	36 cm	14 cm

* Insert product code



Medi Shoe

Product Code: 811

A quality padded upper with return loop buckles accommodates the oedematous or bandaged foot very comfortably. The ridged non-slip sole is semi-flexible to allow normal gait, and lined with closed cell foam for comfort.

Size	Code	Sole Length	Sole Width
XS	*42	20 cm	8 cm
S	*43	23 cm	8 cm
M	*47	27 cm	9 cm
L	*35	29 cm	10 cm
XL	*37	31 cm	12 cm

* Insert product code

Comfortable
padded heel counter

Return Loop
buckles



Durable hard-wearing sole unit

Rigid Post Op Shoe

Product Code: RPO

A padded durable nylon upper with return loop buckles.

The upper will accommodate swelling and dressings.

Padded insole with fully rigid, lightweight hard wearing outer sole unit.

The sole unit is designed for improved stability to aid rehabilitation following foot surgery, trauma or any condition that requires immobilisation to aid healing.

Size	Code	NHSCC Code	Sole Length	Sole Width
Paediatric	*/PAED	ERF769	20 cm	8 cm
S (Ladies)	*/WS	ERF773	24 cm	8 cm
M (Ladies)	*/WM	ERF772	25 cm	9 cm
L (Ladies)	*/WL	ERF771	27 cm	9 cm

Size	Code	NHSCC Code	Sole Length	Sole Width
S (Gents)	*/MS	ERF770	28 cm	10 cm
M (Gents)	*/MM	ERF767	29 cm	10 cm
L (Gents)	*/ML	ERF766	30 cm	10 cm
XL	*/MXL	ERF768	32 cm	11 cm

* Insert product code



Repositionable
hook tabs



Basic Post-Op Shoe

Product Code: BPO

A padded nylon upper with return loop buckles accommodates the oedematous or bandaged foot very comfortably. The rigid non-slip sole is semi flexible to allow normal gait.

Size	Code	Sole Length	Sole Width
XS	*XS	20 cm	8 cm
S	*S	23 cm	8 cm
M	*M	27 cm	9 cm
L	*L	29 cm	10 cm
XL	*XL	31 cm	12 cm

* Insert product code



Open Toe Cast Sandal

Product Code: CS/08

A hard-wearing canvas upper with hook/loop strapping is stitched to a rigid non-slip rocker sole. For use with patients in a weight-bearing cast.

Size	Code	NHSCC Codes	Sole Length	Sole Width
XS	*XS	ERP2079	18 cm	9 cm
S	*S	CVL004	23 cm	10 cm
M	*M	CVL000	27 cm	11.5 cm
L	*L	CVL003	31 cm	13 cm
XL	*XL	ERP2078	36 cm	14 cm

* Insert product code



Open Toe Cast Sandal with Buckles

Product Code: CS/48

Return loop buckles ensure a good fit on the patient. A hard-wearing canvas upper with hook/loop strapping is stitched to a rigid non-slip rocker sole. For use with patients in a weight bearing cast.

Size	Code	Sole Length	Sole Width
XS	*/XS	18 cm	9 cm
S	*/S	23 cm	10 cm
M	*/M	27 cm	11.5 cm
L	*/L	31 cm	13 cm
XL	*/XL	36 cm	14 cm

* Insert product code



Rigid Post-Op Shoe

Product Code: 445

A padded durable nylon upper with return loop buckles. The upper will accommodate swelling and dressings.

Padded insole with fully rigid lightweight composite outer sole unit. The sole unit is rigid to prevent any flexion following foot surgery that requires immobilisation to aid healing and rehabilitation.

Size	Code	Sole Length	Sole Width
S (Ladies)	*/WS	24 cm	8 cm
M (Ladies)	*/WM	25 cm	9 cm
L (Ladies)	*/WL	27 cm	9 cm


Size	Code	Sole Length	Sole Width
S (Gents)	*/MS	28 cm	10 cm
M (Gents)	*/MM	29 cm	10 cm
L (Gents)	*/ML	30 cm	10 cm
XL	*/XL	32 cm	11 cm



Dress Shoe Comfort Inserts

- Visco elastic polymer providing excellent absorption of shock and shear force
- Fits easier into thinner footwear, eg dress shoes
- Unique shaping allows multiple functions
- Scallop pad cushion MT heads
- Heel rim band absorbs and dissipates shock/shear
- Heel cup area aids painful heel and heel spur area


Size	Shoe Size (Ladies)	Code	Shoe Size (Gents)	Code
S	3 – 5	L-DSF1	5 – 7	L-DSM1
M	6 – 7	L-DSF2	8 – 10	L-DSM2
L	8 – 10	L-DSF3	11 – 13	L-DSM3



PQ – 6mm Heel Lift

- Viscoelastic polymer
- Shock absorbing properties
- 6mm raise for shortage accommodation
- sold in single units


Shoe Size	Code
3 – 5	L-PQ61
6 – 8	L-PQ62
9 – 11	L-PQ63



PQ – 9mm Heel Lift

- Viscoelastic polymer
- Shock absorbing properties
- 9mm raise for shortage accommodation
- sold in single units

Shoe Size	Code
3 – 5	L-PQ91
6 – 8	L-PQ92
9 – 11	L-PQ93



Star Heel Pad

- 6mm thick Visco elastic heel pad
- Unique star raise pattern improves dissipation of forces
- Well tapered towards pad edge

Size	Shoe Size	Code
S	3 – 5	L-SHPS
M	6 – 8	L-SHPM
L	9 – 11	L-SHPL



Star Spur Pad

- 6mm thick Viscoelastic heel pad
- Unique star raise pattern improves dissipation of forces
- Well tapered towards pad edge
- Heel relief reduces pain associated with heel spurs and plantar fasciitis

Size	Shoe Size	Code
S	3 – 5	L-SSPS
M	6 – 8	L-SSPM
L	9 – 11	L-SSPL

JURAGel

Jura Gel is a therapeutic range of pure silicone gel containing mineral oil. The products are washable and reusable, manufactured to a high standard to give effective therapy, together with continuous lubrication of the tissue.

Each product is extremely soft and elastic and is comfortable to wear. All products are designed to be worn against the skin and fit easily in the shoe.



Silicone Heel Cup with Soft Heel Spot

- Silicone heel support with stabilising walls
- High viscosity/elasticity cushions loading points when walking/standing
- New design of silicone distributes load more evenly
- Softer silicone area in centre of heel
- Hard-wearing, unbreakable under normal foot loads
- Sold in pairs


Size	Code
S	L-405
M	L-406
L	L-407
XL	L-408



Silicone Pro Insole

- High quality, hard-wearing range for long-term use
- 3/4 length silicone insole with relief pads in heel and metatarsal area
- Sold in pairs

Size	Shoe Size	Code	NHSCC Code
S	3 – 5	L-SRIS	ERP5541
M	6 – 8	L-SRIM	-
L	9 - 12	L-SRIL	-



Epoda Dual Heel Cups

Product Code: EPCB

Dual density design with heel cradle for extra comfort and to aid foot positioning.

Size	Shoe Size	Code	NHSCC Codes
S	4 – 5	*/03	-
M	6 – 7	*/05	ERP2664
L	8 – 9	*/07	ERP2665
XL	10 – 11	*/08	ERP2666

* Insert product code




Epoda Dual Heel Pads

Product Code: EPB

Dual density wedge shaped design with integral extra soft blue area.

Size	Shoe Size	Code (5mm Thick)	Code (10mm Thick)
S	4 – 5	*/53	*/103
M	6 – 7	*/55	*/105
L	8 – 9	*/57	*/107
XL	10 – 11	*/58	*/108


* Insert product code



Cap

- Elasticated protective toecap lined with silicone gel
- Used in the treatment of toe or nail lesions
- Sold singly

Shoe Size	Code	NHSCC Codes
3 – 5	L-455	ERP5534
6 – 8	L-456	-
9 – 12	L-457	ERP5535




Elastic Sheath

- Seamless elasticated tubular protection lined with silicone gel on one side
- Can be cut to size
- Useful in the treatment of callouses, blisters, soft corns and hammer toes
- Sold singly

Shoe Size	Code	NHSCC Codes
3 – 5	L-461	-
6 – 8	L-462	-
9 – 12	L-463	ERP5537

Contains Latex




Tubular Protection A

- Seamless elasticated tubular protection lined with silicone gel on one side
- Can be cut to size
- Useful in the treatment of callouses, blisters, soft corns and hammer toes
- Sold singly

Shoe Size	Code
3 – 5	L-449
6 – 8	L-450
9 – 12	L-451

Contains Latex




Tubular Protection B

- Seamless elasticated tubular protection
- Cut to length lined with silicone gel
- Useful in the treatment of callouses, blisters, soft corns and hammer toes
- Sold singly

Shoe Size	Code
3 – 5	L-443
6 – 8	L-444
9 – 12	L-445

Contains Latex




Silicone Gel Metatarsal Pad with Elastic Band

- Comfortable elastic band with large metatarsal silicone gel pad
- Guaranteed shock absorbing action on the entire plantar areas of the forefoot
- Useful for all types of metatarsal pain
- Sold in single units

Size	Code
S – M	L-473
L – XL	L-474


Contains Latex



Silicone Gel Discs

- Fitted on adhesive fabric for multiple use
- Protects against rubbing or direct pressure on particular areas of the foot
- Sold in pairs

Shoe Size	Code
Universal	L-495



Silicone Gel Retractor / Protector

- Helps to realign the hallux valgus, preventing the first toe from overlapping the second
- Prevents rubbing between toes and reduces pain caused by callouses
- Sold singly


Shoe Size	Code	NHSCC Codes
3 – 5	L-412	ERP5532
6 – 8	L-413	-
9 – 12	L-414	-



Silicone Gel Hallux Valgus Protector

- Pad design helps distribute load about joint
- Reduces and removes pain associated with bursitis of first MP joint
- Sold in single units

Size	Code
S – M	L-439
M – L	L-440



Toe Retractor

- Silicone gel pad avoids overlapping of toes
- Reduces pain caused by callouses between toes
- Sold singly

Shoe Size	Code
3 – 5	L-430
6 – 8	L-431
9 – 12	L-432



Anti-slip Pad

- Self-adhesive silicone gel pad
- Easily attached to inner heel of shoe
- Prevents the shoe from slipping off and protects the heel
- Sold in pairs

Size	Code
Universal	L-478



Maternity Support



- Adaptable to a standing or sitting position and can be worn over clothing. Provides strong circumferential support around the pelvis holding the two sacroiliac joints and the symphysis pubis joint firmly.

- Also available in black

Pelvic Instability Belt

Product Code: SPD

- Designed by obstetric physiotherapists and trialed by expectant mothers, this 'many-tailed bandage' pattern allows variable tension on each set of straps, with strong tension over the pubic area and lighter tension under the abdomen.

Contains Latex

Code (White)	NHSCC Code (White)	Code (Black)	Circumference (At widest part of abdomen)
*1	GRZ1261	*1/BK	90 – 96 cm
*2	GRZ1266	*2/BK	96 – 102 cm
*3	GRZ1267	*3/BK	102 – 108 cm
*4	GRZ1262	*4/BK	108 – 116 cm
*5	GRZ1263	*5/BK	116 – 126 cm
*6	GRZ1264	*6/BK	126 – 136 cm
*7	GRZ1265	*7/BK	136 – 146 cm
*8	GRZ1244	*8/BK	146 – 156 cm

* Insert product code



Wiltshire Maternity Belt

Product Code: WMB

- A pelvic support belt with an elastic back section and two padded side straps, plus two variable tension straps.
- The slimline design of this belt makes it comfortable to wear when sitting.
- Provides strong circumferential support around the pelvis holding the two sacroiliac joints and the symphysis pubis joint firmly.
- Back depth 10cm, front depth 8cm.
- Also available in black.

Contains Latex

Code (White)	NHSCC Code	Code (Black)	NHSCC Code	Circumference
*1	GRZ023	*1/BK	GRZ1240	96 – 116 cm
*2	GRZ024	*2/BK	GRZ1241	116 – 130 cm
*3	GRZ022	*3/BK	GRZ1242	130 – 146 cm
*4	GRZ025	*4/BK	GRZ1243	146 – 160 cm

* Insert product code



Trochanteric Belt

Product Code: TROCH

- A strong webbing strap with foam padding gives firm circumferential support at trochanteric level to ease the pain of separated symphysis.
- Can also be used to treat sport-induced symphysis pubis pain.
- Also available in black.

Size	Code (Navy Blue)	NHSCC Code (Navy Blue)	Code (Black)
One size	*10	GRZ1245	*10/BK

* Insert product code



Maternity Belt Style A

Product Code: MBAS

- Made in quality elastic with oblique-angled hook/loop closure for circumferential support around the pelvis. Depth 20cm back and front. Also available in black.

Contains Latex

Code (White)	NHSCC Code (White)	Code (Black)	Circumference (At widest part of abdomen)
*1	-	*1/BK	90 – 96 cm
*2	-	*2/BK	96 – 102 cm
*3	GPK1508	*3/BK	102 – 108 cm
*4	GPK1509	*4/BK	108 – 116 cm
*5	GPK1510	*5/BK	116 – 126 cm
*6	GPK1524	*6/BK	126 – 136 cm
*7	GPK1525	*7/BK	136 – 146 cm
*8	GPK1526	*8/BK	146 – 156 cm
*9	-	*9/BK	156 – 165 cm
*10	-	*10/BK	165 – 178 cm

* Insert product code



Maternity Belt Style B

Product Code: MBB

- Cut away design at the front gives comfort when sitting and prevents belt from rolling up.
- Angled posterior seams give improved shaping around the sacrum.

Code (White)	NHSCC Code (White)	Code (Black)	Circumference (At widest part of abdomen)
*1	GPK1527	*1/BK	90 – 96 cm
*2	GPK1528	*2/BK	96 – 102 cm
*3	GPK1529	*3/BK	102 – 108 cm
*4	GPK1530	*4/BK	108 – 116 cm
*5	GPK1531	*5/BK	116 – 126 cm
*6	GPK1532	*6/BK	126 – 136 cm
*7	GPK1533	*7/BK	136 – 146 cm
*8	GPK1534	*8/BK	146 – 156 cm
*9	-	*9/BK	156 – 165 cm
*10	-	*10/BK	165 – 178 cm
*11	-	*11/BK	178 – 185 cm
*12	-	*12/BK	185 – 196 cm

* Insert product code

- Made in stronger Leno Weave elastic/cotton mix for extra support.
- Back depth 20cm, front depth 17cm.
- Also available in black.

Contains Latex



Maternity Belt Style C

Product Code: MBC

- The low-cut design of this belt makes it comfortable to wear when sitting.
- Angled posterior seams give improved shaping around the sacrum.

Code (White)	NHSCC Code	Code (Black)	Circumference (At widest part of abdomen)
*1	GRZ027	*1/BK	90 – 96 cm
*2	GRZ028	*2/BK	96 – 102 cm
*3	GRZ029	*3/BK	102 – 108 cm
*4	GRZ030	*4/BK	108 – 116 cm
*5	GRZ026	*5/BK	116 – 126 cm
*6	GRZ031	*6/BK	126 – 136 cm
*7	GRZ032	*7/BK	136 – 146 cm
*8	GRZ033	*8/BK	146 – 156 cm
*9	GRZ1260	*9/BK	156 – 165 cm

* Insert product code

- Made in stronger Leno Weave elastic/cotton mix for extra support.
- Back depth 16cm, front depth 11cm.
- Also available in black.

Contains Latex



Contains Latex

Deluxe Maternity Belt

Product Code: DMB

- Anterior foam laminate panels – low profile shape for improved fit, especially when sitting.
- Hook/loop closure.
- Adjustable stretch fulcrums give additional support.
- Helps relieve lower back pain and associated symptoms.

Code	NHSCC Codes	Circumference (At widest part of abdomen)
*1	-	90 – 96 cm
*2	-	96 – 102 cm
*3	-	102 – 108 cm
*4	GPK1451	108 – 116 cm
*5	GPK1452	116 – 126 cm
*6	GPK1453	126 – 136 cm
*7	GPK1454	136 – 146 cm
*8	GPK1455	146 – 156 cm
*9	GPK1456	156 – 165 cm
*10	GPK1448	165 – 178 cm

* Insert product code



- Additional narrow pelvic wrap provides extra support around the pelvis for improved stability at SI joints.
- Also available in black.

Duofit Maternity Belt

Product Code: DUB

- 20cm deep circumferential elastic wrap with anterior cut-away to improve comfort when sitting.
- Angled hook/loop closure to allow accurate positioning under tummy.

Code (White)	Code (Black)	Circumference (At widest part of abdomen)
*1	*1/BK	90 – 96 cm
*2	*2/BK	96 – 102 cm
*3	*3/BK	102 – 108 cm
*4	*4/BK	108 – 116 cm
*5	*5/BK	116 – 126 cm
*6	*6/BK	126 – 136 cm
*7	*7/BK	136 – 146 cm
*8	*8/BK	146 – 156 cm
*9	*9/BK	156 – 165 cm
*10	*10/BK	165 – 178 cm

* Insert product code



Fracture Bracing & Splinting



Standard Humeral Fracture Brace

Product Code: SHC (full kit) SHS (sleeve only)

A polypropylene humeral sleeve with closed cell foam lining, and a polythene insert for comfort and improved fit. The sleeve can be cut with bandage shears and is heat mouldable. Kit includes two humeral cast socks and a fully adjustable sling.

Size	Code (Full Kit)	Code (Sleeve Only)	Length	Circumference of upper arm	
				Distal	Proximal
S	SHC/03	SHS/03	15 cm	18 – 23 cm	23 – 28 cm
M	SHC/05	SHS/05	17 cm	23 – 28 cm	28 – 33 cm
L	SHC/07	SHS/07	18 cm	28 – 33 cm	33 – 38 cm
XL	SHC/08	SHS/08	23 cm	33 – 38 cm	38 – 43 cm

Sling strap

Sling Strap is a non-stretch padded 50mm wide material that is combined with Hook Y-Tabs to provide an adjustable and versatile sling system.

Sling Strap can be cut to the required length to produce a simple collar n' cuff or similar type of sling. This reduces waste and the adjustability of the material ensures that the sling will fit and can be easily adjusted. Sling Strap is padded for comfort around the neck and the modern black colour improves patient compliance and looks better for longer. Sling Strap is supplied in a pack that contains: 2 x 12m or 2 x 6m rolls and 20 Y-tab closures.



Product Code	Length	NHSCC Code	Width
SS	2 x 12 Metres	GRZ554	5 cm
SS/5.5	2 x 6 Metres	ESM137	5.5cm



Humeral Stabilising Orthosis

Product Code: **HSO/10/K** (full kit)
HSO/10 (sleeve only)

A unique lightweight, comfortable humeral sleeve with a ventilated 'skin friendly' protective lining, which means there is no humeral sock required.

The universal humeral is fully customisable, the thermo formable external sleeve can be altered using standard scissors (A sizing guide is included with sleeve).

An optional kit is available which ensures immobilisation of the arm, eliminating swinging of the arm. Kit includes a Deluxe Shoulder Immobiliser, product code: DSI (See pg. 33).

Code	Description
HSO/10/K	Full Kit
HSO/10	Sleeve Only
HSO/10/E	Sleeve Only (Extended)



Promedics Hip Joint

Product Code: **PHJ/10**

The Promedics Hip Joint is 40cm in length and allows control of abduction (between - 5 to +25°) with variable flexion and extension stops. The abduction axis may be set at various positions while the brace is on the patient, by simply turning the central adjustment screw with the allen key provided.

Universal Hip Abduction Waist Belt

Product Code: **JBB/10**

Universal Hip Abduction Belt with Shoulder Strap

Product Code: **JBBS/10**



Soft Humeral


Product Code: **PSH**

A lightweight preformed brace made from a tri-lamination with outer fabric facing, semi-flexible foam and inner smart-textile lining.

- Inner lining has CoolMax treatment to reduce perspiration and Ion/silver treatment to improve temperature control within the brace.
- The kit is supplied with lightweight humeral sock and fully adjustable hook/loop strap closures.
- If required, the brace can be trimmed.

Size	Code	NHSCC Code	Length	Proximal Circ	Distal Circ
XS	*/02	GMB090	28 cm	20 – 23 cm	15 – 19 cm
S	*/03	GMB091	30 cm	23 – 28 cm	19 – 23 cm
M	*/05	GMB092	32 cm	28 – 33 cm	23 – 28 cm
L	*/07	GMB093	34 cm	33 – 38 cm	28 – 33 cm
XL	*/08	GMB094	35 cm	38 – 43 cm	33 – 38 cm

* Insert product code



- A geared, polycentric hinge which can be set to limit the range of flexion and extension at the knee.
- Available in offsets of 12mm, 18mm and 32mm, with strong aluminium side steels which can be shaped with a bending iron.
- Detachable headplates allow easy removal of the hinges from the cast brace.
- Order hinges in pairs.



Genurange Fracture Bracing Hinges

Product Code: GNF

Order Genurange Range of Motion Hinges in pairs. Our Genurange fracture brace hinges are available in 3 sizes for medial and lateral offset

- The geared polycentric hinge can be set to limit range of flexion and extension at the knee
- Strong aluminium side steels which can be shaped with a bending iron
- Detachable headplates allow easy removal of the hinges from the cast braces
- Supplied in pairs

Size	Medial Offset	Lateral Offset
*12/R	12 mm	12 mm
*12/L	12 mm	12 mm
Size	Medial Offset	Lateral Offset
*18/R	18 mm	12 mm
*18/L	12 mm	18 mm
Size	Medial Offset	Lateral Offset
*32/R	32 mm	12 mm
*32/L	12 mm	32 mm

* Insert product code



Replacement Hinge Cover Screws

Product Code: GEN/CS

Replacement screws for Genurange hinge cover.

Replacement Headplates

Product Code: GNP/10

Replacement headplates for Genurange hinges. Detachable for easy removal of the cast brace; supplied in sets of 4.

Replacement Screws

Product Code: GNS/10

Replacement screws to fix headplates on the Genurange hinge; supplied in sets of 4.





Derby Heel Cup

Product Code: DHC

Made in polypropylene plastic with straight or butterfly stems, to fix the heel cup in plaster or synthetic bandage and prevent cast migration.

The hinged ankle gives ankle stability, allowing plantar/dorsiflexion but limiting medial/lateral movement.

Universal design fits left and right.

Size	Code	To fit shoe size
Universal/Standard (Butterfly Stems)	*1	3 – 6
Universal/Large (Butterfly Stems)	*2	6 – 9
Universal/Standard (Straight Stems)	*3	3 – 6
Universal/Large (Straight Stems)	*4	6 – 9

* Insert product code



Femoral Sock

Product Code: FEM

- Full length femoral sock
- Acrylic/lycra fabric
- Range of sizes for individual fitting
- Available singly
- Wicking properties help reduce perspiration

Size	Code	Dimensions
1	*1	91 x 15 cm
2	*2	89 x 23 cm
3	*3	99 x 23 cm
4	*4	107 x 23 cm
5	*5	117 x 23 cm

* Insert product code

Humeral Sock

Product Code: HS/10 (Universal)

- Acrylic/lycra fabric
- Universal fitting
- Available in pairs only
- Wicking properties help reduce perspiration

Tibial Sock

Product Code: BK/10 (Universal)

- Below knee tibial sock
- Acrylic/lycra fabric
- Universal fitting
- Available singly
- Wicking properties help reduce perspiration

Goniometer

Perfect for precision finger measurements and great for measuring difficult finger angles.

Code
GON/10





Discover more in our
new 2022 Hand Therapy
Catalogue

Available to download today from
www.promedics.co.uk



ORFIT® COLOURS NS
NOW AVAILABLE
IN PRECUTS

Full details on Page 8 of the Hand Therapy Catalogue



Loop Strap

Product Code: LS

Loop strap is a hard-wearing strap that will accommodate multiple hook closures. Loop Strap is hook receptive on both sides and can be cut to the required length. Loop Strap really is a hard-wearing, long lasting strap solution that will have multiple applications in the splinting workshop. Loop Strap will last and last!

Width	Length	Code
20 cm	10 Meter roll	LS/20
25 cm	10 Meter roll	LS/25
30 cm	10 Meter roll	LS/30
40 cm	10 Meter roll	LS/40
50 cm	10 Meter roll	LS/50



Terry Strap

Product Code: TS

Terry Strap is a padded, adhesive backed, terry towel lined material that is available in three widths. 3mm thick Terry Strap is a versatile material that is effective when used to line splints or to pad the edges of casts or splints. The structure of the material allows for a limited amount of stretching and shaping to help when forming it to the edges or inner surface of a splint or cast. This is a really useful material to have in the plaster room and the splinting workshop and it will greatly improve the comfort and finish of any splint or cast.

Width	Length	Code
3 cm	15 Metre roll	TS/30
5 cm	15 Metre roll	TS/50
8 cm	15 Metre roll	TS/80



SL Iron

Product Code: NC/79615

Small, narrow iron is perfect for applying iron on tapes and fasteners. Low-profile, trowel shaped iron is ideal for sealing narrow seam tape and fasteners to neoprene. Heats from 66 degrees to 199 degrees in 90 seconds. Blinking light alerts to when ready. Temperature dial is on the handle for easy reference. Includes a sock cover and a sturdy metal stand for safe storage.

Uses standard 120 volt wall outlet. Total length of the iron measures.

Code	Description
NC/79615	SL Iron



Steinel Heat Gun

Product Code: HT/GUN

A 240 Volt heat gun, with variable temperature control between 300 degrees and 525 degrees, and 2-speed airflow control. Lightweight and free standing, designed for use with all types of thermoplastics.

Code	Description
HT/GUN	240 Volt Heat Gun



Leader in thermoplastic innovations

Based on innovative technologies, Orfit develops and produces the most precise and reliable thermoplastic materials for medical devices that improve patient treatment around the world.

Orthotic fabrication materials from Orfit offer the highest possible performance to both therapists and patients in physical rehabilitation.

Our wide product range consists of engineered thermoplastics that will help you solving your most challenging splinting requirements. Through workshops and courses, our experts provide lots of tips and tricks to make your orthotic fabrications functional and successful.

All of our thermoplastic materials are developed and produced to make sure that they are perfectly suited to their purpose and your needs.



Orfilight (NS)

- Immobilization systems for cancer patients in Radiation Oncology
- Orthotic fabrication materials for patients in Physical Rehabilitation
- Prosthetic socket materials for amputee patients

www.orfit.com



Orficast



Orfit Classic



AquaFit NS



Orfit Colours NS



Orfit Eco (NS)



**Dynasyst -
Dynamic Splinting**

**For further details on any Orfit product,
please contact Customer Services:**

**Call: 01475 746 400 or
E-mail: enquiries@promedics.co.uk**



Standard Hook / Loop

Product Code: HAL

- 25m reels of plain nylon hook and loop materials
- Available in White, Black, Navy Blue, Royal Blue, Beige, Red, Dayglo Pink, Dayglo Yellow, Olive Green
- A stock matrix is available upon request
- When ordering please specify colour and width, hook or loop

Code (Hook)	Code (Loop)	Width
*15/H/RL	*15/L/RL	15 mm
*20/H/RL	*20/L/RL	20 mm
*25/H/RL	*25/L/RL	25 mm
*30/H/RL	*30/L/RL	30 mm
*50/H/RL	*50/L/RL	50 mm

* Insert product code



Self Adhesive Hook / Loop

Product Code: SA

- 25m reels of adhesive backed nylon hook and loop materials
- Available in White, Black, Navy Blue, Royal Blue, Red, Dayglo Pink, Dayglo Yellow
- A stock matrix is available upon request
- When ordering please specify colour and width, hook or loop

Code (Hook)	Code (Loop)	Width
*H/15/RL	*L/15/RL	15 mm
*H/20/RL	*L/20/RL	20 mm
*H/25/RL	*L/25/RL	25 mm
*H/30/RL	*L/30/RL	30 mm
*H/50/RL	*L/50/RL	50 mm

* Insert product code

Hook / Loop Colours



Black



White



Navy Blue



Beige



Red



Olive Green



Royal Blue



Dayglo Yellow



Dayglo Pink



Juraprene Splinting Material

Juraprene is a welcome addition to Promedics splinting range. A latex free neoprene material which can be used to create custom splints. Juraprene offers creative splinting options while providing warmth, compression, comfort and support.

Features of Juraprene include:

- 3mm thickness to provide light, restrictive support and compression as required.
- Plush lined neoprene for improved comfort and to reduce pressure points.
- Latex Free
- Machine Washable
- Sheet size 45 x 61 cm



Code	Description
JPR/BEI	Beige
JPR/PUR	Purple
JPR/RBL	Royal Blue
JPR/RED	Red
JPR/BLK	Black



Orfit Dry Heater

Product Code: 32503/230UK (Dry Heater)

Portable heating device for the activation of thermoplastic orthotic fabrication materials.

Features:

- * Full power at your fingertips.
- * Fast activation of your orthotic fabrication material
- * Smart casing. Portable case that makes your work easy on the go.
- * Effective heating control. Excellent temperature control and temperature uniformity with a digital display and integrated timer.

Voltage: 220-240V

Dimensions: 455 mm L x 320 mm W x 30 mm H

Code	Description
32503/230UK	Dry Heater
32501/4	Individual Teflon sheet



**Leader in
thermoplastic
innovations**

Benefits:

- Short warming up time and very quick set-up.
- Quick activation.
- No water necessary.
- Compatible with different electrical grids.
- Thermoplastic splint cools down more slowly when dry heated.
- Compact device
- Mobile and small device.
- Can be used for various thermoplastic sheet materials and precuts.
- Safe to use.
- Hygienic.



Splinting Baths

Product Code: WT BATH

For working with heat activated thermoplastic splinting material.

- Stainless steel
- Removable lid
- Perforated tray
- Thermostatically controlled
- Over-heat protection
- Illuminated on/off switch
- Insulated handles for moving the bath



5 Litre



9 Litre



28 Litre



56 Litre

Size	Code	Internal (mm)	External (mm)	Weight (Kg)	Capacity (Litres)	Features
XXS	WTBATH/XXS	65 x 235 x 300	135 x 361 x 332	5	5	Single lid, Siphon drainage
XS	WTBATH/XS	150 x 300 x 240	400 x 335 x 275	5.9	9	Single lid, Siphon drainage
S	WTBATH/S	200 x 500 x 300	540 x 335 x 290	10.8	28	Single lid, Drainage tap
L	WTBATH/L	200 x 620 x 500	260 x 660 x 540	18.3	56	Double lid, Built-in tap



Hot & Cold Therapy



Knee Wrap

Shoulder Wrap



Universal Cold Therapy Wraps

Product Code: CTKW/10 (Knee) CTSW/10 (Shoulder)

The Promedics Cold Therapy Wraps feature two reusable 12 x 24 cm Hot/Cold packs that slip into the internal pouches. Hook and loop secures hot/cold packs to inside of wrap. Latex Free, Washable, and Reusable. Universal size, fits left or right.

Code	Description	Size
CTKW/10	Knee Wrap	Universal
CTSW/10	Shoulder Wrap	Universal



Reusable heat pouches



Internal pouches secured by hook and loop

- ⊕ Washable
- ⊕ Latex Free
- ⊕ Universal Sizing
- ⊕ Reusable
- ⊕ Two hot/cold packs included
- ⊕ Hook and loop secures internal pouches

HCP/L
Reusable
Hot/Cold
Pack

Jura Hot/Cold Low Back Wrap

The Jura Hot/Cold Low Back Wrap features a large reusable 18 x 25 cm Hot/Cold pack (Code: HCP/L) that slips into the washable smooth nylon pouch. Hook and loop secures the hot/cold pack to inside of wrap. Dual hook straps run lengthwise edge to edge of the hot/cold pack pouch to prevent movement of the pouch and help easy application to a sore back. Perforated panels ventilate the side sections of the body where heat is not usually required. Latex Free, Washable, and Reusable.

One reusable Hot/Cold pack (Code: HCP/L) included.

Code	Size	Waist Circumference
JHC/BW1	S - M	< 35"
JHC/BW2	L - XL	35" - 54"

HCP/S
Reusable
Hot/Cold
Pack

Jura Hot/Cold Ankle Wrap

The Jura Hot/Cold Ankle Wrap features two reusable 11 x 18cm Hot/Cold packs (Code: HCP/S) that slip into the washable smooth nylon pouches (See diagram below). Hook and loop secures hot/cold packs to inside of wrap. Latex Free, Washable, and Reusable.

2 x Small Hot/Cold packs (Code: HCP/S) included.

Code	Size	Shoe Size
JHC/AW1	S - M	2 - 9
JHC/AW2	L - XL	10 +

HCP/S
Reusable
Hot/Cold
Pack

Jura Hot/Cold Elbow Wrap

The Jura Hot/Cold Elbow Wrap features two reusable 11 x 18cm Hot/Cold packs (Code: HCP/S) that slip into the washable smooth nylon pouches. Hook and loop secures hot/cold packs to inside of wrap. Latex Free, Washable, and Reusable.

2 x small Hot/Cold packs (Code: HCP/S) included.

Code	Size	Arm Circumference
JHC/EW1	XS - S	22 - 25 cm
JHC/EW2	M - L	25 - 35.5 cm



Hot/Cold Reusable Packs

Product Code: HCP

Two different sizes of reusable Hot/Cold packs that slip into the washable smooth nylon pouches. Hook and loop secures hot/cold packs to inside of wrap. Latex Free, Washable, and Reusable.

Code	Size	Dimensions
*/S	S	11 x 18 cm
*/L	L	18 x 25 cm

* Insert product code





Paediatric



Lycrafleece™

Paediatric Thumb Spica

Product Code: PLTS

- Manufactured in Lycrafleece material
- Colourful hook / loop design
- Flexible splint to support thumb and provide greater comfort with thumb movement

Size	Code (Left)	NHSCC Code	Code (Right)	NHSCC Code	Wrist Circumference
XS	*/12	GRZ1810	*/22	GRZ1815	12 – 13 cm
S	*/13	-	*/23	-	13 – 14 cm
M	*/15	-	*/25	-	14 – 15 cm
L	*/17	-	*/27	-	15 – 17 cm
XL	*/18	-	*/28	-	17 – 18 cm

* Insert product code



Lycrafleece™

Paediatric Wrist Splint

Product Code: PLWB

- Manufactured in Lycrafleece material
- Colourful hook / loop design
- 15cm depth
- Sold in single units

Size	Code (Left)	NHSCC Code	Code (Right)	NHSCC Code	Wrist Circumference
XS	*/12	GRZ862	*/22	GRZ863	12 – 13 cm
S	*/13	GRZ864	*/23	GRZ865	13 – 14 cm
M	*/15	GRZ866	*/25	GRZ867	14 – 15 cm
L	*/17	GRZ868	*/27	GRZ869	15 – 17 cm
XL	*/18	GRZ870	*/28	GRZ871	17 – 18 cm

* Insert product code



Lycrafleece™

Paediatric Wrist/Thumb Splint

Product Code: PLWTB

- Manufactured in Lycrafleece material
- Colourful hook / loop design
- 15cm depth
- Sold in single units

Size	Code (Left)	NHSCC Code	Code (Right)	NHSCC Code	Wrist Circumference
XS	*/12	GRZ872	*/22	GRZ873	12 – 13 cm
S	*/13	GRZ874	*/23	GRZ875	13 – 14 cm
M	*/15	GRZ876	*/25	GRZ877	14 – 15 cm
L	*/17	GRZ878	*/27	GRZ879	15 – 17 cm
XL	*/18	GRZ880	*/28	GRZ881	17 – 18 cm

* Insert product code



Contains Latex

Paediatric Elastic Wrist Brace

Product Code: PEW

Manufactured in single compression elastic. Contrasting colourful bar casing and edge binding. 15cm length. Sold in single units.

Size	Code (Left)	Code (Right)	Wrist Circumference
1	*/13	*/23	11.5 cm
2	*/15	*/25	14.5 cm
3	*/17	*/27	16 cm

* Insert product code



Contains Latex

Paediatric Elastic Wrist/Thumb Brace

Product Code: PEWT

Manufactured in single compression blue elastic. Contrasting colourful bar casing and edge binding. 15cm length. Sold in single units.

Size	Code (Left)	Code (Right)	Wrist Circumference
1	*/13	*/23	11.5 cm
2	*/15	*/25	14.5 cm
3	*/17	*/27	16 cm

* Insert product code



Paediatric Neoprene Wrist Brace

Product Code: PNW

- Manufactured in a range of coloured Neoprene.
- Contrasting colourful bar casing and hook/loop fastening
- 15cm length
- Sold in single units

Red = **RD** Black = **BK**

Navy = **NY** Pink = **PK**

Purple = **PU** Royal Blue = **RB**

Selected sizes available NHS Supply Chain

Size	Colour	Code (Left)	Code (Right)	Wrist Circumference
XXS	Black	PNW/11/BK	PNW/21/BK	9 cm
XS	Black	PNW/12/BK	PNW/22/BK	10.5
S	Black	PNW/13/BK	PNW/23/BK	11.5 cm
M	Black	PNW/15/BK	PNW/25/BK	14.5 cm
L	Black	PNW/17/BK	PNW/27/BK	16 cm
XXS	Navy	PNW/11/NY	PNW/21/NY	9 cm
XS	Navy	PNW/12/NY	PNW/22/NY	10.5
S	Navy	PNW/13/NY	PNW/23/NY	11.5 cm
M	Navy	PNW/15/NY	PNW/25/NY	14.5 cm
L	Navy	PNW/17/NY	PNW/27/NY	16 cm
XXS	Pink	PNW/11/PK	PNW/21/PK	9 cm
XS	Pink	PNW/12/PK	PNW/22/PK	10.5
S	Pink	PNW/13/PK	PNW/23/PK	11.5 cm
M	Pink	PNW/15/PK	PNW/25/PK	14.5 cm
L	Pink	PNW/17/PK	PNW/27/PK	16 cm
XXS	Purple	PNW/11/PU	PNW/21/PU	9 cm
XS	Purple	PNW/12/PU	PNW/22/PU	10.5
S	Purple	PNW/13/PU	PNW/23/PU	11.5 cm
M	Purple	PNW/15/PU	PNW/25/PU	14.5 cm
L	Purple	PNW/17/PU	PNW/27/PU	16 cm
XXS	Royal Blue	PNW/11/RB	PNW/21/RB	9 cm
XS	Royal Blue	PNW/12/RB	PNW/22/RB	10.5
S	Royal Blue	PNW/13/RB	PNW/23/RB	11.5 cm
M	Royal Blue	PNW/15/RB	PNW/25/RB	14.5 cm
L	Royal Blue	PNW/17/RB	PNW/27/RB	16 cm
XXS	Red	PNW/11/RD	PNW/21/RD	9 cm
XS	Red	PNW/12/RD	PNW/22/RD	10.5
S	Red	PNW/13/RD	PNW/23/RD	11.5 cm
M	Red	PNW/15/RD	PNW/25/RD	14.5 cm
L	Red	PNW/17/RD	PNW/27/RD	16 cm



Neoprene Wrist/Thumb Brace

Product Code: PNWT

- Manufactured in a range of coloured Neoprene.
- Contrasting colourful bar casing and hook / loop fastening
- 15cm length
- Sold in single units

Black = **BK** Navy = **NY** Pink = **PK** Purple = **PU** Royal Blue = **RB** Red = **RD**

Size	Colour	Code (Left)	Code (Right)	Wrist Circumference
1	Black	PNWT/13/BK	PNWT/23/BK	11.5 cm
2	Black	PNWT/15/BK	PNWT/25/BK	14.5 cm
3	Black	PNWT/17/BK	PNWT/27/BK	16 cm
1	Navy	PNWT/13/NY	PNWT/23/NY	11.5 cm
2	Navy	PNWT/15/NY	PNWT/25/NY	14.5 cm
3	Navy	PNWT/17/NY	PNWT/27/NY	16 cm
1	Pink	PNWT/13/PK	PNWT/23/PK	11.5 cm
2	Pink	PNWT/15/PK	PNWT/25/PK	14.5 cm
3	Pink	PNWT/17/PK	PNWT/27/PK	16 cm
1	Purple	PNWT/13/PU	PNWT/23/PU	11.5 cm
2	Purple	PNWT/15/PU	PNWT/25/PU	14.5 cm
3	Purple	PNWT/17/PU	PNWT/27/PU	16 cm
1	Royal Blue	PNWT/13/RB	PNWT/23/RB	11.5 cm
2	Royal Blue	PNWT/15/RB	PNWT/25/RB	14.5 cm
3	Royal Blue	PNWT/17/RB	PNWT/27/RB	16 cm
1	Red	PNWT/13/RD	PNWT/23/RD	11.5 cm
2	Red	PNWT/15/RD	PNWT/25/RD	14.5 cm
3	Red	PNWT/17/RD	PNWT/27/RD	16 cm



Junior Wrist/Thumb Brace

Product Code: JWTO

A wrist brace made of Neoprene which can be cut without fraying to fit around swollen joints, giving compression and excellent conformity.

- Incorporating a malleable metal thumb stay to control an unstable CMC joint; the thumb section can be trimmed to allow IP flexion
- Use for: Juvenile chronic arthritis, soft tissue injury.

Size	Code (Left)	Code (Right)	Wrist Circumference	Palmar Bar Length
S	*/13	*/23	12 cm	14 cm
M	*/15	*/25	14 cm	14 cm
L	*/17	*/27	16 cm	16 cm

* Insert product code



Tracy Brace

Product Code: TWB

A wrist brace made of quality woven elastic for comfort and durability, enclosing an aluminium palmar bar which can be easily removed.

- The design is cut below the palmar crease and will not impinge on the MCP joints allowing for a fully functional hand
- Use for: Colles' fracture, wrist fracture and arthritis
- Length of wrist brace is 14-15cm

Size	Code (Left)	Code (Right)	Wrist Circumference
S	*/13	*/23	12 cm
M	*/15	*/25	14 cm
L	*/17	*/27	16 cm

* Insert product code

Contains Latex

Paediatric Wrist Braces

Our new Paediatric Wrist Braces are now available for a wide range of injuries including Buckle Fractures, Green Stick Fractures and post injury.



Product Code: **CDRS**

Product Code: **CDRL**

Product Code: **BFB**



Children's D Ring Brace Short

Product Code: **CDRS** (Available Patterns - C, H, I & J)

Size	Code (Left)	Code (Right)	Length	Wrist Circ.
S	CDRS/13*	CDRS/23*	12.5 cm	12 - 13 cm
M	CDRS/15*	CDRS/25*	15 cm	14 - 15 cm
L	CDRS/17*	CDRS/27*	18 cm	16 - 17 cm
XL	CDRS/18*	CDRS/28*	18 cm	18 - 19 cm

*Add /C, /H, /I or /J after above code

Fabrics available:



Code

C



Code

E



Code


F

Children's D Ring Brace Long

Product Code: **CDRL** (Available Patterns - E, F, K & L)


Size	Code (Left)	Code (Right)	Length	Wrist Circ.
S	CDRL/13*	CDRL/23*	12.5 cm	12 - 13 cm
M	CDRL/15*	CDRL/25*	15 cm	14 - 15 cm
L	CDRL/17*	CDRL/27*	18 cm	16 - 17 cm
XL	CDRL/18*	CDRL/28*	18 cm	18 - 19 cm

*Add /E, /F, /K or /L after above code



Code

G



Code

H




Code

I



Code

J



Code

K



Code

L

Children's Buckle Fracture Brace

Product Code: **BFB**

Size	Code (Left)	Code (Right)	Wrist Circ.
S	BFB/13	BFB/23	16 - 18 cm
M	BFB/15	BFB/25	18 - 20 cm
L	BFB/17	BFB/27	20 - 22 cm

Please note that pattern codes C, G, H, I & J may vary



Paediatric Ulnar Deviation Support

Product Code: PUDS

- Manufactured in colourful Neoprene.
- Contrasting multicoloured hook/loop closure straps.
- Hard wearing palmar bar casing with adjustable aluminium palmar bar.
- Adjustable ulnar support can be shaped to help control ulnar drift.

Size	Code (Left)	Code (Right)	Wrist Circumference
Baby	*/11	*/21	6.5 cm
XS	*/12	*/22	10 cm
S	*/13	*/23	12 cm
M	*/15	*/25	14 cm
L	*/17	*/27	16 cm


* Insert product code



Paediatric Thumb Splint

Product Code: PTS

Paediatric thumb abduction and opposition for patients with cerebral palsy, hemiplegia, indwelling thumb and other conditions such as congenital clasped thumb.



Paediatric sizes are suitable for children aged approx 3 months to 10 years. Older children may require adult sizes. Measure the circumference of the distal thumb joint (A). Measure the distance on the palm from the distal thumb joint to the centre of the wrist crease (B). Refer to chart provided.

Size	A (cm)	B (cm)
*/1	3 - 3.8	3 - 4
*/2	3.4 - 4	3.6 - 5
*/3	3.7 - 4.5	4 - 6.3
*/4	4 - 5	4.5 - 6.9
*/5	4.3 - 5.5	5.5 - 7.9
*/6	4.8 - 6	6.2 - 9

* Insert product code



Paediatric Lancaster Shoulder Brace

Product Code: LSSB/10 (baby) LSSP (Paediatric)

A universal adjustable sling with shoulder strap and body wrap for shoulder immobilisation.

- Arm envelope can be easily trimmed for custom-fit.
- Fully adjustable shoulder strap with padded sleeve for comfort around neck.
- This sling can be easily adjusted to fit most paediatric patients.

Size	Code	Age
Baby	LSSB/10	0 - 18 Months
Paediatric	LSSP/10	2 - 10 Years



Paediatric Trauma Sling

Product Code: TAU/P

The Trauma Sling is a high quality universal sling that has multiple uses within the healthcare field. The Trauma Sling can be used as a traditional sling or for arm elevation following trauma, burns or post-operative management of the shoulder or upper arm.

The design and quality of materials used offer protection, comfort and ease of access to the damaged limb. The Trauma Sling is easy to apply and is supplied vacuum packed to reduce storage volume.

Size
One size



Fun Pattern

Arm & Knee Immobiliser

Product Code: FPI

- Hard wearing fun print material
- Wrap-around with hook / loop closures
- Three rigid stays for effective immobilisation
- Please quote required length of brace after product code
- After each code please specify choice of fabric from below.
- Patterns may vary, but will be age/gender appropriate.

Size	Code	Length	Circumference (Top)	Circumference (Bottom)
1	FPI8/02/*	20 cm	17 – 24 cm	12 – 19 cm
2	FPI8/03/*	20 cm	23 – 30 cm	14 – 21 cm
3	FPI10/02/*	20 cm	17 – 24 cm	12 – 19 cm
4	FPI10/03/*	25 cm	23 – 30 cm	14 – 21 cm
5	FPI12/03/*	30 cm	23 – 30 cm	14 – 21 cm
6	FPI12/05/*	30 cm	30 – 37 cm	17 – 24 cm
7	FPI14/05/*	35 cm	30 – 37 cm	17 – 24 cm
8	FPI16/07/*	40 cm	37 – 44 cm	20 – 27 cm
9	FPI18/07/*	45 cm	37 – 44 cm	20 – 27 cm
10	FPI18/08/*	45 cm	44 – 51 cm	23 – 30 cm
11	FPI20/07/*	50 cm	37 – 44 cm	20 – 27 cm
12	FPI20/08/*	50 cm	44 – 51 cm	23 – 30 cm

* Add one of the codes below for fabric selection:
(e.g. Size 2, with fabric code 'K' = FPI8/03/K)

Fabrics available:



Please note that patterns may vary. Check on our website:
www.promedics.co.uk





Paediatric Hi-Sided Cast Sandal

Product Code: HSCSL

Designed with hard wearing, non-stretch nylon faced padded upper for long term wear. Loop-lock D-ring fastening for secure and easy fastening. Non-slip rocker sole.

Size	Code	NHSCC Codes	Length	Width
S	*/CS	-	14 cm	5.5 cm
M	*/CM	ERP5518	17 cm	6.5 cm
L	*/CL	ERP5505	20 cm	7.5 cm

* Insert product code



Paediatric Cast Sandal

Product Code: CS

- Strong and durable
- Canvas upper with jolly print
- Grooved non-slip rocker sole unit
- Hook / loop fastening

Size	Code	Length	Width
XXXS	*/98	12.5 cm	9 cm
XXS	*/99	14 cm	9 cm

* Insert product code



Paediatric Walker

Product Code: PJWF

A comfortable paediatric walker boot with adjustable uprights for custom fit on lower limb, providing improved stability and immobilisation.

New lightweight design sole unit to greatly improve shock absorption and gait pattern. Robust foot section for improved wear.

Simple, secure hook / loop fastening. Comfortable detachable liners with replacements available.

Product Features

- Lightweight & durable design
- Improved stability and support
- No-slip, rocker sole

Size	Code	NHSCC Code	Child UK Shoe Size
S	PJWF/03	GAA104	2 – 6
M	PJWF/05	GAA105	7 – 10
L	PJWF/07	GAA106	11 – 13

No-slip tread



Rocker Sole



Paediatric Soft Foam Collar

Product Code: JBC

Made of fully washable, closed cell foam with a soft stockinette covering. Use to stabilise the neck after trauma, or for long term use as a night resting collar.

Size	Code	NHSCC Codes	Depth x Length
S	*1	GQB1208	5 x 30 cm
M	*2	GQB1209	5 x 35 cm
L	*3	GQB1210	6.5 x 37.5 cm

* Insert product code





Protem Belt

Product Code: TBP

This paediatric lumbar support has been used effectively in the management of abnormal tone by providing proximal stability and improving proprioceptive feedback.

Contains Latex

Size	Code	Waist Circumference
XXS	*/01	46 – 54 cm
XS	*/02	56 – 64 cm

* Insert product code



Pavlik Harness

Product Code: PAV

To control hip flexion and abduction in the primary treatment of congenital dislocation of the hips. Colourful strap identification with adjustable hook and loop closures, a padded harness with a unique reversible toe cover which can be worn open or closed toe.

Size	Code	NHSCC Code	Chest Circumference	Age of Baby
XXS	*/01	GRQ022	28 – 34 cm	Prem
XS	*/02	GRZ2048	30 – 36 cm	0 – 3 Months
S	*/03	GRZ2049	36 – 42 cm	3 – 6 Months
M	*/05	GRZ2050	39 – 45 cm	6 – 9 Months
L	*/07	GRZ2051	46 – 55 cm	9 – 12 Months
XL	*/08	-	56 – 59 cm	12 Months +

* Insert product code



CDH Splint

Product Code: FCDH (Standard) ACDH (Adjustable)

A lightweight, rigid polypropylene frame is lined with a hygienic, non-absorbent, closed cell foam padding.

Use to maintain legs in 90 degrees hip flexion and up to 60 degrees of abduction.

Size	Code	Code	Age of Baby
XS	FCDH/02	ACDH/02	0 – 3 Months
S	FCDH/03	ACDH/03	3 – 6 Months
M	FCDH/05	ACDH/05	6 – 12 Months
L	FCDH/07	ACDH/07	12 – 18 Months
XL	FCDH/08	-	18 + Months



Index

Hand Therapy

Code	Name	Pg
ATW	Abducted Thumb Wrap.....	021
ASX	Air SX Wrist Brace	032
LPS	Air X Low Profile Spica.....	031
AWB	Air X One Piece Wrist.....	015
ATRS	Air X Thumb Restriction Splint	014
ATS	Air X Thumb Spica.....	015
AFS	Air X Two Piece Wrist.....	015
AWT	Air X Wrist Thumb.....	014
100302P	Aluminium Finger Splints.....	038
BWTB	Bedford Wrist Thumb Brace	020
BWB	Beta Wrist Brace	008
SSX	Black Seal SX Wrist Brace.....	032
CAP/5	Cambridge Finger Splint.....	036
CSX	Carbon SX Wrist Brace	032
COT	Carpotec Carpal Brace	034
DFP	Deltaform Plus Wrist Brace	009
DFS	Deltaform Wrist Brace	009
CAP/3	Double Armchair.....	036
DBF	Double Bedford Finger Stall	037
ESNBW	Economy Standard Neobrace (3mm thick).....	010
ESNB	Economy Standard Neobrace.....	010
ETSMLF	Elastic Thumb Spica - Modified	022
ETSLF	Elastic Thumb Spica.....	022
WTE	Elastic Wrist Thumb Brace	019
FSX	Freedom SX Wrist Brace	010
FWB	Freedom Wrist Brace.....	010
JRS	Functional Resting Splint (Padded).....	033
JRSE	Functional Restraining Splint (Easy Liner).....	033
JRS/TP	Functional Restraining Splint (Thermoplastic Insert).....	033
IWB	Intrinsic Wrist Brace	017
JBPM	Jura Basic Mallet Finger Splints.....	038
JBWB	Jura Basic Wrist Brace Elastic	034
JRNWB	Jura Basic Wrist Brace Neoprene	033
JBWTN	Jura Basic Wrist Thumb Neoprene.....	031
JRWB	Jura Black Wrist Brace (Regular)	028
JRWBL	Jura Black Wrist Brace (Long).....	028
JRWTB	Jura Black Wrist/Thumb Brace.....	028
JURACMC	Jura CMC Thumb Brace	035
JLFW	Jura Latex Free Wrist Brace.....	013
JPC	Jura Palm Protector.....	017
JP01	Rainbow Putty.....	031
JUDSE	Jura Ulnar Deviation Support Extended.....	024
JUDS	Jura Ulnar Deviation Support.....	024
JUWBU	Jura Universal Wrist Brace Unbound	030
JUWB	Jura Universal Wrist Brace.....	030
JUWTB	Jura Universal Wrist Thumb Brace.....	033
JWB	Jura Wrist Brace.....	016
JWTL	Jura Wrist Thumb (Long).....	016
JWT	Jura Wrist Thumb.....	016
JLWW	Juraprene Long Wrist Wrap - Sizes.....	025
JSWWU	Juraprene Short Wrist Wrap - Universal	025
JTWLW	Juraprene Wrist Thumb Wrap - Long	025
JTWWS	Juraprene Wrist Thumb Wrap	024
BLTS	Lycrafleece Thumb Spica	012
BLWB	Lycrafleece Wrist Brace	012
MTW	Mini Thumb Wrap	030
NNWB	Neobrace PLUS Wrist Brace.....	030
LNB	Neobrace Wrist Brace (Long)	013
SNB	Neobrace Wrist Brace (Standard).....	013
NTS	Neoprene Thumb Spica.....	034
WTO	Neoprene Wrist Thumb Brace	020
NUF	Nuform Wrist Brace.....	009
P1008-K	Oval Eight Finger Splint (Complete Kit)	038
P1008-S	Oval Eight Finger Splint (Sizing Kit)	038
MFS	Plastic Mallet Splint	037
LPS	Pollo Form Spica - Long.....	023
SPF	Pollo Form Spica - Standard	023
PFTB	Poroflex Pull & Thumb.....	011
PFWB	Poroflex Pull	011
PRHS	Pro-rheuma Thumb Spica	023
PRHW	Pro-rheuma Wrist Brace (Beige)	012
PRHWB	Pro-rheuma Wrist Brace (Black).....	012
PRHT	Pro-rheuma Wrist Thumb Brace (Beige).....	013
PRHTB	Pro-rheuma Wrist Thumb Brace (Black).....	013
PDRF	Procool D-ring Wrist Splint (Firm)	027
PDR	Procool D-ring Wrist Splint (Standard).....	027
PDWT	Procool Deluxe Wrist Thumb Brace.....	027
PTRSB	Procool Thumb Restriction Splint (Beige)	026
PTRS	Procool Thumb Restriction Splint (Black)	026
PWS	Procool Wrist Support	026
PWTRS	Procool Wrist Thumb Restriction	027
PWTS	Procool Wrist Thumb Support.....	032
TS	Promedics Thumb Shield	022
CMC	Push CMC.....	035
RWB	Reversible Wrist.....	029
CAP/4	Simpson Modified Armchair	037
CAP/2	Single Armchair.....	036
SWB	Standard Wrist Brace.....	008
TUJ	Thumbs Up.....	017
CAP/8	Triple Armchair	036
TPNS	Two Piece Neoprene Spica	030
TRSV	Ventilated Thumb Restriction Splint	029
WTRSV	Ventilated Thumb Restriction Splint	029
VNB	Ventoprene Neobrace	018
VWNB	Ventoprene Wrist Brace (Beige)	018
R/VWNB	Ventoprene Wrist Brace (Blue).....	018

WFS	Whitechapel Finger Splints.....	037
WNB	Wrightington Neobrace.....	018
WWB	Wrightington Wrist Brace.....	008
WTN	Wrist Thumb Brace.....	020
WTW	Wrist Thumb Wrap.....	021
WW	Wrist Wrap.....	021
UDSA	Wrist/Ulnar Deviation Support.....	019
WL	Wristlet.....	019

Spinal Support

Code	Name	Pg
983119	Aspen Cervical Collar Large Back Panel.....	049
983110	Aspen Cervical Collar Regular.....	049
983108	Aspen Cervical Collar Short.....	049
983118	Aspen Cervical Collar Small Back Panel.....	049
983120	Aspen Cervical Collar Standard Back Panel.....	049
983112	Aspen Cervical Collar Tall.....	049
983114	Aspen Cervical Collar X-Tall.....	049
983130	Aspen Cervical Regular Collar Set with Replacement Pads.....	049
983128	Aspen Cervical Short Collar Set with Replacement Pads.....	049
983132	Aspen Cervical Tall Collar Set with Replacement Pads.....	049
983134	Aspen Cervical X-Tall Collar Set with Replacement Pads.....	049
983159	Aspen Collar Replacement Pads.....	049
983161	Aspen Collar Small Back Panel Replacement Pad.....	049
983330	Aspen CTO Control Clamps (2).....	053
983321	Aspen CTO Large Back Panel.....	053
983311	Aspen CTO Regular Front Panel.....	053
983310	Aspen CTO Regular.....	053
983320	Aspen CTO Replacement Pads.....	053
983309	Aspen CTO Short Front Panel.....	053
983308	Aspen CTO Short.....	053
983318	Aspen CTO Small Back Panel.....	053
983319	Aspen CTO Standard Back Panel.....	053
983313	Aspen CTO Tall Front Panel.....	053
983312	Aspen CTO Tall.....	053
983315	Aspen CTO X-Tall Front Panel.....	053
983314	Aspen CTO X-Tall.....	053
992505	Aspen Extension Panel.....	057
992505	Aspen Extension Panel.....	059
992505	Aspen Extension Panel.....	075
991035	Aspen LSO / TLSO Replacement Pads.....	071
990021	Aspen LSO LoPro Large.....	071
990025	Aspen LSO LoPro Replacement Pads.....	071
990020	Aspen LSO LoPro Small.....	071
990027	Aspen LSO LoPro X-Large Replacement Pads.....	071
990022	Aspen LSO LoPro X-Large.....	071
990041	Aspen LSO Short / Large.....	071
990040	Aspen LSO Short / Small.....	071
990061	Aspen LSO Tall / Large.....	071
990060	Aspen LSO Tall / Small.....	071
580104	Aspen Lumbar Large.....	073
580103	Aspen Lumbar Medium.....	073
580102	Aspen Lumbar Small.....	073
580105	Aspen Lumbar X-Large.....	073
580101	Aspen Lumbar X-Small.....	073
580106	Aspen Lumbar XX-Large.....	073
983199	Aspen Pediatric Collar Replacement Pads.....	051
501322	Aspen Sierra Adjustable Extension Panel.....	069

501111	Aspen Sierra Adjustable LSOI.....	069
501221	Aspen Sierra Adjustable TLSOI.....	069
501321	Aspen Sierra Adjustable Upgrade Kit.....	069
983125	Aspen Sierra Universal Collar Adult.....	055
983150	Aspen Sierra Universal Collar Paediatric.....	055
991031	Aspen TLSO Short / Large.....	071
991030	Aspen TLSO Short / Small.....	071
991071	Aspen TLSO Tall / Large.....	071
991070	Aspen TLSO Tall / Small.....	071
993905	Aspen Tri-Point FSO Dual Strut Configuration Kit.....	075
993911	Aspen Tri-Point FSO One Size Adjustable.....	075
993061	Black Extension Panel.....	065
993201	Black Large RAP (Rigid Anterior Panel).....	065
993200	Black Regular RAP (Rigid Anterior Panel).....	065
BC	Budget Collar.....	080
CBB	Chairback Brace.....	076
992212	Contour Back Panel 12" Back Panel.....	067
992215	Contour Back Panel 15" Back Panel.....	067
994352	Contour TLSO SPK Complete Large.....	066
994351	Contour TLSO SPK Complete Medium.....	066
994353	Contour TLSO SPK Complete X-Large.....	066
994354	Contour TLSO SPK Complete XX-Large.....	066
994465	Contour TLSO SPK Upgrades Large.....	066
994460	Contour TLSO SPK Upgrades Medium.....	066
994455	Contour TLSO SPK Upgrades Small.....	066
994260	Contour TLSO Upgrades Medium Chest Panel.....	066
994151	Contour TLSO Upgrades Short Shoulder Strap.....	066
994152	Contour TLSO Upgrades Tall Shoulder Strap.....	066
LBF/12	Deep Lumbarform Belt.....	078
36-2106	Evergreen 621 SI Belt Large.....	063
36-2104	Evergreen 621 SI Belt Medium.....	063
36-2102	Evergreen 621 SI Belt Small.....	063
FC	Halton Collar.....	080
993740	Horizon 456 TLSO.....	061
993710	Horizon 627 Lumbar.....	061
993730	Horizon 637 LSO.....	061
993745	Horizon Extension Panels.....	061
LBS/9	Lumbarform Belt with Straps.....	078
TO1303N	MP45 Plus Three Point Hyperextension Brace.....	076
TO410	Nhyco Cruciform Hyperextension Orthosis.....	076
983301	PD CTO with PD4 and PD5 Front.....	053
983192	PD1 Set with Replacement Pads.....	051
983190	PD1.....	051
983197	PD2 Set with Replacement Pads.....	051
983195	PD2.....	051
983126	PD3 Set with Replacement Pads.....	051
983116	PD3, PD4, PD5 Back Panel.....	051
983106	PD3.....	051
983122	PD4 Set with Replacement Pads.....	051
983102	PD4.....	051
983124	PD5 Set with Replacement Pads.....	051
983104	PD5.....	051
PCC	Philly Collar.....	080
TB	Pro Tem Belt.....	077
PF/8	Proflex Belt (Flexible Panel).....	078
PFN/8	Proflex Belt (No Panel).....	078
PFM	Proform Belt.....	077
993031	QuikDraw PRO Black Large.....	065
993021	QuikDraw PRO Black Medium.....	065
993011	QuikDraw PRO Black Small.....	065
993041	QuikDraw PRO Black X-Large.....	065
993081	QuikDraw PRO Black X-Small.....	065
993091	QuikDraw PRO Black XX-Large.....	065

993030	QuikDraw PRO White Large	065
993020	QuikDraw PRO White Medium.....	065
993010	QuikDraw PRO White Small.....	065
993040	QuikDraw PRO White X-Large	065
993080	QuikDraw PRO White X-Small.....	065
993090	QuikDraw PRO White XX-Large	065
SLS	Sacro Lumbar Support	079
LBF/9	Standard Lumbarform Belt.....	078
992710	Summit 456 TLSO One Size Adjustable	059
992310	Summit 637 LSO One Size Adjustable	059
60002	TLSO Back Brace	079
ABW	Universal Abdominal Binder.....	079
UTP/10	Universal Two Piece Collar.....	080
993644	Vista 464 TLSO Upgrade Kit	057
993640	Vista 464 TLSO	057
984028	Vista Collar Large Back Panel Replacement Pad	045
984020	Vista Collar Replacement Pads.....	045
984002	Vista Collar Set with Replacement Pads	045
984000	Vista Collar	045
984504	Vista CTO Upgrade Kit • Anterior.....	047
984510	Vista CTO.....	047
984520	Vista CTO4 / Vista CTO Replacement Pads.....	047
984554	Vista CTO4 Upgrade Kit • Posterior	047
984550	Vista CTO4	047
984620	Vista ICU Back Panel Replacement Pad	045
984600	Vista ICU Back Panel.....	045
984220	Vista MultiPost Collar Replacement Pads.....	045
984202	Vista MultiPost Collar Set with Replacement Pads	045
984200	Vista MultiPost Collar.....	045
984252	Vista MultiPost Therapy Collar Replacement Pads	045
984257	Vista MultiPost Therapy Collar Replacement Therapy Pack.....	045
984250	Vista MultiPost Therapy Collar.....	045
984120	Vista TX Collar Replacement Pads	045
984102	Vista TX Collar Set with Replacement Pads.....	045
984100	Vista TX Collar.....	045
984008	Vista Upper Spine Large Back Panel	045

Upper Limb & Shoulder

Code	Name	Pg
AES	Arm Elevator Sling.....	083
ASCE	Ascender Telescopic ROM Elbow Brace.....	089
BOS	Boscombe Sling.....	084
AESB	Bradford Arm Elevator Sling	083
CLS	Clavicle Splint	088
CTS/10	Cubital Tunnel Splint.....	089
DSI	Deluxe Shoulder Immobiliser	083
ESI	Economy Shoulder Immobiliser.....	083
EAS	Evans Abduction Sling	087
HEW	Hinged Elbow Wrap	090
JASSP/10	Jura Abduction Sling PLUS.....	086
JASS	Jura Abduction Sling	087
JUC	Jura Clasp.....	090
JCS/10	Jura Clavicle Splint.....	088
JSI	Jura Shoulder Immobiliser	086
JSL/10	Jura Shoulder Lock.....	086
NSB	Neoprene Shoulder Brace	088
OES	Orthoprene Elbow.....	090
TAU/P	Paediatric Trauma Sling	085
RVP	Replay Tennis Elbow Support.....	091
SS	Sling Strap	091

SSW	Sling & Swathe	084
SCF	Subluxation Cuff	088
TEP	Tennis Elbow Strap.....	090
TAU/10	Trauma Sling	085
UCS/10	Universal Combo Sling	087
LSS/25	Universal Lancaster Shoulder Sling.....	082
UPS	Universal Pro Sling.....	084
VPS	Vari Pad Tennis Elbow Support.....	091

Lower Limb

Code	Name	Pg
ARK	Air X Long ROM Knee Brace	98
AMP	Air X Polycentric Knee Brace	98
ASK	Air X SK Knee Brace	98
AVK	Air X VK Long Knee Brace	99
ASC/10	Ascender Telescopic ROM knee Brace	93
ASCFP/10	Ascender Telescopic ROM knee Brace with Foam ... Pads	93
BN	Buttress Knee Sleeve.....	103
CNI	Canvas Knee Immobiliser.....	97
EN	Economy Knee Sleeve	102
GN/08	Genurange Lite	95
GN/0510	Genurange Universal	94
GN/04	Genuranger	94
GN/0310	Genurehab	95
HN	Hinged Knee Sleeve	103
JTI	Jura 3 Panel Immobiliser	99
BTPK	Jura Knee Immobiliser	96
OPK	Jura One Piece Immobiliser	97
BLPL	Lycrafleece Patella Lux	101
NS	Neosplint	97
NSK	Nuffield SK Knee Brace	100
PKS/20	Post-OP Knee Splint	96
SK	SK Knee Brace	100
SP	Spiral Knee Sleeve	102
ST	Stabilised Knee Sleeve	103
SNO	StandardKnee Sleeve – Open Patella	102
TPK	Three Panel Knee Immobiliser.....	96
VSKP	Ventilated SK Knee Plus	100
VK	VK Knee Brace	101

Ankle & Footwear

Code	Name	Pg
KTAF01/M	AFO-1 Classic Medial/KT AFO-1 Classic Medial EasyFit.....	112
KTAF0/L	AFO-2 Classic Lateral/KT AFO-1 Classic Lateral EasyFit.....	112
PRAFKT3	AFO-3 Flexible/KT AFO-3 Flexible EasyFit	113
PREGFIX/AFO	AFO-4 Anterior Carbon AFO	113
AGA/10	Air/Gel Ankle.....	114
LEB	Ankle Contracture AFO	111
ANK	Anklet.....	116
L-478	Anti-Slip Pad.....	127
EWB	Barnsley Ankle Brace.....	114
BPO	Basic Post-Op Shoe.....	122
BFS	Boot Foot Support.....	109
L-45	Cap	126
L-DFS	Dress Shoe Comfort Inserts.....	123
L-46	Elastic Sheath.....	126
EPCB	Epoda Dual Heel Cups.....	125
EPB	Epoda Dual Heel Pads	125
F8LA	Figure of 8 Ankle Brace.....	115
821	Fixed Ankle	114

HVN	Hallux Valgus Night Splint.....	115
HUP/10	Heal-Up Protection Boot	110
HRS	Heel Relief Shoe	119
HSCS	Hi-Sided Cast Sandal.....	120
JWW	Jura Achilles Wedges.....	106
JWA-ES	Jura Air Shield Walker Extension Straps	106
JWAS	Jura Air Shield Walker Short.....	106
JWA	Jura Air Shield Walker	106
JWFS*	Jura Fixed Walker Short with Air Liner.....	107
JWFS	Jura Fixed Walker Short.....	107
JWF*	Jura Fixed Walker with Air Liner	107
JWF	Jura Fixed Walker	107
HP-OP	Jura High Post-Op Shoe.....	117
P-OP*	Jura Post-Op Shoe.....	117
825	Leaf Spring Orthosis	111
LAB	Lincoln Ankle Brace	116
M.A1028	Malleostong.....	115
811	Medi Shoe	120
OWH	Off-Loader Heel Shoe	118
-	Off-Loader Insoles.....	118
CS/30	Open Toe Cast Sandal with Buckles	120
CS/48	Open Toe Cast Sandal with Buckles	122
CS/08	Open Toe Cast Sandal	122
PB	Podus Boot.....	108
L-PQ6	PQ - 6mm Heel Lift.....	124
L-PQ9	PQ - 9mm Heel Lift.....	124
ASSW	Promedics Airstep Walker (Short).....	105
ASW	Promedics Airstep Walker (Tall)	105
ASAW	Promedics Airstep Walker Achilles Wedges	105
445	Rigid Post-Op Shoe.....	123
JWR	ROM Walker	108
L-495	Silicone Gel Discs.....	127
L-43*	Silicone Gel Hallux Valgus Protector.....	127
L-47*	Silicone Gel Metatarsal Pad with Elastic Band	126
L-41*	Silicone Gel Retractor/Protector.....	127
L40*	Silicone Heel Cup with Soft Heel Spot.....	125
L-SRI*	Silicone Pro Insole.....	125
L-SHPS	Star Heel Pad.....	124
L-SSPS	Star Spur Pad.....	124
RPO	Superstep Post-Op Shoe.....	121
M.9050	Swedish Ankle Foot Orthosis.....	111
L-43*	Toe Retractor	127
L-44*	Tubular Protection A.....	126
L-44*	Tubular Protection B.....	126
UFA	Universal Foam Ankle	114
UFO	Universal Plantar Fasciitis Orthosis	109
JSHV	Universal Soft Hallux Valgus	115
ASWSS	Walker Boot Socks Short.....	106
ASWS	Walker Boot Socks Tall.....	106

Maternity Support

Code	Name	Pg
DMB	Deluxe Maternity Belt.....	131
DUB	Duofit Maternity Belt.....	131
MBAS	Maternity Belt Style A.....	130
MBB	Maternity Belt Style B.....	130
MBC	Maternity Belt Style C.....	130
SPD	Pelvic Instability Belt.....	129
TROCH	Trochanteric Belt.....	129
WMB	Wiltshire Maternity Belt.....	129

Fracture Bracing & Splinting

Code	Name	Pg
DHC	Derby Heel Cup	136
EHS	Extended Humeral Fracture Brace Sleeve Only.....	133
FEM	Femoral Sock.....	136
GNF	Genurange Fracture Bracing Hinges	135
GNP/10	Genurange Replacement Headplates	135
GEN/CS	Genurange Replacement Hinge Cover Screws	135
GNS/10	Genurange Replacement Screws	135
GON/10	Goniometer	136
HS/10	Humeral Sock.....	136
HSO/10/K	Humeral Stabilising Orthosis Full Kit.....	134
HSO/10	Humeral Stabilising Orthosis Sleeve Only.....	134
JPR*	Juraprene Splinting Material	140
LS	Loop Strap	138
32503*	Orfit Dry Heater.....	141
PHJ/10	Promedics Hip Joint.....	134
SA	Self Adhesive Hook/Loop	140
NC/79615	SL Iron	138
SS	Sling Strap	133
PSH	Soft Humeral.....	134
WTBATH	Splinting Baths	141
HAL	Standard Hook/Loop	140
SHC	Standard Humeral Fracture Brace Full Kit	133
SHS	Standard Humeral Fracture Brace Sleeve Only.....	133
HT/GUN	Steinel Heat Gun	138
TS	Terry Strap	138
BK/10	Tibial Sock	136
JBBS/10	Universal Hip Abduction Belt with Shoulder Straps	134
JBB/10	Universal Hip Abduction Waist Belt.....	134

Hot & Cold Therapy

Code	Name	Pg
HCP	Hot/Cold Reusable Packs.....	144
JHC*	Jura Hot/Cold Ankle Wrap	144
JHC*	Jura Hot/Cold Elbow Wrap.....	144
JHC*	Jura Hot/Cold Low Back Wrap.....	144
CTKW/10	Universal Cold Therapy Wraps (Knee).....	143
CTSW/10	Universal Cold Therapy Wraps (Shoulder).....	143

Paediatric

Code	Name	Pg
ACDH	CDH Splint (Adjustable)	154
FCDH	CDH Splint (Standard)	154
BFB	Children's Buckle Fracture Brace.....	149
CDRL	Children's D Ring Brace Long	149
CDRS	Children's D Ring Brace Short.....	149
FPI	Fun Pattern Arm & Knee Immobiliser	151
JWTO	Junior Wrist/Thumb Brace.....	148
CS	Paediatric Cast Sandal	152
PEW	Paediatric Elastic Wrist Brace	147
PEWT	Paediatric Elastic Wrist/Thumb Brace	147
HSCSL	Paediatric Hi-Sided Cast Sandal.....	152
LSSB/10	Paediatric Lancaster Shoulder Brace (Baby).....	150
LSSB	Paediatric Lancaster Shoulder Brace (Paediatric).....	150
PNW	Paediatric Neoprene Wrist Brace.....	147
PNWT	Paediatric Neoprene Wrist/Thumb Brace	148
JBC	Paediatric Soft Foam Collar	153
PLTS	Paediatric Thumb Spica.....	146
PUDS	Paediatric Thumb Splint.....	150
TAU/P	Paediatric Trauma Sling	151
PUDS	Paediatric Ulnar Deviation Support	150

PJWF	Paediatric Walker.....	153
PLWB	Paediatric Wrist Splint.....	146
PLWTB	Paediatric Wrist/Thumb Splint.....	146
PAV	Pavlik Harness.....	154
TBP	Protem Belt	154
TWB	Tracy Brace	148





Terms & Conditions

Promedics Orthopaedics Limited Terms and Conditions for the supply of goods to trade customers

1. Definitions

1.1 "Contract" means the contract for the sale and purchase of Goods in accordance with these Terms. "Customer" means the person, firm or company whose order for Goods is accepted by Promedics in accordance with these Terms. "Goods" means the goods (including any instalment of the goods) which Promedics is to supply in accordance with these Terms. "Promedics" means Promedics Orthopaedics Limited, (company number 06455477) whose registered office is at Suite A, 6 Honduras Street, London, EC1Y 0TH and whose main trading address is Port Glasgow Industrial Estate, Gareloch Road, Port Glasgow, Scotland, PA14 5XH, and whose registered VAT number is 927148511. "Terms" means the standard terms of sale set out in this document as amended from time to time and (unless the context otherwise requires) includes any special terms agreed in writing between the Customer and Promedics. "Website" means the website operated by Promedics with web address <https://promedics-healthcare-professionals.co.uk/>.

1.2 A reference in these Terms to a provision of a statute shall be construed as a reference to that provision as amended, re-enacted or extended at the relevant time. The headings in these Terms are for convenience only and shall not affect their interpretation.

2. Basis of the sale

2.1 Promedics shall sell and the Customer shall purchase the Goods in accordance with these Terms. These Terms apply to the Contract to the exclusion of any other terms that the Customer seeks to impose or incorporate, or which are implied by trade, custom, practice or course of dealing.

2.2 No variation to these Terms shall be binding unless agreed in writing between an authorised representative of Promedics and the Customer.

2.3 Any typographical, clerical or other error or omission in any sales literature, quotation, price list, acceptance of offer, invoice or other document or information issued by Promedics (including information on the Website) shall be subject to correction without any liability on the part of Promedics.

3. Orders and acceptance

3.1 Prior to placing orders, the Customer shall be required to complete an application form provided by Promedics. The Customer shall ensure that all information provided to Promedics on such application form is complete and accurate in all respects.

3.2 The Customer may place an order by telephone or email or, where agreed between Promedics and the Customer, via Electronic Data Interchange (EDI). Each order is an offer by the Customer to buy the Goods specified in the order subject to these Terms.

3.3 Where Promedics acknowledges receipt of an order, such acknowledgement shall not constitute acceptance of the order. No order submitted by the Customer shall be deemed to be accepted by Promedics, and the Contract shall not come into existence, unless and until:

3.3.1 in the case of custom-made Goods, acceptance is confirmed in writing by Promedics's authorised representative by issuing an order confirmation; or

3.3.2 in the case of Goods that are not custom-made, acceptance is confirmed by Promedics despatching the Goods, or by Promedics informing the Customer that the Goods are ready for collection (where applicable), or by Promedics's authorised representative issuing

an order confirmation, whichever is the earlier.

3.4 Each order accepted by Promedics shall constitute a separate Contract, formed in accordance with these Terms.

3.5 The Customer shall be responsible for ensuring that any order and any specification or design (where applicable) submitted by the Customer is complete and accurate.

3.6 Where the Customer's order refers to Goods which are unavailable, or refers to prices of Goods or delivery charges which are incorrect or not then current, Promedics shall contact the Customer and agree an amended order, after which Promedics's authorised representative will confirm acceptance of such order in writing by issuing an order confirmation, at which point the Contract shall come into existence.

3.7 The Customer shall be responsible to Promedics for ensuring the accuracy of the terms of any order confirmation, and shall carefully check any order confirmation provided by Promedics on receipt thereof and promptly inform Promedics of any discrepancies therein. Failure to inform Promedics promptly (and in any event prior to despatch of the Goods or any instalment thereof) of any discrepancies within an order confirmation shall be deemed to be acceptance by the Customer of the accuracy of such order confirmation.

3.8 The quantity, quality and description of the Goods shall be as set out in Promedics's despatch note.

3.9 Promedics reserves the right to make any changes in the specification of the Goods, which do not materially affect their quality or performance, or which are required to ensure they comply with applicable safety, statutory or regulatory requirements.

3.10 No order which has been accepted by Promedics may be cancelled or amended by the Customer except with the agreement in writing of Promedics and on terms that the Customer shall indemnify Promedics in full against all loss (including loss of profit), costs (including the cost of all labour and materials used), damages, charges and expenses incurred by Promedics as a result of such cancellation or amendment.

3.11 The Customer acknowledges that Promedics is not obliged to accept orders submitted by the Customer.

4. Price of the goods

4.1 Subject to clause 2.3, the price of the Goods shall be the prices shown on the current price list agreed between the Customer and Promedics. Where no price list has been agreed, or where a price list has expired or been cancelled by Promedics, the prices set out in Promedics's then-current standard price list shall apply. Notwithstanding the foregoing, the price for custom-made Goods shall be the price set out on Promedics's order confirmation.

4.2 In the case of custom-made Goods, the price for the Goods shall be the price agreed between the Customer and Promedics in relation to such custom-made Goods.

4.3 Prices for Goods may change from time to time without notice.

4.4 Notwithstanding condition 4.1, Promedics reserves the right, by giving written notice to the Customer at any time before despatch of the Goods, to increase the price of the Goods or the delivery charges to reflect any increase in the cost to Promedics which is due to any increase in the cost of delivery charges, a change in delivery dates, quantities or specifications for the Goods which is requested by the Customer, or any delay caused by any of the Customer's instructions or failure by the Customer to give Promedics adequate information or instructions.

4.5 Unless otherwise agreed between the Customer and Promedics in writing, the Customer shall be liable to pay the delivery charges (which may include the cost of insurance where applicable) in respect of each shipment of the Goods. Such delivery charges shall be as set out on the current price list agreed between the Customer and Promedics, or where such charges are not set out on such price list, as set out on the order confirmation, or where such charges are not set out on the order confirmation, in accordance with Promedics standard delivery charges in force from time to time. Any delivery information set out on the Website (including details of delivery charges) is indicative only and shall not be binding on Promedics.

4.6 If the Customer orders Goods for delivery to a destination outside of the UK, such order may be subject to import duties and taxes which are applied when the delivery reaches that destination. Please note that Promedics has no control over these charges and cannot predict their amount.

4.7 Unless otherwise agreed between the Customer and Promedics in writing, the Customer is responsible for payment of any such import duties and taxes. Please contact your local customs office for further information before placing an order.

4.8 The Customer must comply with all applicable laws and regulations of the country for which the Goods are destined. Promedics will not be liable or responsible if the Customer breaks any such laws or regulations.

4.9 The price of the Goods is exclusive of any applicable value added tax, which, where applicable, the Customer shall be additionally liable to pay to Promedics.

5. Terms of payment

5.1 Subject to any special terms agreed in writing between the Customer and Promedics, Promedics will raise an invoice at or after the time of despatch of the Goods, or the time at which Promedics notifies the Customer that the Goods are ready for collection, as the case may be.

5.2 Unless otherwise agreed between the Customer and Promedics in writing, the Customer shall pay the full amount invoiced to it by Promedics within 30 days of the date of Promedics's invoice and in the currency shown on the invoice. Promedics shall be entitled to recover the price of the Goods and delivery charges (where applicable) notwithstanding that delivery of the Goods may not have taken place.

5.3 Notwithstanding condition 5.1 and condition 5.2, Promedics may require a deposit or payment in full from Customers prior to the despatch or making available of Goods. In such cases, Promedics will issue an invoice to the Customer which will be payable immediately on receipt.

5.4 All payments shall be made by the Customer without any deduction or set off by electronic transfer into a bank account nominated by Promedics from time to time. The time of payment of the price shall be of the essence of the Contract.

5.5 If the Customer fails to make any payment on the due date then, without limiting any other right or remedy available to Promedics, Promedics may:

5.5.1 terminate the Contract or suspend any further deliveries to the Customer;

5.5.2 appropriate any payment made by the Customer to such of the Goods (or the Goods supplied under any other contract between the Customer and Promedics) as Promedics may think fit (notwithstanding any purported appropriation by the Customer);

5.5.3 withhold and/or offset amounts due from Promedics to the Customer against amounts due from the Customer to Promedics under any other contract between Promedics and the Customer (notwithstanding any terms of the Customer to the contrary); and

5.5.4 charge the Customer interest (both before and after any judgment) on the amount unpaid, at the rate of 3 per cent per annum above the Bank of England's base rate from time to time, but at 3 per cent per annum for any period when that base rate is below zero per cent, until payment in full is made (a part of a month being treated as a full month for the purpose of calculating interest). The parties agree that this constitutes a substantial remedy in terms of the Late Payment of Commercial Debts (Interest) Act 1998.

5.6 Where Promedics terminates the Contract in accordance with condition 5.5.1, the Customer shall immediately pay to Promedics all of Promedics's outstanding unpaid invoices and interest and, in respect of Goods supplied but for which no invoice has been submitted, Promedics shall submit an invoice, which shall be payable by the Customer immediately on receipt.

5.7 Where Promedics has agreed to provide the Customer with a limited credit facility, Promedics may refuse to accept any further orders from the Customer where the Customer has exceeded the agreed credit facility at the time of such order(s), or where such order(s) may increase the sums due from the Customer beyond the agreed credit facility.

5.8 Promedics may withdraw or reduce credit facilities at any time and without notice.

6. Delivery

6.1 Delivery of the Goods shall be made by the Customer or a carrier on behalf of the Customer collecting the Goods at Promedics's premises at any time after Promedics has notified the Customer that the Goods are ready for collection, or, if some other place for delivery is agreed in writing by Promedics, by Promedics delivering the Goods to that place.

6.2 Any dates quoted for delivery of the Goods are approximate only and Promedics shall not be liable for a delay however caused. Time for delivery shall not be of the essence of the Contract. The Customer may specify a preferred delivery date when placing an order, however this preferred delivery date shall not be binding on Promedics. The Goods may be delivered in advance of the quoted or requested delivery date on giving reasonable notice to the Customer.

6.3 Where the Goods are to be delivered in instalments, each delivery shall constitute a separate Contract and failure by Promedics to deliver any one or more of the instalments in accordance with these Terms or any claim by the Customer in respect of any one or more instalments shall not entitle the Customer to treat the Contract as a whole as repudiated.

6.4 If, following acceptance of an order by Promedics, Promedics fails to deliver the Goods (or any instalment) for any reason other than any cause beyond Promedics's reasonable control or the Customer's fault, and Promedics is accordingly liable to the Customer, Promedics's liability shall be limited to the cost to the Customer (in the cheapest available market) of similar goods to replace those Goods not delivered, less the price of those Goods not delivered, where such liability shall be limited to a sum equivalent to the price of the Goods not delivered.

6.5 If the Customer fails to take delivery of the Goods or fails to give Promedics adequate delivery instructions for delivery (otherwise than by reason of any cause beyond the Customer's reasonable control or by reason of Promedics's fault) then, without limiting any other right or remedy available to Promedics, Promedics may:

6.5.1 store the Goods until actual delivery and charge the Customer for the reasonable costs (including insurance) of storage; or

6.5.2 sell the Goods at the best price readily obtainable and (after deducting all reasonable storage and selling expenses) charge the Customer for any shortfall below the price under the Contract.

6.6 Damage on delivery ascertainable on reasonable inspection, and errors in the quantity or specifications of the Goods delivered against the applicable despatch note, must be reported to Promedics and the carrier within 21 days of delivery. Failure to report damage/errors within this period (specifying invoice or despatch note number and date of delivery) may lead to the Customer being charged the full amount of the corresponding invoice or the full price of the Goods (whichever is the higher), and the provisions of condition 6.7 will not apply. Promedics may require photographic evidence of damaged goods, and the Customer shall provide all reasonable assistance in this regard.

6.7 Subject to condition 6.6, Promedics will, at its option, where Goods are accepted by Promedics as damaged, repair or replace the damaged Goods or refund the price of the damaged Goods in full, save to the extent that Promedics holds that such damage is caused by the carrier.

6.8 If requested by Promedics, the Customer shall return the damaged

or excess Goods or those Goods that do not match the specification to Promedics at the Customer's cost. Customer shall return such Goods in a saleable condition, in their original boxes with seals and shrink-wrap intact, and Promedics shall be entitled to charge the Customer the full price of such Goods if they are not returned in accordance with the provisions of this condition 6.8.

- 6.9 Promedics will, at its sole discretion, accept the return of Goods (excluding custom-made Goods) that have been ordered by the Customer in error provided that the Customer notifies Promedics within 21 days of delivery of its wish to return such Goods. Customer shall return the Goods at the Customer's cost, within 14 days of notifying Promedics of its wish to return such Goods, in a saleable condition, in their original boxes with seals and shrink-wrap intact. Promedics shall be entitled to charge the Customer a restocking fee.

7. Risk and property

- 7.1 Except as otherwise agreed between Promedics and the Customer in writing, risk of damage to or loss of the Goods shall pass to the Customer:

7.1.1 in the case of Goods to be delivered at Promedics's premises, at the time when Promedics notifies the Customer that the Goods are available for collection; or

7.1.2 in the case of Goods to be delivered otherwise than at Promedics's premises, at the time of delivery or, if the Customer wrongfully fails to take delivery of the Goods, at the time when Promedics has tendered delivery of the Goods.

- 7.2 Notwithstanding delivery and the passing of risk in the Goods, or any other provision of these Terms, the property in the Goods shall not pass to the Customer until Promedics has received in cleared funds payment in full of the price of the Goods and the delivery charges (if applicable) and of all other Goods agreed to be sold by Promedics to the Customer for which payment is then due.

- 7.3 Until such time as the property in the Goods passes to the Customer, the Customer shall hold the Goods as Promedics's fiduciary agent and bailee, and shall keep the Goods separate from those of the Customer and third parties and properly stored protected and insured and identified as Promedics's property.

- 7.4 Until such time as the property in the Goods passes to the Customer (and provided the Goods are still in existence and have not been resold), Promedics may at any time require the Customer to deliver up the Goods to Promedics and, if the Customer fails to do so forthwith, enter on any premises of the Customer or any third party where the Goods are stored and repossess the Goods. Promedics shall not be liable to the Customer for any damage as a result of the removal or retrieval of the Goods.

- 7.5 The Customer shall not be entitled to pledge or in any way charge by way of security for any indebtedness any of the Goods which remain the property of Promedics, but if the Customer does so all monies owing by the Customer to Promedics shall (without limiting any other right or remedy of Promedics) forthwith become due and payable.

8. Warranties and liability

- 8.1 Promedics provides a warranty that on delivery and for a period of 12 months from delivery, the Goods shall:

8.1.1 subject to condition 9.1, conform in all material respects with their description in the despatch note; and

8.1.2 be free from material defects in design, material and workmanship.

- 8.2 Subject to condition 8.3, if:

8.2.1 the Customer gives Promedics notice in writing within a reasonable time of discovery (specifying invoice or despatch note number and date of delivery) that some or all of the Goods do not comply with the warranty set out in condition 8.1;

8.2.2 Promedics is given a reasonable opportunity of examining the Goods; and

8.2.3 if requested by Promedics, the Customer returns the Goods to Promedics at the Customer's cost,

Promedics will, at its option, where Goods are accepted as defective, repair or replace the defective Goods, or refund the price of the defective Goods in full.

- 8.3 Promedics shall not be liable for breach of the warranty set out in condition 8.1 if:

8.3.1 the Customer makes any further use of the Goods after giving notice to Promedics under condition 8.2, including selling the Goods to a third party;

8.3.2 the defect arises as a result of Promedics following any drawing, design, specification or other request supplied by the Customer;

8.3.3 the Customer or a third party alters or repairs the Goods without Promedics's written consent;

8.3.4 the defect arises because the Customer failed to follow Promedics's oral or written instructions as to the storage, use or maintenance of the Goods or (if there are none) good trade practice;

8.3.5 the defect arises as a result of fair wear and tear, wilful damage, negligence, failure to use the Goods in accordance with instructions provided, or abnormal storage or working conditions;

8.3.6 the defect arises as a result of excessive use or use beyond the normal lifespan or intended period of use of the Goods (for example where Goods are single-use items) or where the Goods are used outside their intended purpose; or

8.3.7 the Goods differ from their description or specification as a result of changes made to ensure they comply with applicable statutory or regulatory requirements.

- 8.4 Except as provided in this condition 8, Promedics shall have no liability to the Customer in respect of the Goods' failure to comply with the warranty set out in condition 8.1.

- 8.5 The terms implied by sections 13 to 15 of the Sale of Goods Act 1979 are, to the fullest extent permitted by law, excluded from the Contract. Except as expressly provided in these Terms, and except where the Goods are sold to a person dealing as a consumer (within the meaning of the Unfair Contract Terms Act 1977), all warranties, conditions or other terms implied by statute or common law are excluded to the fullest extent permitted by law.

- 8.6 These Terms also apply to any repaired or replacement Goods supplied by Promedics to the Customer.

- 8.7 Any Goods replaced will belong to Promedics.

- 8.8 Promedics shall be under no liability under the warranty in condition 8.1 (or any other warranty, condition or guarantee) if the total price for the Goods has not been paid by the due date for payment.

9. Quality/Description

- 9.1 All drawings, descriptive matter, specifications and advertising issued by Promedics and any descriptions or illustrations contained in Promedics's catalogues/brochures or on Promedics's Website are issued or published for the sole purpose of giving an approximate idea of the Goods described in them. They will not form part of the Contract and the Customer accepts that colours reproduced in any electronic image or brochure may differ from the actual Goods.

- 9.2 The Customer shall indemnify Promedics and keep Promedics indemnified in full against all liabilities, costs, expenses, damages and losses (including any direct, indirect or consequential losses, loss of profit, loss of reputation and all interest, penalties and legal and other professional costs and expenses) suffered or incurred by Promedics as a result of or in connection with: (i) any defect in Goods caused in whole or in part by Promedics supplying such Goods in accordance with the Customer's instructions or specification; and (ii) any infringement or alleged infringement of any intellectual property rights or any breach of applicable law or regulation caused in whole or in part by Promedics complying or attempting to comply with the Customer's instructions or specification or the use, manufacture or supply of Goods in accordance with the Customer's instructions or specification. This condition shall survive the termination of the Contract.

10. Limitation of Liability

- 10.1 Subject to condition 6.4 and condition 8, the following provisions set out the entire financial liability of Promedics (including any liability for

the acts or omissions of its employees, agents and sub-contractors) to the Customer in respect of any breach of these Terms and any representation, statement or tortious act or omission including negligence arising under or in connection with the Contract.

10.2 Nothing in these Terms excludes or limits the liability of Promedics for death or personal injury caused by Promedics's negligence, or liability for defective products under the Consumer Protection Act 1987, or for fraudulent misrepresentation.

10.3 Subject to conditions 10.1 and 10.2:

10.3.1 Promedics's total liability in contract, tort (including negligence or breach of statutory duty), misrepresentation, restitution or otherwise, arising in connection with the performance or contemplated performance of the Contract shall be limited to the price of the Goods, except as expressly provided in these Terms; and

10.3.2 Promedics shall not be liable to the Customer for any indirect or consequential loss or damage (whether for loss of profit, loss of business, depletion of goodwill or otherwise), costs, expenses or other claims for consequential compensation whatsoever (howsoever caused) which arise out of or in connection with the Contract.

10.4 Promedics shall not be liable to the Customer or be deemed to be in breach of the Contract by reason of any delay in performing, or any failure to perform any of Promedics's obligations in relation to the Goods, if the delay or failure was due to any cause beyond Promedics's reasonable control.

11. Insolvency of the Customer

11.1 This condition 11 applies if:

11.1.1 the Customer makes a voluntary arrangement with its creditors or (being an individual or firm) becomes bankrupt or (being a company) becomes subject to an administration order or goes into liquidation (otherwise than for the purposes of amalgamation or reconstruction); or

11.1.2 an encumbrancer takes possession of, or a receiver is appointed to, any of the property or assets of the Customer; or

11.1.3 the Customer ceases, or threatens to cease, to carry on business; or

11.1.4 the Customer's financial position deteriorates to such an extent that in Promedics's opinion the Customer's capability to adequately fulfil its obligations under the Contract has been placed in jeopardy; or

11.1.5 Promedics reasonably apprehends that any of the events mentioned above is about to occur in relation to the Customer and notifies the Customer accordingly.

11.2 If this condition applies then, without limiting any other right or remedy available to Promedics, Promedics may cancel the Contract or suspend any further deliveries under the Contract without any liability to the Customer, and if the Goods have been delivered but not paid for the price shall become immediately due and payable notwithstanding any previous agreement or arrangement to the contrary.

12. General

12.1 A notice required or permitted to be given by either party to the other under these Terms shall be in writing addressed to that other party at its main trading address or such other address (including email address) as may at the relevant time have been notified pursuant to this provision to the party giving the notice.

12.2 Any notice shall be deemed to have been received:

12.2.1 if delivered by hand, at the time the notice is left at the proper address;

12.2.2 if sent by pre-paid first-class post or other next working day delivery service, at 9.00 am on the second business day after posting; or

12.2.3 if sent by email, at the time of transmission, or, if this time falls outside business hours in the place of receipt, when business hours resume. In this condition 12.2.3, business hours means 9.00am to 5.00pm Monday to Friday on a day that is not a public holiday in the place of receipt.

12.3 Condition 12.2 does not apply to the service of any proceedings or other documents in any legal action or, where applicable, any arbitration or other method of dispute resolution.

12.4 The Contract is between Promedics and the Customer. No other person, firm or company has any rights to enforce any of its terms.

12.5 No waiver by Promedics of any breach of the Contract by the Customer shall be considered as a waiver of any subsequent breach of the same or any other provision.

12.6 If any provision of the Contract is held by a court or other competent authority to be invalid or unenforceable in whole or in part the validity of the other provisions of the Contract and the remainder of the provision in question shall not be affected.

12.7 Any dispute arising under or in connection with the Contract shall be governed by the laws of England, and the Customer agrees to submit to the exclusive jurisdiction of the English courts.

12.8 Promedics shall be entitled to assign its rights and obligations under these Terms and/or the Contract.

T: +44 (0) 1475 746 400
E: enquiries@promedics.co.uk

www.promedics.co.uk

Promedics Orthopaedics Ltd.
Port Glasgow Industrial Estate
Gareloch Road, Port Glasgow
Scotland PA14 5XH



Life Changing Innovations

